

SFWMD

SFWMD
COMPREHENSIVE
QUALITY ASSURANCE
PLAN

September 24, 1996

WRE INVENTORY #346


SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT

SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT

SOUTH FLORIDA WATER MANAGEMENT DISTRICT
SOUTH FLORIDA WATER MANAGEMENT DISTRICT


Comprehensive Quality Assurance Plan #870166G
for
South Florida Water Management District
1480-9 Skees Road, West Palm Beach, FL 33411
(561) 686-8800

Prepared by
South Florida Water Management District
1480-9 Skees Road, West Palm Beach, FL 33411
(561) 686-8800




Dr. Leslie A. Wedderburn, Department Director
Department of Water Resources Evaluation

10/15/96
Date




Maxine Cheesman, Division Director
Water Quality Monitoring Division

10/16/96
Date



Maria Loucraft-Manzano, Division QA Officer
Water Quality Monitoring Division

10-4-96
Date




Linda Lindstrom, Division Director
Resource Assessment Division

10-15-96
Date




Tim Bechtel, Division QA Officer
Resource Assessment Division

10-15-96
Date




Keith Smith, Division Director
Hydrogeology Division

10/16/96
Date



Carole Milliman, Division QA Officer
Hydrogeology Division


10/16/96
Date



Allan Hall, Department Director
Ecosystem Restoration and Implementation Dept.

10-15-96


Date



Allan Steinman, Division Director
Okeechobee Systems Research

10/15/96


Date



Chuck Hanlon, Division QA Officer
Okeechobee Systems Research Division

10-16-96

Date



Patricia Strayer, Division Director
Kissimmee River Basin Division

10/15/96

Date



Brad Jones, Division QA Officer
Kissimmee River Basin Division

10/15/96


Date



Thomas Fontaine, Division Director
Everglades Systems Research Division

10/15/96

Date



Steve Kelly, Division QA Officer
Everglades Systems Research Division

Date



Nick Aumen, Division Director
Ecologically Engineered Systems Research Division

10/11/96

Date



Kathy Pietro, Division QA Officer
Ecologically Engineered Systems Research Division

10/11/96

Date

TABLE OF CONTENTS

QA Plan Elements

Section	Title	# of Pages	Revision Date
1.0	Title Page	3	09/24/96
2.0	Table of Contents	3	09/24/96
3.0	Statement of Policy	1	09/24/96
4.0	Organization and Responsibility	9	09/24/96
5.0	QA Targets For Precision, Accuracy and Method Detection Limits	7	09/24/96
6.0	Sampling Procedures	34	09/24/96
7.0	Sample Custody	22	09/24/96
8.0	Analytical Procedures	3	09/24/96
9.0	Calibration Procedures and Frequency	8	09/24/96
10.0	Preventative Maintenance	6	09/24/96
11.0	Quality Control Checks, Routines to Assess Precision and Accuracy, and Calculation of Method Detection Limits	9	09/24/96
12.0	Data Reduction, Validation and Reporting	7	09/24/96
13.0	Corrective Action	3	09/24/96
14.0	Performance and Systems Audits	12	09/24/96
15.0	Quality Assurance Reports	2	09/24/96

FIGURES

Figure #	Description	Page #	Revision Date
4.1	SFWMD Organization Chart	4.4	09/24/96
4.2	Department of Water Resources Evaluation Organization Chart	4.5	09/24/96
4.3	Water Quality Monitoring Division Organization Chart	4.6	09/24/96
4.4	Hydrogeology Organization Chart	4.7	09/24/96
4.4	Resource Assessment Organization Chart	4.8	09/24/96
4.6	Departments of Ecosystem Restoration and Regulation Organization Charts	4.9	09/24/96
6.1	Inorganic Surface Water Sampling Procedure	6.17	09/24/96

Figure #	Description	Page #	Revision Date
7.1	Field Trip Check List	7.2	09/24/96
7.2	Sample Submission Diagram for Unfiltered Water	7.3	09/24/96
7.3	Sample Submission Diagram for Filtered Water	7.4	09/24/96
7.4	Chemistry Field Data Log	7.5	09/24/96
7.5	Sample Log-In to LIMS	7.6	09/24/96
7.6	Sample Product List	7.7	09/24/96
7.7	Atomic Absorption Analysis Log	7.9	09/24/96
7.8	Atomic Absorption Sample Log	7.10	09/24/96
7.9	RFA Log	7.11	09/24/96
7.10	Physical Parameters Log	7.12	09/24/96
7.11	Ion Chromatography Log	7.13	09/24/96
7.12	Carbon Analyzer Log	7.14	09/24/96
7.13	Total Coliform Log Sheet	7.15	09/24/96
7.14	Fecal Coliform Log Sheet	7.16	09/24/96
7.15	MPN Log Sheet	7.17	09/24/96
7.16	Fecal Strep Log Sheet	7.18	09/24/96
12.1	Example of Final Report	12.7	09/24/96
14.1	Systems Audit Checklist	14.3	09/24/96

TABLES

Table #	Description	Page #	Revision Date
5.1	Laboratory Quality Assurance Objectives	5.2	09/24/96
5.2	Microbiology Parameters QA Objectives	5.7	09/24/96
5.3	Field Quality Assurance Objectives	5.7	09/24/96
6.1	Sampling Capabilities by Major Category	6.1	09/24/96
6.2	Inorganic Surface Water Sampling Equip.	6.3	09/24/96
6.3	Groundwater Sampling Equipment	6.5	09/24/96
6.4	Soil/Sediment Sampling Equipment	6.7	09/24/96
6.5	Biological Tissues Sampling Equipment	6.7	09/24/96
6.6	Atmospheric Deposition Equipment	6.8	09/24/96
6.7	Autosampler Equipment	6.8	09/24/96
6.8	Low Level Mercury Sampling Equipment	6.9	09/24/96
6.9	Macrobenthic Invertebrate Sampling Equip	6.10	09/24/96
6.10	Holding Time and Preservation for Water Inorganics	6.31	09/24/96
6.11	Holding Time and Preservation for Water Organics	6.32	09/24/96

TABLES

Table #	Description	Page #	Revision Date
6.12	Holding Times and Preservation for Sediments, Soils and Tissues	6.32	09/24/96
6.13	Holding Times and Preservation for Microbiologicals	6.33	09/24/96
6.14	Holding Times and Preservation for Biologicals	6.33	09/24/96
6.15	Reagent and Standard Storage	6.34	09/24/96
8.1	Reagent and Chemical Storage	8.2	09/24/96
9.1	Standard Sources and Preparation	9.4	09/24/96
9.2	Solutions Requiring Standardization	9.5	09/24/96
9.3	Laboratory Instrument Calibration	9.6	09/24/96
9.4	Laboratory Equipment Calibration	9.7	09/24/96
9.5	Field Instrument Calibration	9.8	09/24/96
10.1	Field Equipment Maintenance Schedule	10.2	09/24/96
10.2	Laboratory Equipment Maintenance Schedule	10.3	09/24/96
11.1	Field Quality Control Checks	11.2	09/24/96
11.2	Laboratory Quality Control Checks	11.5	09/24/96
11.3	Microbiological Checks	11.5	09/24/96
11.4	Procedures Used to Determine Precision and Accuracy	11.9	09/24/96
12.1	Formulas Used for Calculations	12.5	09/24/96
13.1	Corrective Actions for the Laboratory	13.2	09/24/96
13.2	Corrective Actions for the Field	13.3	09/24/96
13.3	Corrective Actions Resulting from FQC	13.3	09/24/96

APPENDICES

Appendix ID	Description	# of Pages	Revision Date
A	Color Method SOP	5	01/18/94
B	Silica Method SOP	9	09/26/94
C	Soluble Reactive Phosphate	2	07/14/95
D	Alkaline Phosphatase Method SOP	8	10/01/94
E	Benthic Macroinvertebrate	18	07/14/95
F	Total Inorganic Carbon SOP & Analysis Sheet	7	07/14/95
G	Porewater SOP	24	08/03/95

3.0 STATEMENT OF POLICY

This is a comprehensive quality assurance plan (CQAP) for the sampling and analytical programs of the South Florida Water Management District (SFWMD). These programs encompass surface, estuarine, ground water, atmospheric deposition, biological tissue and sediment monitoring programs. The Water Quality Monitoring Division Laboratory is certified by the Florida Department of Health and Rehabilitative Services as an environmental laboratory capable of performing metals, nutrients, microbiology and general parameter 1 & 2 categories.

This CQAP is intended to be used as a reference and guide for those personnel collecting and analyzing samples and evaluating the results obtained for those analyses. It sets forth the minimum standards to be complied with for these activities and provides a reference for evaluating the procedures used during the time this CQAP is in effect.

The SFWMD is committed to the use of good QA/QC management practices to produce data of a verifiable quality.

4.0 Organization and Responsibility

4.1 Capabilities

The South Florida Water Management District conducts field sampling for surface water, ground water, atmospheric deposition, sediments, tissues and soils. The laboratory is capable of performing analyses for inorganic anions, metals, physical properties, organics (TOC), microbiology and other tests such as chlorophyll and macrobenthic invertebrates.

4.2 Key Personnel

The following are key personnel associated with the collection of samples and the analysis of the samples.

Department Director, Department of Water Resources Evaluation: Responsible for the allocation of resources throughout the department to meet the needs of the SFWMD for sampling and analytical services and championing the quality assurance program.

Department Directors, Department of Ecosystem Restoration and Department of Regulation: Responsible for the allocation of resources throughout the department to meet the needs of the SFWMD for sampling and supporting the quality assurance program within their respective departments.

Division Directors, Hydrogeology, Resource Assessment, Okeechobee Systems Research, Kissimmee Systems Research, Everglades Systems Research, Ecologically Engineered Systems Research, Surface Water Management, Okeechobee Service Center and Field Engineering: Responsible for the allocation of resources, training of personnel, collection of samples and supporting the quality assurance program within their respective division to meet the needs of the SFWMD for sampling services.

Division Director, Water Quality Monitoring (WQM) Division: Responsible for the allocation of resources, training of personnel, collection and analysis of samples, and reporting of results to meet the needs of the SFWMD for analytical services.

Supervising Professional, WQM Division Quality Assurance Officer: Responsible for assuring that the laboratory and field personnel adhere to the approved methods of sample collection, analysis, maintenance of the Comprehensive Quality Assurance Plan, method validation studies, issuance of new methods, the administration of the internal and external laboratory audits, field audits and the review of legislation pertaining to laboratory quality assurance.

Supervising Professional, Chemist: Responsible for the supervision of the laboratory, review of quality control results, review of data, release of samples, training of personnel and adherence to required quality control procedures.

Senior Chemists: Responsible for the supervision of the assigned shift, review of quality control

results, review of data, release of samples, training of personnel and adherence to required quality control procedures.

Senior Environmental Scientist-Special Projects, WQM: Responsible for overseeing organics monitoring for the division.

Staff Environmental Scientists, Water Quality Monitoring Division: Duties include project management, report generation, data review, and collecting samples.

Hydrogeology, Okeechobee Systems Research, Kissimmee Systems Research, Everglades Systems Research, Ecologically Engineered Systems Research, Surface Water Management, Okeechobee Field Station and Field Engineering Division Quality Assurance Officers: Responsible for coordination of all project quality assurance plans and quality assurance (QA) reports for the Division, review of quality control results for the projects, conducting field audits and training of division personnel in quality control procedures. The role of QA officer is separate and distinct from all other responsibilities for any specific project.

Resource Assessment Division Quality Assurance Officer: Responsible for coordination of all Department of Ecosystem Restoration and Resources Assessment Division project quality assurance plans and QA reports, review of quality control results for the projects, conducting field audits and training of Ecosystem Research and Implementation department personnel in quality control procedures. The role of QA officer is separate and distinct from all other responsibilities for any specific project.

Supervising Professional-Environmentalists and Hydrogeologists: Responsible for the supervision of the project managers and field operations supervisors, training of personnel, coordination with divisional QA officers for QA/QC issues, development of sampling networks, review of quality control data and analytical results, and development of research projects.

Supervising Professional-Okeechobee and West Palm Beach Units, Water Quality Monitoring Division: Responsible for supervision of the Water Quality Monitoring Division's Okeechobee and West Palm Beach Sample Collection Units, training of personnel in sampling and quality control procedures, design of sampling networks and research projects, review of quality control data and analytical results, and coordination with divisional QA officer.

Field Operations Supervisor, Hydrogeology Division: Responsible for the allocation of personnel for ground water sampling, drilling and logging of monitor wells, and training of personnel.

Senior Technician Supervisors-Scientific, Water Quality Monitoring Division: Responsible for the allocation of personnel and equipment for surface water sampling, training of personnel in sampling and quality control procedures, review of quality control data, review of analytical results for specified projects.

Staff Programmer Analyst, Water Quality Monitoring Division: Responsible for the maintenance

and integrity of the laboratory information management system, programming computers, testing software packages, and operation of the LIMS.

4.3 Organization Charts

The following charts show the organization of the South Florida Water Management District.

Figure 4.1 is the overall organization of the SFWMD showing the Governing Board, Executive and Deputy Directors and the major departments and offices.

Figure 4.2 is the organization of the Department of Water Resources Evaluation showing the reporting relationships to the Department Director.


Figure 4.3 shows the organization of the Water Quality Monitoring Division through the supervisory levels.

Figure 4.4 shows the organization of the Hydrogeology Division through the supervisory levels.

Figure 4.5 shows the organization of the Resource Assessment Division through the supervisory levels.

Figure 4.6 shows the organization of the Department of Ecosystem Restoration and the Department of Regulation through the supervisory levels for divisions conducting sample collection.

SOUTH FLORIDA WATER MANAGEMENT DISTRICT



Planning	Regulation	Water Resource Evaluation	Ecosystem Restoration Department	Operations & Maintenance	Management Services	Construction & Lead Mgmt.
Director Don Cary Deputy Director Dana Powell	Director Tami Bates Deputy Director Ken Ammon	Director Linda Waddellburn Deputy Director Trevor Campbell	Director Alan Hall Deputy Director Sharon Trout	Director Joseph Schwagert Deputy Director Jennie Hall Deputy Director Victor Powell	Director James Yager Deputy Director Barry Atwood Deputy Director Sandra C. Turnquist Office of Information Systems Director John Lynch	Director William Malone Deputy Director Chuck Ronald Deputy Director Construction Albert Bando

Approvals:
 Executive Office: Michael Slayton

Date: April 1, 1996

Department of Water Resources Evaluation




Figure 4.2 Department of Water Resources Evaluation

Water Quality Monitoring Division




Figure 4.3 Water Quality Monitoring Division Organizational Chart

Hydrogeology Division




Figure 4.4 Hydrogeology Division Organizational Chart




Figure 4.5 Resource Assessment Division Organizational Chart




Figure 4.6 Department of Ecosystem Restoration and Department of Regulation Organizational Chart

5.0 QA TARGETS FOR PRECISION, ACCURACY AND METHOD DETECTION LIMITS

Tables 5.1 and 5.2 present the laboratory quality assurance objectives and Table 5.3 presents the field quality assurance objectives used by SFWMD. The data are generated from historical data collected in the laboratory and field.

Table 5.1 - QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Precision % RSD	Accuracy % Recv	MDL mg/L
Alkalinity	Surface H ₂ O Ground H ₂ O	EPA 310.1	0 - 5.0	NA	1.0
Ammonia	Surface H ₂ O Ground H ₂ O	SM 4500-NH3E	0 - 3.0	94.2 - 106 L	0.015
Ammonia, Un-ionized	Surface H ₂ O Ground H ₂ O	DEP SOP 10/3/83	NA	NA	0.015
Calcium, dissolved	Surface H ₂ O Ground H ₂ O	SM 3120B	0 - 10	80 - 120 M	0.01
Chloride	Surface H ₂ O Ground H ₂ O	EPA 300.A	0 - 2.6	95.4 - 105 M	0.5
Fluoride	Surface H ₂ O Ground H ₂ O	SM 4500F.C	0 - 5.0	92.9 - 103L	0.01
Nitrogen, Organic	Surface H ₂ O Ground H ₂ O	EPA 351.2 - SM 4500NH3H	NA	NA	0.5
Nitrogen, Total Kjeldahl	Surface H ₂ O Ground H ₂ O	EPA 351.2	0 - 8.2	86.9 - 114 M	0.50
Magnesium, dissolved	Surface H ₂ O Ground H ₂ O	SM 3120B	0 - 10	80 - 120 M	0.01
Nitrate + nitrite	Surface H ₂ O Ground H ₂ O	SM 4500NO3F	0 - 3.4	94.6 - 108 M	0.015
Nitrite	Surface H ₂ O Ground H ₂ O	SM 4500NO2B	0 - 2.9	90.2 - 101 M	0.004
Nitrate	Surface H ₂ O Ground H ₂ O	SM 4500NO3F - SM 4500NO2B	NA	NA	0.015

Table S.1 (cont'd) - QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Precision % RSD	Accuracy % Recv	MDL mg/L
Total phosphorus	Surface H ₂ O Ground H ₂ O	SM 4500PF	0 - 10.8	88.1 - 110 M	0.004
Inorganic carbon, total and dissolved, MVP	Surface H ₂ O Ground H ₂ O	SEWMD 3150.2	0 - 15	93 - 117 M	1.0
Organic carbon, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 415.1	0 - 10	87.0 - 114 M	1.0
Orthophosphate	Surface H ₂ O Ground H ₂ O	SM 4500PF	0 - 5.0	92.9 - 109 M	0.004
Potassium, dissolved	Surface H ₂ O Ground H ₂ O	SM 3120B	0 - 10	80 - 120 M	0.01
Residue, filterable	Surface H ₂ O Ground H ₂ O	SM 2540C	0 - 10	NA	13
Residue, nonfilterable	Surface H ₂ O Ground H ₂ O	EPA 160.2	0 - 10	NA	2
Residue, nonfilterable low level	Surface H ₂ O Ground H ₂ O	EPA 160.2	0 - 10	NA	0.4
Residue, volatile	Surface H ₂ O Ground H ₂ O	EPA 160.4	0 - 10	NA	2
Silica, dissolved	Surface H ₂ O Ground H ₂ O	SM 4500SID Modified	0 - 2.8	95.1 - 104 M	1.0
Sodium, dissolved	Surface H ₂ O Ground H ₂ O	SM 3120B	0 - 10	80 - 120 M	0.01
Sulfate	Surface H ₂ O Ground H ₂ O	EPA 300.A	0 - 2.7	94.0 - 104 M	1.0

Table 5.1 (cont'd) - QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Digestion Method #	Precision % RSD	Accuracy % Recv	MDL ug/L
Aluminum, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 202.2	SW - 3015	0 - 10	81 - 126 M	2.5
		SM 3120B	SW - 3015	0 - 3.1	84 - 121 M	4.0
Antimony, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 204.2	SW - 3015	0 - 10	80 - 120 M	2.2
		SM 3120B	SW - 3015	0 - 17	82 - 121 M	2.2
Arsenic, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 206.2	SW - 3015	0 - 5.4	89.6 - 110 M	1.5
		SM 3120B	SW - 3015	0 - 10	80 - 120 M	4.0
Barium, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 208.2	SW - 3015	0 - 10	80 - 120 M	2.0
		SM 3120B	SW - 3015	0 - 5.3	78 - 116 M	0.2
Beryllium, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 210.2	SW - 3015	0 - 10	80 - 120 M	0.2
		SM 3120B	SW - 3015	0 - 2.6	91 - 110 M	0.1
Cadmium, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 213.2	SW - 3015	0 - 4.9	87.0 - 109 M	0.1
		SM 3120B	SW - 3015	0 - 3.8	93.4 - 105 M	0.3
Chromium, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 218.2	SW - 3015	0 - 5.0	80 - 113 M	0.5
		SM 3120B	SW - 3015	0 - 4.4	93.8 - 120 M	1.2
Chromium VI, total and dissolved	Surface H ₂ O Ground H ₂ O	SM 3500CR-D	SW - 3015 SW - 3015	0 - 10	80 - 120 M	3
Copper, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 220.2	SW - 3015	0 - 4.5	83.9 - 111 M	0.5
		SM 3120B	SW - 3015	0 - 5.6	91.4 - 110 M	1.2
Iron, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 236.1	SW - 3015	0 - 9	87.3 - 110 M	4.0
		SM 3120B	SW - 3015	0 - 3.1	93.5 - 115 M	3.0
Lead, total and dissolved	Surface H ₂ O Ground H ₂ O	EPA 239.2	SW - 3015	0 - 6.2	84.3 - 109 M	0.8
		SM 3120B	SW - 3015	0 - 10	90.9 - 111 M	4.0

Table 5.1 (cont'd) - QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Digestion Method #	Precision % RSD	Accuracy % Recv	MDL (ug/L)
Manganese, total and dissolved	Surface H ₂ O	EPA 243.2	SW - 3015	0 - 10	80 - 120 M	1.0
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 2.7	87.9 - 112 M	0.2
Mercury, total and dissolved	Surface H ₂ O	EPA 245.1	SW - 3015	0 - 9.4	91.2 - 119 L	0.2
	Ground H ₂ O		SW - 3015			
Nickel, total and dissolved	Surface H ₂ O	EPA 249.2	SW - 3015	0 - 10	80 - 120 M	0.5
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 3.0	88 - 120 M	0.5
Selenium, total and dissolved	Surface H ₂ O	EPA 270.2	SW - 3015	0 - 10	80 - 120 M	1.0
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 10	80 - 120 M	6.0
Silver, total and dissolved	Surface H ₂ O	EPA 272.2	SW - 3015	0 - 10	80 - 120 M	0.5
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 10	80 - 120 M	1.0
Strontium, total and dissolved	Surface H ₂ O	SM 3111B	SW - 3015	0 - 6.7	94.2 - 110 M	0.2 mg/L
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 10	80 - 120 M	0.3 mg/L
Thallium, total and dissolved	Surface H ₂ O	EPA 279.2	SW - 3015	0 - 10	80 - 120 M	0.5
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 10	80 - 120 M	4.0
Tin, total and dissolved	Surface H ₂ O	EPA 282.2	EPA 282.2	0 - 10	80 - 120 M	0.7
	Ground H ₂ O	SM 3120B	SM 3120B	0 - 10	80 - 120 M	10.0*
Zinc, total and dissolved	Surface H ₂ O	EPA 289.1	SW - 3015	0 - 5.0	95.4 - 106 M	R
	Ground H ₂ O	SM 3120B	SW - 3015	0 - 4.9	82.3 - 122 M	4.0

Table 5.1 (cont'd) - QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Precision % RSD	Accuracy % Recv	MDL
Hardness	Surface H ₂ O Ground H ₂ O	SM 2340B by calculation	NA	NA	5 mg eq CaCO ₃ /L
Alkaline Phosphatase MVP	Surface H ₂ O Ground H ₂ O	SFWMD 3160.1	0 - 9	NA	1 nM/min-ml
Biochemical Oxygen Demand	Surface H ₂ O Ground H ₂ O	EPA 405.1	0 - 10	80 - 120 M	2 mg/L
Carbonaceous Biochemical Oxygen Demand	Surface H ₂ O Ground H ₂ O	AM 5210B	0 - 10	80 - 120 M	2 mg/L
Chemical Oxygen Demand	Surface H ₂ O Ground H ₂ O	EPA 410.1	0 - 10	80 - 120 M	2 mg/L*
Chlorine Residual, total	Surface H ₂ O Ground H ₂ O	EPA 330.4	0 - 10	NA	0.02 mg/l
Chlorophyll	Surface H ₂ O	SM 10200H	0 - 15	NA	1 mg/m ³
Color	Surface H ₂ O Ground H ₂ O	SM 2120B modified	0 - 5.0	NA	1 Pt. Co unit
Macrobenthic Invertebrates	Surface H ₂ O Sediments	SM 10500	NA	NA	NA
pH	Surface H ₂ O Ground H ₂ O	SM 4500H-B	0 - 5.0	NA	0.5 pH unit
Specific Conductance	Surface H ₂ O Ground H ₂ O	SM 2510B	0 - 5.0	NA	0.4 uS/cm
Turbidity	Surface H ₂ O Ground H ₂ O	SM 2130B	0 - 6.5	NA	0.1 NTU

*These are the targeted MDLs based on the instrument manufacturer's claims. Calculated MDLs will be reported in the next revision of the Comprehensive Quality Assurance Plan.

Table 5.2 - MICROBIOLOGY PARAMETERS QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Precision %RSD	Accuracy % Recovery ¹	MDL ²
Total Coliform	Surface H ₂ O Ground H ₂ O	SM 909A/9222B SM 908A/9221B	0 - 8.9 0 - 20	NA	NA
Fecal Coliform	Surface H ₂ O Ground H ₂ O	SM 909C/9222D SM 908C/9221C	0 - 15 0 - 20	NA	NA
Heterotrophic Plate Count	Surface H ₂ O Ground H ₂ O	SM 907A/9215B	0 - 40	NA	NA
Fecal Streptococci	Surface H ₂ O Ground H ₂ O	SM 910B/9230C	0 - 37	NA	NA

¹ There are currently no established accuracy data for relating analytical results to field populations.

² Since the ratio of culturable to unculturable cells depends on the sample, a generic MDL can not be determined.

Table 5.3 - FIELD QUALITY ASSURANCE OBJECTIVES

Component	Matrix	Analytical Method #	Precision % RSD	Accuracy % Recy	MDL
pH	Surface H ₂ O Ground H ₂ O	EPA 150.1	0 - 5	NA	0.5 pH unit
Oxygen, dissolved	Surface H ₂ O Ground H ₂ O	EPA 360.1	0 - 20	NA	0.1 mg/L
Specific conductance	Surface H ₂ O Ground H ₂ O	EPA 120.1	0 - 5	NA	50 uS/cm
Temperature	Surface H ₂ O Ground H ₂ O	EPA 170.1	0 - 3	NA	5°C
Salinity	Surface H ₂ O Ground H ₂ O	SM 2520B	0 - 5	NA	0.1 ppt
Turbidity	Surface H ₂ O Ground H ₂ O	EPA 180.1	0 - 20	NA	1 NTU

6.0 SAMPLING PROCEDURES

6.1 Sampling Capabilities

The sampling capabilities of SFWMD are shown in Table 6.1. The reference used in developing sampling procedures is the EPA Region IV Engineering Support Branch Standard Operating Procedures and Quality Assurance Manual (1991). This reference is available in the field and is referred to in this document as EPA SOP & QAM.

Table 6.1 Sampling Capabilities by Major Category

<u>Matrix</u>	<u>Parameters</u>
Surface Water	Inorganic Anions Metals Physical Properties Organics Extractable Organics Volatile Organics Microbiology Other (Macrobenthic Invertebrates/Chlorophyll)
Ground Water/ Porewater	Inorganic Anions Metals Physical Properties Organics Extractable Organics Volatile Organics Microbiological
Soil/sediments	Inorganic Anions Metals Physical Properties Organics Extractable Organics Volatile Organics
Atmospheric Deposition	Inorganic Anions Metals Physical Properties
Biological Tissues	Inorganic Anions Metals Physical Properties Organics

Samples are collected from the least to the most contaminated areas whenever possible. Time limitations and distance between sites may make this impossible. However, the majority of the District sampling sites are ambient water sources with little variation in concentration levels within a sampling trip. For ground water sampling, if a well is suspected to have free product, then it is not sampled. A new pair of disposable latex/PVC gloves is used at each sampling point for all types of sampling. The preferred order of sample collection is: 1) VOC, 2) POX, 3) TOX 4) TOC, 5) extractable organics, 6) total metals, 7) dissolved metals, 8) microbiological, 8) inorganics, 9) turbidity and 10) macrobenthic invertebrates.

6.2 Field Equipment

6.2.1 Surface Water Sampling Equipment

Table 6.2 lists the equipment used for inorganic surface water quality sampling.

Table 6.2 Inorganic Surface Water Quality Sampling Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>	<u>Notes</u>
Wildco 2.2 liter vertical sampling bottle	acrylic	collection	1
Wildco 2.2 liter horizontal sampling bottle	PVC	collection	1
Wildco 3.2 liter horizontal sampling bottle	PVC	collection	1
5.0 liter bucket	Polyethylene	collection	1
Nalgene sample bottles	Polyethylene	collection	1
Millipore Swinnex filter holder	acrylic/polypropylene	filtration	1
Poretics prefilter	glass fiber	filtration	1
0.45 micron filter	polycarbonate	filtration	1
Millipore 60/150 ml syringe	plastic	filtration	1
Masterflex pump	Not Applicable	collection	1
Pump tubing	C-FLEX	collection	1
Sampling boom	PVC pipe with LDPE sample bottle attached	collection	1,2

¹Not suitable for the collection of organics, extractable organics and VOCs

²Used to hold tubing at desired depth while collecting sample with the pump. Does not contact sample

Field Instrumentation

1. Hydrolab Model 4031
2. Hydrolab Model 4041
3. Hydrolab Surveyor II
4. Hydrolab Surveyor III
5. Solomat WP803
6. YSI Model 58 DO Meter
7. YSI Model 6000UPG Multiprobe - Temperature, pH, Spec. Conductance, DO, Turbidity, EH
8. YSI Model 600XL Multiprobe - Temperature, pH, Spec. Conductance, DO, EH
9. HF Scientific DRT-15CE - Turbidity
10. Secchi Depth Disc
11. Licor spherical quantum sensors

Navigational Aids

1. Apelco DXL 6000 Loran with antenna and battery- for site location
2. USGS Quadrangle maps- for site location
3. Project location maps- for site location
4. WMD low band radio for communication
5. Mark Hurd Areal Photographs-for site location
6. Global Positioning Systems, Pathfinder Basic Plus- for site location

Boating Supplies

1. Safety equipment: life vests, horn, flare kit, oars, fire extinguisher and throwable cushion
2. Spare tire for trailer
3. Engine oil for boat

6.2.2 Ground Water Sampling Equipment

The following is a list of the equipment used for ground water quality sampling.

Table 6.3. Ground Water Sampling Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>	<u>Parameters</u>
Purging Equipment			
1. Gorman Rupp/Centrifugal Pump	Iron/Steel/Rubber	Purge Only	Inorganic/Organic
2. Honda/Centrifugal Pump	Iron/Steel/Rubber	Purge Only	Inorganic/Organic
3. Grundfos Submersible Pump	Stainless Steel/Teflon	Purge Only	Inorganic
4. Suction Hose	Flex PVC	Purge Only	Inorganic
5. Suction Hose	Polypropylene	Purge Only	Inorganic/Organic
6. Drop Pipe	Teflon/Stainless Steel	Purge Only	Inorganic/Organic
7. Check Valve	Teflon	Purge Only	Inorganic/Organic
8. 2 Gallon Bucket	Polyethylene	Purge Only	Inorganic/Organic
9. Lufkin Tapes	Stainless Steel	Prior to Purge	Inorganic/Organic
10. Solinst Electronic Water Level Indicator	Teflon/Stainless Steel	Prior to Purge	Inorganic/Organic

Table 6.3. Ground Water Sampling Equipment (Continued)

Equipment Description	Type of Material	Use	Parameters
1. Sampling Equipment Norwell Bailer & Top/bottom/filter adapter /control-flow bottom	Teflon	Sample	Inorganic/Organic
2. Lanyard	Teflon Coated Stainless Steel	Sample	Inorganic/Organic
Field Filtration Units			
1. QED FF-8200	Filter: Acrylic Copolymer with Polypropylene Body 0.45 micron. 1.0 micron is used for metals collected for permit compliance	Filtration	Inorganic
2. Milyvac Hand Held Vacuum Pump	Tygon Tubing/PVC	Filtration	Inorganic

6.2.3 Soil/Sediment Equipment

The following equipment is used to collect soil/sediment samples.

Table 6.4 Soil/Sediment Sampling Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>
Petite Ponar™	Stainless Steel	Sampling
Scoop/spoon	Stainless Steel	Sampling/ Compositing
Bowl/Tray	Stainless Steel	Compositing/ Homogenizing (not VOCs)
Core	Stainless Steel	Sampling
Core	PVC	Sampling (not for Organics)
Core attachment	PVC	Sampling (not for Organics)
Core stoppers	Rubber	Sampling (not for Organics)

6.2.4 Biological Tissues Sampling Equipment

The following equipment is used to collect biological tissue samples.

Table 6.5 Biological Tissues Sampling Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>
Shears	Stainless Steel	Sampling (above ground biota)
Knife	Stainless Steel	Sampling (above ground biota)
Core	PVC	Sampling (below ground biota)
Opaque Bags	Plastic	Sample storage (biota)
Soil Sieve, 5mm mesh size)	Stainless Steel	Sample processing

6.2.5 Atmospheric Deposition Equipment

The following equipment is used to collect atmospheric deposition samples.

Table 6.6 Atmospheric Deposition Sampling and Processing Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>	<u>Notes</u>
Aerochem Metric Wet/Dry Precipitation Collector	Aluminum	Hold sample collection buckets	1
Aerochem dryfall bucket and lid	PVC	Sample collection	2
Aerochem wetfall bucket and lid	PVC	Sample collection	2
5 liter bucket	Polyethylene	Sample processing	2
Rubber spatula/tweezers/scoop	Plastic	Sample processing	2

Notes:

¹ Equipment will not contact sample.

² For inorganic ions, physical properties and Aluminum and Iron collection/processing only.

6.2.6 Autosampler Equipment

The following types of autosamplers and equipment are used:

Table 6.7 Autosampler Collection and Processing Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>	<u>Notes</u>
American Sigma STREAMLINE Model 700/800 Autosampler	Polyethylene	Hold sample containers	1,4
American Sigma Model 6201 Autosampler	Polyethylene	Hold sample containers	1,3
Sample bottle 5 gal.	Polypropylene	Sample collection	2
Sample bottles 1 liter	Polypropylene	Sample collection	2
Pump tubing/Intake tubing	Silicon/PVC	Sample collection	2

Notes:

¹ Does not contact sample ² For Inorganic Anions, and physical parameters with long holding times only

³ Cooled by a refrigeration unit ⁴ Autosampler bottles usually pre-preserved, unit not cooled

6.27 Low Level Mercury Sampling Equipment

The following equipment is used to collect low level mercury samples.

Table 6.8 Low Level Mercury Sampling Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>
gloves (long & short)	plastic	sample collection
gloves (short)	latex	sample collection
bags (small & large)	plastic	wrapping bottles/equipment
bottles	FEP	sample collection
Masterflex pump	Not Applicable	sample collection
pump tubing	Teflon	sample collection
pump tubing	C-FLEX	sample collection
battery (12 volt)	Not Applicable	pump operation
filter units, high capacity FS8200 from QED (0.45 micron)	nylon	filtrations
screen (100 micron)	nitex	pre-screening
filter holder	Teflon	hold screen
sampling boom	PVC pipe with LDPE sample bottle attached & plastic ties	hold tubing at desired depth while collecting sample, when using the pump only.

6.2.8 Macrobenthic Invertebrate Sampling Equipment

The following equipment is used to collect macrobenthic samples. The equipment used will depend on the habitats available in the area.

Table 6.9 Macrobenthic Invertebrate Sampling Equipment

<u>Equipment Description</u>	<u>Type of Material</u>	<u>Use</u>
Hester-Dendy artificial substrates	Tempered hardboard	Sample collection
dip nets	Nylon	Sample collection
Petite Ponar	Steel	Sample collection
Ekman grab sampler	Steel	Sample collection
Sieve (U.S. Standard No. 30)	Steel	Processing
10% buffered formalin		Preservation
70-80% ethyl alcohol		Preservation
Glass jars	Glass	Sample Storage
Sorting trays	Stainless Steel	Processing

6.2.9 Miscellaneous Equipment

Sample bottles

60 ml Nalgene polyethylene sample bottles - for inorganics only

175 ml Nalgene polyethylene sample bottles - for inorganics only

250 mL Nalgene polyethylene sample bottles - for inorganics only

500 mL Nalgene polyethylene sample bottles - for inorganics only

1000 mL amber Nalgene polyethylene sample bottles - for inorganics and chlorophyll

125 mL amber Nalgene polyethylene rectangular sample bottles (mercury only)

250 mL Nalgene polyethylene sample bottles (trace metals only)

Note: Sample containers for VOCs organics & extractable organics are provided by contract laboratories.

Microbiology sample containers
125 ml Presterilized Whirlpaks™

Sample preservation supplies

1. 50% H₂SO₄ in 60 mL dropping bottle (fresh weekly, ACS reagent grade)
2. 50% HNO₃ in 60 mL dropping bottle (fresh weekly, Trace Metal grade)
3. pH strips, 0 - 3 range
4. Safety goggles
5. Neutralizing buffer solution, pH 6.9, for acid spill clean up

Note: Preservatives for VOCs, organics and extractable organics are provided by contract laboratories.

Hydrolab field calibration kit

1. Certified pH buffers - 4, 7, and 10
2. Certified conductivity standards - range appropriate for sampling sites for particular trip
3. Ring stand
4. Screwdrivers - straight edge and Phillips
5. Calibration cup with both hard and soft end caps
6. DO membranes and O-rings
7. DO and pH electrolyte
8. Turbidity Standards
9. pH cell caps
10. Technical manual and tools

Coolers with wet ice

1. Igloo brand or equivalent (48, 84 or 196 quart)
2. Gott brand or equivalent (60 or 80 quart)

QA/QC supplies

1. Trip Spike/Trip Blanks
2. Analyte free water - for field blanks and rinsing equipment

Miscellaneous supplies

1. Five (5) liter polyethylene bucket - for inorganic sample processing
2. Disposable latex/PVC gloves, powder free (PVC not for organics or VOCs)
3. Polyethylene trays - for providing clean working areas (does not contact sample)
4. Sample tags
5. #16 rubber bands
6. Waterproof pens (Sharpies)
7. Clipboard
8. Field notebook (bound, waterproof)
9. Chemistry field data log sheets
10. First aid kit
11. Watch (with seconds hand or timer)

12. Personal protective equipment

6.3 Decontamination Procedures

All sampling equipment is transported to the field pre-cleaned and ready to use. Cleaning procedures are described in Section 6.3.1 through 6.3.10 and are as required by EPA SOP & QAM, Appendix B. All sample collection equipment and unpreserved containers are rinsed three times with sample water before the sample is collected with the following exceptions: VOCs, bacteriological samples (Whirlpaks™), Total Recoverable Petroleum Hydrocarbons (TRPHs) and Oil and Grease.

Analyte free water is obtained from the laboratory or field prep areas which have water purification systems. Analyte free water is generated to provide a source of water in which all interferences and analytes are below detection limits. Field personnel use analyte free water to prepare field blanks, equipment blanks, and for the final decontamination rinse of field equipment. The reliability and purity of the analyte free water is monitored through routine conductivity measurements and the results obtained from the equipment and field blanks. The District does not provide a decontamination service to its clients.

6.3.1 Sampling Equipment Cleaning - Laboratory

In the laboratory the sampling equipment is cleaned using the following procedure:

1. Wash all surfaces thoroughly with a hot tap water and phosphate free soap (such as Liquinox™) solution. Use a brush to contact all surfaces and remove stubborn debris. Heavily contaminated equipment is disposed of properly.
2. Rinse thoroughly with hot tap water.
3. Rinse with 10% hydrochloric acid. Use 10% nitric acid for trace metal equipment only.
4. Rinse thoroughly with analyte free water.
5. Rinse thoroughly with pesticide grade isopropyl alcohol (equipment used for organic sampling only).
6. Rinse thoroughly with analyte free water.
7. Allow to air dry completely.
8. Wrap in aluminum foil with the shiny side out (if applicable) or untreated butcher paper, for storage and transportation.

Equipment is properly disposed of if decontamination is not effective. The cleaning procedures used for the field equipment are documented in the equipment logbook and include which equipment was cleaned, the procedure used, the date and initials of the person performing the procedure.

6.3.2 Surface Water Equipment (Inorganics Only)

Sampling equipment for nutrients, major ions and physical parameters that is reused from site to site is rinsed twice with analyte free water at the site and then rinsed three times with sample at the next sample site. Sampling equipment for both trace metals and nutrients is rinsed with 10 % hydrochloric acid, and DI rinsed. When collecting metals only, 20% nitric acid is used instead of the hydrochloric acid. Before the next sampling event, all equipment must be cleaned as stated in Section 6.3.1. When using peristaltic pumps, the exterior must be cleaned with a detergent wash followed by tap and DI water rinses. Hydrolab systems are to be rinsed in a light soap solution and stored in tap water.

6.3.3 Surface Water Equipment (Organics)

A subsurface grab sampler is used to collect organics in surface water. Samples are collected by immersing the sampler and bottle upstream into the surface water body until it is full. The subsurface sampler does not come into contact with the sample.

6.3.4 Ground Water Equipment

All sampling equipment is transported to the field pre-cleaned and ready to use. Cleaning procedures are described in sections 6.3.1 and 6.3.4. The equipment is used only once in the field. The equipment ID numbers are recorded in the field notes at each well. All sample collection equipment and unpreserved sample containers are rinsed once with sample water before the sample is collected with the following exceptions: VOCs, bacteriological samples, Oil & Grease, Total Recoverable Petroleum Hydrocarbons, Bacteriological samples, or any sample bottle containing premeasured preservative. After use, the equipment is rinsed with analyte free water and returned to the sample preparation area for thorough in-house cleaning.

6.3.4.1 Teflon Drop Pipe Decontamination

Teflon drop pipes are decontaminated in-house according to section 6.3.1.

6.3.4.2 PVC Tubing/Hose

A PVC hose connects the purge pump to the Teflon drop pipe when purging groundwater from a well. This hose is decontaminated by rinsing with dilute Liquinox[®] followed by a DI water rinse and wiping dry with clean Kim Wipe[®] towels before purging the next well. This hose does not come in contact with the sample. The drop pipe contains a check valve and is removed while the purge pump is still running to reduce the possibility of water draining back into the well from the inside of the suction hose.

6.3.4.3 Filtration Units

Filtration units (QED FF 8200) are disposable and individually packed in plastic bags to prevent contamination prior to use.

6.3.4.4 Teflon Lanyards

Teflon coated, stainless steel lanyards are decontaminated in-house according to section 6.3.1. In the field, the lanyards are rinsed with dilute Liquinox^R and DI water before using on the next well. All equipment is rinsed three times with sample water before sample collection.

6.3.4.5 Submersible Pumps

Submersible pumps are used for purging only and are decontaminated in-house according to section 6.3.1, with the exception of the acid rinse (due to the stainless steel construction). The exterior of the pump is cleaned between sites per section 6.3.1 (the solvent rinse is optional if organics are not sampled). The interior is cleaned between sites by flushing thoroughly with DI water.

6.3.5 Sediment/Soils Equipment

The field sampling equipment is cleaned prior to being taken to the field and in-between sites by the following procedure: The dredge, scoop, bowl, and core (stainless steel) are washed with LiquinoxTM detergent, rinsed three times with tap water, rinsed with analyte free water, rinsed twice with pesticide grade isopropyl alcohol, and rinsed three times with analyte free water. The isopropyl alcohol is containerized for proper disposal. PVC cores are washed with LiquinoxTM detergent, rinsed with tap water, dipped in 10% HCl, rinsed with tap water, then rinsed with DI water and allowed to dry. The cores are bagged in polypropylene bags until taken into the field.

6.3.6 Sample Bottles and Filtering Equipment

In the laboratory, the reusable sample bottles (except for trace metals) and filter holders are washed with LiquinoxTM detergent, rinsed with hot tap water then with 10% hydrochloric acid, tap water, analyte free water, and finally air dried. Dry bottles are stored with the caps on and and filter holders are stored in a closed container. Unpreserved sample bottles are rinsed with sample three times before the final sample is collected. Disposable bottles, certified precleaned from the vendor are also used. The decontamination procedure of both the reusable and disposable bottles is evaluated by checking 5% of each day's washing or purchased lot number.

Polyethylene trace metal bottles are washed with LiquinoxTM detergent, rinsed with hot tap water, soaked in 20% nitric acid overnight, rinsed with analyte free water, and stored with deionized water containing 1% nitric acid until sample collection. At the time of sampling, the acid is dumped and the bottle rinsed three times with sample water. The decontamination procedures are evaluated by checking 5% of the sample bottles. This is accomplished by filling the selected bottles with DI water and analyzing.

The reusable filter holders are soaked in a weak Liquinox^R solution before washing to soften any residues which may be in the filter holders from the previous sampling process. Reusable filter

holders are then washed following the same procedure as inorganic sample bottles. Disposable filter holders are purchased from Poretics and are not cleaned before use. At least three 10 ml portions of sample water are passed through the filter holder and filters and discarded before filling sample containers.

Filtering syringes are washed with Liquinox™, then rinsed with hot tap water, followed by analyte free water and allowed to air dry.

Clean bottles/equipment are collected from the lab by field personnel and are not shipped. Bottles for the collection of samples sent to contract laboratories are provided by the laboratory performing the analyses. These bottles have been cleaned by the contract laboratory according to that laboratory's procedures prior to shipment to the District. The laboratory must have approved cleaning procedures in their Comprehensive Quality Assurance Plan.

6.3.7 Autosamplers

All autosampler bottles are cleaned in the laboratory using the standard bottle cleaning procedure described in Section 6.3.1 as per EPA SOP & QAM guidelines, Appendix B, Sections B.5 and B.6. Tubing is not cleaned, but rather, replaced. The sampler and sampler tubing is dedicated to the site. The sampler is programmed to rinse the tubing twice with sample before collection and then purge the tubing after the sample is taken. For the American Sigma Model 6201 autosamplers located at pump stations, the sampler and pump tubing are replaced with clean tubing every four storm events, or, at a minimum, quarterly in the absence of storm events. The pump tubing for American Sigma Models 700 and 800 autosamplers is replaced at a minimum quarterly, or preferably, monthly. The entire sampler (American Sigma Model 6201 700, or 800) is cleaned when it is removed from the site for repair or transferred to another site. In an emergency, pump tubing may be cleaned by soaking in hot water and Liquinox® and scrubbing outside and ends with a small bottle brush. The tubing is then rinsed with hot tap water then DI water. If only metals are required the tubing will be rinsed with 10% Nitric Acid. Finally, the tubing is rinsed with DI water and both ends capped with aluminum foil until ready for installation.

6.3.8 Atmospheric Deposition Equipment

Atmospheric Deposition dry and wetfall buckets are cleaned according to section 6.3.1. Once dry, the buckets are capped and taken to the site. Each lot of cleaned buckets is tested for proper decontamination by pouring one liter of DI water into the bucket, swirling, allowing to steep overnight, and analyzing for the parameters of interest.

6.3.9 Low Level Mercury Equipment

All low level mercury sample collection equipment (including sample bottles) is precleaned before use by the contract laboratory using approved methods. The filtration units are disposable. The

SFWMD does not clean any low level mercury equipment. Sample bottles are rinsed three times in the field with sample before the sample is collected. When using the pump, a minimum of three sample hose volumes (200 ml) of sample water are flushed through the tubing or filter before the sample is collected. The contract laboratories ship all equipment clean and in bags. All equipment is taken into the field in plastic bags.

6.3.10 Analyte Free Water Containers

The analyte free water containers are cleaned in the laboratory using the standard bottle cleaning procedures described in Section 6.3.1 as per EPA SOP & QAM guidelines, Appendix B.

6.4 Sampling Procedures


Water sampling locations in the SFWMD are frequently established at water control structures where known flow rates can be combined with chemical concentrations to determine loadings. However, in lakes, rivers, marshes, estuarine systems, storm water runoff, and agricultural point sources, sampling sites are chosen based on other criteria as described in the appropriate Quality Assurance Project Plans. In any case, the sample must be collected in such a manner as to ensure that it is representative of the water body being studied. A new pair of disposable latex/PVC gloves are used at each sampling point for all types of sampling.

6.4.1 Surface Water (except organics and autosamplers)

The surface water sampling procedures (except organics and autosamplers) are given in the flow chart in Figure 6-1. The following special considerations are observed when applicable:

1. If a boat is used, the sample is taken from the bow, and/or upwind and upstream from the motor.
2. When wading, the sample is collected upstream from the boat or body.
3. Care is taken not to disturb the sediment in the immediate sampling area.
4. Pre-preserved containers are not used as collection containers.
5. Intermediate containers are inverted, immersed to the appropriate depth, and turned upright pointed in the direction of flow, if applicable.
6. Preserve samples according to Table 6.10:

Figure 6.1 - Inorganic Surface Water Sampling Procedure



*Preservation/Filtration occurs within 15 minutes of collection

6.4.2 Surface Water - Organics

Surface water samples for organic analysis are collected using a subsurface grab sampler equipped with a precleaned glass bottle provided by the contract laboratory. All unpreserved containers (except for VOCs) are rinsed three times before sample collection. The sample containers are immersed inverted to 0.5 meter below the surface, turned upright, pointing in the direction of flow, until full.

For Volatile organic samples, the water sample is poured slowly down the edge of the 40 mL Teflon-lined septum glass vial from the precleaned glass bottle to minimize aeration. The vial is filled to the point of creating a convex meniscus. The septum is placed, teflon side down, on the meniscus and sealed with the screw cap. The vial is inverted and lightly tapped on the lid to dislodge any entrapped air bubbles. The absence of air bubbles indicates a proper seal. If air bubbles are present, the bottle is opened, additional sample is added and the vial is resealed. Additional sample is added a maximum of three times. If a seal cannot be obtained, the vial, sample, and septum are discarded, a new vial is used and the sampling procedure is repeated.

Each sample is identified by the project code and sequential sample number. During sample collection, the date, time, location, water movement, weather and site conditions are recorded. The preservation technique and holding times are shown in Table 6.10. Following collection of the sample, the bottles are sealed, tagged, and placed in wet ice.

The samples are sent to the certified contract laboratory having a DEP approved Comp QAP with appropriate sample identification and chain of custody form. The samples are packed in coolers with bubble-wrap or other appropriate packing material to avoid breakage. Samples are kept at 4°C with wet ice. Samples for VOC analysis are placed into separate bubble-pack bags for each station. At least one trip blank must be included in each separate cooler. Coolers are taped shut with shipping tape, labeled appropriately and shipped to the laboratory using common carrier overnight delivery.

6.4.3 Surface Water - Autosamplers

Water quality autosamplers are used for the collection of daily composite or discrete samples at selected sample sites. The sampling procedure is controlled by the automatic sampler programming. The automatic sampler, American Sigma Model Streamline 700 or Model 800, is programmed to rinse the sample collection tubing twice prior to sample collection and to purge the tubing following collection. The samplers are programmed depending on the project requirements. For example, the programming may include collection of daily composites at a rate of 80 mL of sample at 144 minute intervals and addition of the sample to the correct composite sample bottle for a total of 10 samples per bottle. When samples are analyzed for total Kjeldahl nitrogen, total phosphorus, ammonia and nitrate plus nitrite only, one ml of 50% sulfuric acid is pre-added to the discrete autosampler bottles before sample collection. This amount of acid has been determined to maintain pH <2 after sample collection.

After sample collection the bottles are capped, thoroughly mixed by inverting five times, labeled, placed on ice immediately and transported to the field laboratory and stored at 4°C. Within 24 hours, aliquots of the samples are transferred to 175 mL bottles, tagged, and stored at 4°C until transported to the analytical laboratory. Samples are submitted to the analytical laboratory within two days of the day the sampler was serviced. If possible, samples are poured into sample bottles at the collection site, preservation is checked and samples are immediately placed on wet ice for transport to the analytical laboratory. These procedures are only followed for long holding time parameters.

Short holding time parameters are not routinely collected using an automatic sampler. However, if short holding time parameters are collected, autosampler bottles are not preacidified. The samples are collected within 24 hours after sample uptake or compositing by the autosampler has been completed, and analyzed within holding times.

Sampler intakes at water control structures are located 1 to 2 feet below the surface and 1 foot off the wing wall.

American Sigma Model 6201 autosamplers located at pump stations are dedicated to the site and work in synchrony with each of the pumps in the pumping battery. The autosampler is activated once any of the pumps becomes operational, and is usually programmed to collect a 50 ml aliquot for every 10,000 RPM on each operating pump. The aliquot is dispensed into a refrigerated five gallon polyethylene bottle. The refrigerated autosampler is deployed for up to one week and samples picked up within 24 hours of after sample collection has been completed. Deployment times may vary depending on project requirements as specified in the project specific Quality Assurance Plan. The refrigerated composite sample is homogenized by shaking the capped bottle. An aliquot is transferred unfiltered to a 250 mL sample bottle for TKN and TPO₄, and another aliquot is filtered into a 60 mL sample bottle, for NO_x and NH₄ testing. Both aliquots are preserved as indicated in Table 6.10.

Autosamplers are dedicated to a specific location as long as they are working properly. If an autosampler fails to perform, the sampler head including the pump hose, electronics, and flow sensor are replaced as required. Sampler intake tubing is dedicated to a collection site.

6.4.4 Sediment/Soil

Sediment samples for organics are collected by hand grab using a stainless steel scoop/spoon, stainless steel petite Ponar™ dredge, or stainless steel core. The stainless steel scoop/spoon is utilized only in quiescent shallow waters. The petite Ponar™ dredge is effective over the wide range of circumstances encountered during the collection of sediment samples. The petite Ponar™ dredge can be used on soft or hard bottoms under a variety of flow and depth conditions. Sediment samples from canals are collected mid-stream using a small boat. The dredge is lowered by rope until contact with the bottom is established. Then the dredge is slowly retrieved and emptied into a stainless steel bowl/tray. Samples are not composited or sectioned into subsamples. The stainless steel core is used when collecting shallow sediment

samples. Each sample (except VOC) is thoroughly mixed in the stainless steel bowl with a precleaned stainless steel spoon. The sediment in the bowl should be scraped from the walls and bottom, rolled to the middle and initially mixed. The mass is quartered, moved to the opposite sides of the bowl, and each quarter individually mixed. Each quarter is then rolled to the center of the bowl and the entire mass is mixed again. This process is continued until the mass is as homogeneous as possible. A precleaned stainless steel spoon or small scoop is used to transfer the required amount of final homogenized material from near the center of the mass into the appropriate sample bottle. This is done quickly to prevent oxidation of metal ions or volatilization of organic compounds. The bottle is filled to its rim.

Duplicates are collected to measure the variability inherent in the sampling process. Duplicates for sediments are collected from the same sampling device. Since a true split sediment sample is almost impossible under field conditions, split soil samples are considered duplicates. If analyses for pesticides, extractable organic compounds, or VOC are to be performed, the sample containers must be glass with teflon lined lids. VOCs bottles must have a teflon-lined septum. For other classes of analytes, glass or plastic jars may be used as sample containers. The containers are completely filled with sample so that air bubbles are not trapped in the container. This is necessary in order to minimize sample oxidation that could influence certain test results. The best way to minimize trapped air is to very gently pack the sample into the container with the spoon or scoop. Preserve according to Table 6.12.

Sediment samples are collected using PVC cores if the sample will not be analyzed for organics. Intact sediment cores are obtained by driving a PVC coring tube to a depth of approximately 50 cm into the soil. Under shallow conditions (water depth < 30 cm) this is achieved using a block of wood to protect the neck of the core tube and then striking the block several times with a hammer such that the core penetrates the soil with minimal compaction. If the water is greater than 30 cm, a PVC coupling should be attached to the coring device. This attachment is comprised of a ball valve and a length of PVC pipe. This is fitted to the neck of the PVC coring tube and a closed ended piece of metal pipe with a diameter larger than the upright PVC pipe is used to pound the core into the sediment. During the insertion, the ball valve should be open.

6.4.5 Ground Water

Ground water samples are collected from the least to most contaminated wells. The order of collection is: 1) VOC, 2) Extractable Organics, 3) Total Metals, 4) Dissolved Metals, 5) Inorganics and 6) Radionuclides. The SFWMD does not use temporary well points or dedicated equipment for ground water monitoring. Samples are not collected from wells with free product showing.

A protective covering of visqueen plastic is placed on the ground at all sampling sites to reduce the potential for contamination. A fresh pair of disposable non-powdered latex gloves are worn at each sample site while purging and another fresh pair worn for collecting samples.

All samples are collected using a Teflon bailer which is suspended by a Teflon coated stainless steel lanyard and carefully lowered into the well. Disposable lanyards are not used. Prior to the collection of the sample, the bailer is rinsed three times with sample. The lanyard and bailer are not allowed to touch the ground during sampling or purging.

The bailer, filled with sample, is retrieved from the well and placed on a bailer stand. VOC samples are collected first via a controlled-flow bailer bottom and are poured slowly down the side of the sample vial to minimize sample aeration. The vial is filled to the point of creating a convex meniscus. The cap is secured with the Teflon side of the septum contacting the sample and no headspace in the sample container. The vial is inverted and gently tapped to locate bubbles. If bubbles are present, the sample in the vial is discarded and sample is recollected and checked for bubbles. Additional sample is added a maximum of three times. If unsuccessful, the sample, vial and septum are discarded, a new vial and septum are used, and the collection procedure begins again.

To filter a sample, a one piece, molded in-line disposable filter is attached to the bottom of the bailer with a Teflon adaptor. The top of the bailer is attached to a vacuum pump by a Teflon attachment so the bailer can be pressurized to force water through the filter. The filters are 0.45 micron unless a dissolved metals sample is collected for permit compliance, in which case the filter pore size is required to be 1.0 micron. The first 100 mls of sample to pass through the filter is discarded as rinse water, and the last 100 mls of sample water in the bailer is not used since it has been in contact with the air at the top of the bailer.

Bottles are supplied by the SFWMD laboratory. If the analyses are to be performed by a contract laboratory, the bottles are cleaned by the contract lab according to that laboratory's procedures prior to shipment to SFWMD. The contract laboratory must have approved cleaning procedures in their Comprehensive Quality Assurance Plan, and perform a 5% analytical check on the washing procedure. Certified disposable VOC vials are purchased from I-Chem.

Duplicates are collected by sampling from consecutive bailers.

Splits are not routinely collected in the groundwater program. When splits are required to measure the performance between two or more laboratories, the following procedure is used:

Splits are collected from the same bailer. For large volume samples that may require more than one bailer load, the first half volume of the first bailer load is poured into the first set of containers and the second half in the second set of containers. Then the first half-volume of the second bailer load is poured into the second set of containers and the second half in the first set of containers, etc., until both sets are full.

The sample bottles are labeled with the date, sample number and project name before sampling at a site begins. After collection the samples are preserved, and immediately placed on ice in a closed container. Exposure of the organic samples to sunlight is kept to a minimum. When VOCs are collected, trip blanks are kept with the organic sample bottles to make certain that the samples have not become contaminated. These trip blanks accompany the samples from the time the empty sample bottles are shipped from the contract laboratory until the samples are analyzed. Trip blanks account for at least 5% of the samples that are analyzed. At least one trip blank must be included in each separate cooler.

6.4.5.1 Wells with In Place Plumbing

These wells are purged for a minimum of 15 minutes, until three bore volumes have been removed, or until the well has chemically stabilized, whichever is greater. The sample is taken from the faucet closest to the source and before any screens, aerators, or filters. The flow rate is reduced to less than 500 ml/min when collecting samples to avoid any undue disturbance. Unfiltered samples are collected directly into the sample containers from the spigot. Filtered samples are collected by filling a bailer from the spigot. The bailer is rinsed three times with sample, then filled and pumped through a filter attached to the bottom of the bailer with a Teflon adaptor. Sample bottles are rinsed once with sample (unless sample bottles are prepreserved, collection is for Oil & Grease, bacteriologicals, Total Recoverable Petroleum Hydrocarbons or VOCs).

6.4.5.2 Flowing Artesian Wells

These wells are purged until three bore volumes have been evacuated, or until the well has chemically stabilized, whichever is greater. The flow of water from the well is adjusted to minimize the aeration and disturbance of samples. Unfiltered and filtered samples are collected directly from the discharging water into a bailer.

6.4.5.3 Monitoring Wells

The depth to water in each well, relative to a measuring point is measured using a graduated stainless steel tape and chalk or an electronic water level indicator. The depth to water is measured twice. Both values, which must be accurate to within 1/10th of a foot, are recorded in the sampler's field log. More stringent measurements may be required for specific project programs. These values will later be used to calculate water elevation relative to mean sea level.

The same measuring point is used every time the well is sampled. When the wells are surveyed, the survey point is used for the measuring point. Water level measuring devices are rinsed with dilute Liquinox[®] and DI water and wiped dry with a Kim Wipe[®] towel before measuring each well.

Prior to the collection of samples, each well is then purged by a purge pump, submersible pump or centrifugal pump. Purge pumps are gas powered and extreme care must be used when handling and locating these pumps to minimize on-site contamination (placed away and downwind from any sampling activities). The purge pump consists of a centrifugal suction pump that does not contact the water within the well. A flexible PVC suction hose connects the pump to a three foot length of rigid Teflon drop pipe. The portion of the drop pipe that comes in contact with the water column is cleaned prior to use as per the procedures outlined for other sampling equipment in section 6.3. The drop pipe is equipped with a Teflon check valve at its bottom end to prevent the back flow of purged water into the well. The drop pipe is lowered to the top of the water column so that the purging process removes all of the standing water. If the water level is drawn down during the purging process and air begins to enter the drop pipe, the drop pipe is lowered to the new water level. If the water level is drawn down too severely, the purge rate is decreased. The SFWMD does not collect samples at any well with a depth to water of more than 25 feet. Submersible pumps are used for purging only. They are decontaminated as specified in section 6.3.4.5.

The well is purged until a minimum of three standing water volumes are removed. The volume of water to be purged is calculated using the following formula:

$$\text{Minimum Purge Volume in gallons} = (D^2) * (DW - DTW) * 0.1224^{(1)}$$

D = Casing Diameter in inches

DTW = Depth to Water in feet

DW = Depth of Well in feet

(1) The constant 0.1224 is an units conversion factor

The volume of water removed from the well must be calculated to avoid excess purging. The flow rate is estimated by measuring the amount of time required to fill a bucket of known volume. The required purge volume is then divided by the estimated flow rate to find the estimated time necessary to purge the well. Flow is measured several times during purging to be certain that it does not change. Though a minimum of three water column volumes must be purged, the well is not considered purged until the well has chemically stabilized. Temperature, pH and conductivity are monitored and readings are recorded at 5 minute time intervals, or at shorter intervals if necessary. Chemical stability readings of the well are taken at intervals required to purge one half the bore volume. These measurements are made in a flow through chamber to minimize atmospheric contact with the sample. The well is considered to be chemically stable when two consecutive readings of temperature, pH and conductivity are within 5% (at least five minutes apart). The purge volume is noted in the field log. All measurements that were recorded while waiting for the well to stabilize are recorded in the field log. Only the final readings taken after the well has stabilized are input into the sample

results database. The drop pipe is slowly raised out of the well while the purge pump is still on to make certain that all of the water above the drop pipe inlet is purged. This procedure also minimizes the possibility of the back flow of water from the drop pipe or suction hose.

Wells with low hydraulic conductivity sediments that can be purged dry at one liter per minute are not sampled because the water quality is not deemed to be representative of the aquifer water quality. The SFWMD does not sample wells, dry purged or slow recovery, for which recovery cannot be accomplished in four hours or less.

All samples are collected using a Teflon bailer. When handling bailers or sample containers, disposable latex gloves are worn. Prior to the collection of the samples, a Teflon bailer is lowered into the well via a lanyard. The bailer is filled and evacuated three times. This procedure is done as a final rinse of the pre-cleaned Teflon bailer and associated Teflon parts before actual sampling begins. The bailer is suspended by a Teflon coated stainless steel lanyard that is carefully lowered into the well. The lanyard and bailer are not allowed to touch the ground during sampling or purging.

The bailer, filled with sample, is retrieved from the well and placed on a bailer stand. VOC samples are collected first via a controlled-flow bailer bottom to avoid sample aeration. The first 100 mls of sample in the bailer are discarded. The sample bottles are filled from the bailer after discarding the first rinse, according to section 6.3.4.

6.4.5.4 Porewater Wells and Peepers

Porewater is collected using porewater wells and peepers (Appendix G) for research projects only. Once collected, samples are processed and handled as routine groundwater samples.

6.4.6 Biological Tissue

6.4.6.1 Above Ground Macrophyte Biomass

Above ground vegetation is clipped at the sediment-water interface with a pair of sheers or knife, and separated into subsamples by species. Each biomass subsample is placed into an individual large opaque plastic bag, labeled, and transported to the laboratory for processing.

6.4.6.2 Below Ground Macrophyte Biomass

After the above ground vegetation has been removed from the sampling quadrant, a soil core is collected from the center of the quadrant using a PVC soil corer. Each soil core is placed into an individual large opaque plastic bag, labeled, and transported to the laboratory for processing.

6.4.7 Hazardous Wastes/Drums

The SFWMD does not sample hazardous waste or drums.

6.4.8 Waste Water

The SFWMD does not sample waste water.

6.4.9 Microbiological

Microbiological samples are collected with extra care in order to prevent sample contamination. The samples are not composited. The personnel collecting the samples do not touch the rims/top of the Whirlpaks™. A Whirlpak™ must be discarded if it is suspected that the tip portion has been touched. Microbiological sampling is closely coordinated with the laboratory performing the analyses to ensure that samples are analyzed within holding times.

6.4.10 Oil and Grease

Oil and Grease samples are collected directly into the sample container without rinsing with sample. Sample containers are not pre-acidified. Sample bottles for Oil and Grease are provided by the contract laboratory.

6.4.11 Low Level Mercury

The SFWMD follows the clean hands (CH), dirty hands (DH) technique as specified in *EPA Method 1669: Sampling Ambient Water for Trace Metals at EPA Water Quality Criteria Levels* when collecting low level total and methyl mercury samples in surface water. CH, using shoulder length plastic gloves, followed by wrist length plastic or latex gloves, touches only clean surfaces (such as new gloves, new plastic bags, or surfaces that have been cleaned in an ultra-clean facility). DH, using wrist length latex or plastic gloves, touches the cooler and the sample bottle covered with two plastic bags. DH opens the outer plastic bag, CH reaches in, pulls up bottle covered with inner plastic bag, and removes the bottle from the inner bag. CH places the inner bag into the outer bag. DH seals the outer bag.

CH reaches into the water or slowly wades (so as not to suspend sediment) and collects 100 ml of surface water by immersing the sample bottle about 10 cm under the surface. The bottle is rinsed by shaking the contents (cap on) and dumping the contents away and downstream from the sampling area. This rinse is repeated two more times. The bottle is filled a final time with more than 100 mL and the cap secured tightly by hand. CH then wades out of the water, puts the bottle in the inside bag, seals it, and places the inside bag into the outside bag. DH seals the outside bag and places it back inside the container. DH records the sample bottle number.

time, site and other information on the field data log sheets, and later transfers the information to the field notebook. Relevant site conditions are recorded into the field notebook.

Samples may also be collected using a pump. The CH/DH procedure is used. DH removes the sampling train (filter holder with nitex screen, Teflon tubing and C-FLEX tubing--all connected) from the cooler and opens the outside bag. CH opens the inside bag and removes the sampling train. While CH holds the tubing, DH secures it to the sampling boom with tie wraps. DH opens the pump head, CH places the pump tubing inside, and DH closes the pump head. DH opens the ring stand clamp, CH places the pump tubing inside and DH closes the clamp.

To begin sample collection, DH positions the end of the sampling train about 10 cm below the water surface and about two meters from shore and secures it to the sampling platform (i.e. place boom end under battery handle). DH starts the pump and begins flushing. DH removes a sample bottle from the cooler and opens the outside bag. CH opens the inside bag and removes the sample bottle. CH places the inner bag into the outer bag. DH seals the outer bag. CH opens the sample bottle. The contents are dumped into a waste container if they are dilute HCL or on the ground (away from the sampling area) if contents are deionized water. After flushing with three sampling-hose volumes of sample water (approximately 400 mL), CH fills the bottle with about 50 mL from the sample water stream. Note: keep cap faced downwards. The bottle is rinsed by shaking the contents (cap on) and dumping the contents away from the sampling area. This rinse is repeated two more times. The bottle is filled a final time with more than 100 mL and the cap secured tightly by hand. CH then puts the bottle in the inside bag, seals it, and stuffs the inside bag down inside the outside bag. DH seals the outside bag and places it back inside the cooler.

To prepare a filtered surface water sample, DH stops the pump, removes a filter cartridge from a cooler and opens the outside bag. CH opens the inside bag, removes the filter and connects it to the end of the sampling train. DH starts the pump and the filtered sample is collected as above. Equipment is rinsed with sample water at the next site before collecting a sample. To control mercury contamination, CH secures a clean plastic bag to each end of the sampling train for transport between sites. The equipment and supplies are stored inside coolers and bags in a place relatively free of mercury contamination. Dust and human breath (from dental amalgams) contain large amounts of mercury and should be avoided during sample collection. Samples should not be collected in the rain.

6.4.12 Benthic Macroinvertebrates

The SFWMD follows macroinvertebrate sampling methods as detailed in the Florida Department of Environmental Protection (FDEP) Biology Section Standard Operating Procedures (Appendix E) and in EPA/600/0-90/000, *Macroinvertebrate Field and Laboratory Methods for Evaluating the Biological Integrity of Surface Waters*. The field methods used to study benthic macroinvertebrate communities relate to qualitative, semi-quantitative and quantitative sampling, and are dependent on the data quality objectives of the study. Generally, there are four categories of macrobenthic invertebrate samples collected: cores, grabs, artificial

substrates, and miscellaneous. Field equipment is cleaned prior to the field trip and after each sample according to the procedures described in Section 6.3.5.

a. Corers

Coring devices may be used at various depths in any substrate that is sufficiently compacted so that an undisturbed sample is retained; however they are best suited for sampling homogeneous soft sediments, such as silt, clay or sand. Macrofauna are sampled using single and multiple-head corers, tubular inverting samplers, open-ended stovepipe devices, hand corers and/or box corers. Devices are lowered slowly into the substrate to ensure good penetration and to prevent organism escapement. Visual inspection of each sample is necessary to ensure an adequate amount of sample is obtained.

b. Grab samplers

Grab samplers are designed to penetrate the substrate by gravity and have spring- or gravity-activated closing mechanisms. They are used to sample a unit area of the habitat. The habitat and substrate type sampled, depth of penetration, angle and completeness of jaw closure, loss of sample during retrieval, disturbance at the water-sediment interface, and effect of high flow velocities all affect the quantity and species of macroinvertebrates collected by a particular grab. Petite Ponar™ or Ekman grab samplers are typically used; however the type and size of the selected device depends on the substrate composition, water depth, and hoisting gear available. The sampler is lowered slowly to avoid a disturbance of the surface sediment and to ensure the device bites vertically. Upon tripping of the closing mechanism, the sediment is contained and the device is retrieved. The sample is then placed into a suitable container for transporting to the laboratory or is placed directly into a sieving device and processed.

c. Artificial substrate samplers

Artificial substrate samplers, such as the multiplate (modified Hester-Dendy) sampler, consist of natural or artificial materials of various composition and configuration. They are placed in the waterbody for a predetermined period of exposure and depth for the colonization of indigenous macroinvertebrates. This type of sampling is used to augment bottom substrate sampling because many of the physical variables encountered in bottom sampling are minimized, e.g. light attenuation, temperature changes, and substrate variation. The samplers are deployed in the euphotic zone of the water column to promote maximum colonization. They are exposed for a preset period of time and retrieved vertically from the water. The samplers are placed in a preservative container and transported to the laboratory for processing.

d. Miscellaneous qualitative devices

Many devices such as dip nets, bare hands, tongs, and forceps may be used to collect benthic macroinvertebrates. Dip nets are used by sweeping the net through the water or by holding the net stationary against the bottom and disturbing the substrate, causing the benthos to be swept into the net.

Samples collected by any of the aforementioned devices contain varying amounts of fine materials which can be removed by sieving immediately in the field, or in the laboratory. If

laboratory sorting can be performed within 24 hours, place samples immediately on ice. If sorting will be delayed, preserve with 10% formalin. After organisms have been removed from detritus, they should be placed into 70% ethanol.

6.4.13 Atmospheric Deposition

Atmospheric deposition samples are collected according to the National Atmospheric Deposition Program (NADP) guidelines.

The District uses Aerochem Metric Wet/Dry Precipitation collectors. The Aerochem is an aluminum structure with two 3.5 gallon PVC buckets which serve as receptacles for wet and dry deposition. An aluminum lid covers the wet side bucket under dry conditions. When the humidity sensor detects rain, the lid slides over the dry side bucket, uncovering the wet side bucket. All buckets are collected on Tuesdays only, according to NADP guidelines. The buckets are capped (while using PVC/Latex gloves) and taken to the laboratory for processing within 24 hours. Clean buckets are placed into the Aerochem collectors. Before clean buckets are placed into the collectors, one bucket from each decontamination lot is checked by pouring one liter of DI water into the bucket, allowing the water to sit overnight (capped), and processing as an equipment blank.

All wet side buckets are weighed using a top loading balance and subtracting the clean bucket tare weight (determined before deployment in the field). Live frogs, insects and nonrepresentative foreign matter are removed with a clean scoop or tweezers. Any removed contaminants are documented in the field notebook. The sample is then processed as a routine surface water sample, including filtration and preservation as required. Processing equipment blanks, splits and field blanks are submitted according to routine QA procedures.

Dry side buckets are inspected for contamination as above. Only loose contamination is removed (with clean tweezers/spatula). Approximately 600 ml of DI water is added to the dry bucket contents rinsing the sides of the bucket. A precleaned spatula is used to rub the inside bottom and sides of the bucket. The total amount of DI water added should not exceed one liter and the volume must be recorded in the field notebook. The sample is then processed as a routine surface water sample, including filtration and preservation as required. Processing equipment blanks, splits and field blanks are submitted according to routine QA procedures.

6.4.14 Duplicates/Split Samples

Duplicates are collected to measure the variability inherent in the sampling process, and should be obtained by duplicating in rapid succession during the same sampling event, the entire acquisition technique used to obtain the first sample. Duplicates for water are collected by sampling from successively collected volumes. Duplicates for soils are collected from the same sample source (i.e., soil obtained from the same soil sampling device). Split samples are collected to measure the variability between laboratories, and should be obtained as subsamples

taken from the same parent sample (a true split of soil, sediment or sludge is almost impossible to accomplish under field conditions; a split soil sample should be considered a duplicate and not a split sample). Split samples (for water) are taken from consecutive sample volumes from the same sampling device (i.e. from same bailer), or by mixing in a large intermediate vessel.

Note: for large volume samples that may require more than one bailer load, the first half volume of the first bailer load is poured into the first set of containers (the second half into the second set of containers), then first half volume of second bailer load is poured into the second set of containers (second half into first set of containers) etc. until both sets are full.

6.5 Documentation

The following is a list of the field records that are maintained.

1. Field Trip Checklist (surface water/groundwater)
2. Chemistry Field Data Log
3. Bound field notebook (project specific)

6.6 Preservation

Holding times and preservation techniques for each parameter are given in Tables 6-10, 6-11, 6-12, 6-13 and 6-14.

Samples are preserved in the field at the time of sample collection either by using preservatives provided by the SFWMD laboratory or the contract laboratory. If sample containers are received with preservative from the contract laboratory (VOCs), additional preservative is requested so it is available if needed. ACS reagent grade or better preservatives are used.

For unfiltered samples, the appropriate preservative is added to the sample bottle after the bottle is filled. The bottle is capped and shaken after which a small amount of the sample is poured onto the pH (0 - 3 pH units) test strip to ensure $\text{pH} < 2$. If pH is not < 2 , additional acid is added drop-wise, the bottle is capped and shaken, and the pH is tested again. This procedure is followed until $\text{pH} < 2$. The pH is checked on all samples requiring preservative and this amount of acid is added to the equipment blank and field blank. The amount of acid added is recorded in the field notebook.

For filtered samples, the acid is added after filtration following the procedure outlined for unfiltered samples.

Fresh preservatives are obtained from stocks weekly for all sampling trips. Preservatives are taken into the field in polyethylene dropper bottles. One mL nitric acid (70%) glass ampoules are used in groundwater sampling for the preservation of metals samples.

6.7 Sample Dispatch

Samples to be analyzed by the SFWMD laboratory are submitted to the laboratory by field personnel. The majority of samples are submitted the same day they are collected but all samples are submitted as soon as possible after collection in order to meet recommended holding times.

Samples to be analyzed by contract laboratories are shipped to the laboratory by common carrier overnight delivery the same or next day they are collected. All samples are carefully packed with appropriate material to prevent breakage. Insulated coolers are used for sample shipment and are sealed with shipping tape to avoid tampering. If samples must be kept at 4°C, wet ice is used during shipping.

6.8 Field Waste Disposal

All field generated wastes and contaminated purge waters are segregated and containerized for proper disposal including the use of the services of a commercial disposal company, if required. The SFWMD does not sample hazardous waste sites so the only field generated wastes are acids and isopropyl alcohol. The isopropyl alcohol is containerized for proper disposal. Acids are diluted and/or neutralized and flushed down the sanitary sewer. The calibration standards for field parameters are flushed into the sanitary sewer.

Table 6.10 Holding Time and Preservation for Water Inorganics¹

Parameter	Holding Time	Preservation	Container & Size
Alkaline Phosphatase	24 hours	Cool, 4°C	Plastic, 175 ml
Alkalinity	14 days	Cool, 4°C	Plastic, 175 ml
Ammonia	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Plastic, 60 ml
Biochemical Oxygen Demand	48 hours	Cool, 4°C	Plastic, 500 ml
Carbonaceous Biochemical Oxygen Demand	48 hours	Cool, 4°C	Plastic, 500 ml
Chemical Oxygen Demand	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Chloride	28 days	None required	Plastic, 60 ml
Chlorine, Total Residual	Analyze immediately	None required	Plastic/glass, 250 ml
Chromium VI	24 hours	Cool, 4°C	Plastic, 250 ml
Color	48 hours	Cool, 4°C	Plastic, 60 ml
Fluoride	28 days	None required	Plastic, 500 ml
Inorganic Carbon	14 days	Cool, 4°C	Plastic, 175 ml ²
pH	Analyze immediately	None required	Plastic, 175 ml ³
Kjeldahl nitrogen	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Mercury	28 days	pH<2(HNO ₃)	Plastic, 125 ml
Metals	6 months	pH<2(HNO ₃)	Plastic, 250 ml
Nitrate + nitrite	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Plastic, 60 ml
Nitrite	48 hours	Cool, 4°C	Plastic, 60 ml
Oil and Grease	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Glass, 1 Liter
Organic Carbon	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Orthophosphate	48 hours	Filter immediately, cool, 4°C	Plastic, 60 ml
Total phosphorus	28 days	Cool, 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Residue, filterable, nonfilterable, volatile	7 days	Cool, 4°C	Plastic, 500 ml ³
Silica	28 days	Cool, 4°C	Plastic, 60 ml
Specific conductance	28 days	Cool, 4°C	Plastic, 500 ml ³
Sulfate	28 days	Cool, 4°C	Plastic, 60 ml
Turbidity	48 hours	Cool, 4°C	Plastic, 500 ml ³

(1) From 40 CFR, Part 136, Table II (7-1-90)

(2) Not listed in 40 CFR, Part 136, Table II (7-1-90) (3) 250 ml plastic bottle may be used if BOD is not requested

Table 6.11 Holding Time and Preservation for Water Organics

Parameter	Holding Time	Preservation ¹	Container & Size
Volatile (Purgeable) Organics	14 days (preserved with 1:1 HCl), 7 days (unpreserved)	Cool, 4°C	Glass, 40 ml. Teflon lined septum
Base neutral acid extractable compounds	7 days until extraction, 40 days after extraction	Cool, 4°C	Amber Glass, 1 L. Teflon lined cap
Organochlorine Pesticides and PCB's	7 days until extraction, 40 days after extraction	Cool, 4°C	Amber Glass, 1 L. Teflon lined cap

(1) If Residual chlorine is present, Sodium Thiosulfate is added to the sample vial first. The vial is then filled to almost full volume with sample, acid is added, and finally the vial is filled as per procedure.

Note: It is not recommended to mix the two preservatives (and sample) together in an intermediate vessel.

Table 6.12 Holding Times and Preservation for Sediments, Soils and Tissues¹

Parameter	Holding Time	Preservation	Container & Size
Volatile organics	14 days	Cool, 4°C	Glass, 4 oz. widemouth with Teflon/silicone septum
Semivolatile organics	14 days until extraction, 40 days after extraction	Cool, 4°C	Glass, 8 oz. widemouth with Teflon/silicone septum (50 grams)
Total metals	6 months	Cool, 4°C	Glass or plastic, 8 oz. widemouth (200 grams)
Mercury	28 days	Cool, 4°C	Glass or plastic, 8 oz. widemouth (200 grams)
Nutrients & Inorganics ²	Not Specified	Cool, 4°C	Glass, 500 ml or plastic, 8 oz. widemouth with Teflon lined closure

(1) From Table 5, Chapter 62-160, F. A. C.

(2) From USEPA Standard Operating Procedures and Quality Assurance Manual (Appendix A).

Table 6.13 Holding Times and Preservation for Microbiologicals⁴

Parameter	Holding Time	Preservation ²	Container & Size ³
Microbiologicals ¹	6 hours	Cool. 4°C/Na ₂ S ₂ O ₃	Glass or Plastic > 125 ml

(1) Parameters included are: Fecal Coliform, Total Coliform and Fecal Streptococci.

(2) Addition of sodium thiosulfate is only required if the sample has a detectable amount of residual chlorine, as indicated by a field test using EPA Method 330.4 or 330.2 or equivalent.

(3) Presterilized Whirlpak bags (or equivalent) are typically used for sampling.

(4) From Table 8, chapter 62-160, F.A.C.

Table 6.14 Holding Times and Preservation for Biologicals

Parameter	Holding Time	Preservation	Container & Size
Chlorophyll	(i) 24 hours to filter (ii) 21 days after filtration	(i) Cool. 4°C, dark, filter with MgCO ₃ . (ii) Filter frozen (until testing)	Amber Plastic, 1 L
Macrobenthic Invertebrates	Preserved in the field with formalin, then ethanol in the laboratory or placed on ice and preserved with formalin within 8 hours, then preserved with ethanol	10% buffered formalin, then 70% ethanol	Glass or Plastic

6.9 Field Reagent and Standard Storage

Table 6-15 lists the storage procedures for reagents, standards and solvents.

Table 6-15. Reagent, Solvent and Standard Storage

<u>Chemical</u>	<u>Method of Storage</u>
Sulfuric Acid	Stored in original containers in vented acid storage cabinet. Note: Each acid is stored in a separate cabinet. Acid is transported in the field in polyethylene dropper bottles.
Nitric Acid	See above.
Phosphoric Acid	See above.
Hydrochloric Acid	See above.
Isopropyl Alcohol	Stored in original containers in vented solvent storage cabinet. No other solvents are stored in the same cabinet. Taken into the field in glass containers carried in a safety carrier.
pH Standards	Stored in cabinet designated for standard and reagent storage. Cabinet is in air conditioned laboratory. Standards taken into the field in polyethylene bottles.
Conductivity Standards	Stored in cabinet designated for standard and reagent storage. Cabinet is in air conditioned laboratory. Standards taken into the field in polyethylene bottles.
Formalin (10%)	Stored in cabinet designated for standard and reagent storage. Cabinet is in air conditioned laboratory. Taken into the field in polyethylene bottles.
Ethyl Alcohol	Stored in original containers in vented solvent storage cabinet. No other solvents are stored in the same cabinet. Taken into the field in approved non-combustible containers.

7.0 SAMPLE CUSTODY

A verifiable trail of documentation for each sample must be maintained from the time of sample collection through the analytical laboratory to the final reporting or archiving of data. SFWMD does not require legal custody procedures be used as the only samples analyzed in the District laboratory are not used in criminal court cases at this time.

The purpose of sample custody is to provide a clear description of sample and container traceability from sample collection to final sample disposition and to identify those persons responsible for collection and analysis.

7.1 Documentation

Custody starts with the person who prepares for the field trip. This person, who may be a Staff Scientific Associate, Senior Scientific Associate, Specialist Scientific Associate or Senior Supervising Scientific Associate, obtains the appropriate size precleaned bottles from the bottle storage area. This person signs the trip checklist (Figure 7-1) when all items necessary have been prepared and loaded into the vehicle.

The field sheets become a part of the project records maintained by the project manager with copies retained by the laboratory. Entries on all records, laboratory and field, are made in waterproof ink with errors deleted by crossing one line through them. All field documentation is signed or initialed by the field personnel.

Samples of the documentation are given in Figures 7-1 through 7-4.

Figures 7.2 and 7.3 show the sample tags which include the preservation instructions, bottle size, and the parameters to be tested. Information on the tag includes project code, the date the sample was taken and the sample number. The sample number is a unique sequential number for that project. Sample seals are not used by SFWMD.

Figure 7-4 shows the Chemistry Field Data Log used to record the field information. This data is manually entered into the Laboratory Information Management System (LIMS) through the Sample Log-In screens shown in Figure 7-5 and 7-6. This data entry process automatically generates a unique sequential number consisting of the login group hyphenated with the number of the sample for that group. A login group is defined as a group of samples for a given project on a given day. The sample number is used to track the progress of the sample through the laboratory.

The Chemistry Field Data Log includes the site name or station code, the date and time of sample collection, the signature of the person relinquishing the samples, the field ID number, the number of samples collected, the intended analyses and preservation requirements, a comment section, and a place for the person who receives and logs-in the sample to sign, date, and record the corresponding laboratory sample numbers. Common carrier is not used for

Figure 7.1

FIELD TRIP PREPARATION AND CHECK-OFF LIST

PROGRAM: _____ DATE: _____ PROGRAM CODE: _____
(Routine / Quarterly / Bi-Annual)

- | | |
|---|-------------------------------------|
| _____ * Cellular Phone | _____ * Sign Out Vehicle |
| _____ * Sign Out Board | _____ * Gas Vehicle |
| _____ * Maps | _____ * Maintenance Check |
| _____ * Vehicle Packet/Credit Cards | _____ * Coolers/Ice |
| _____ * Pre-Cleaned Buckets | _____ * Syringes |
| _____ * Loaded Filters _____ | _____ * Pre-cleaned Sampling Bottle |
| _____ * Processing Tray | |
| _____ * Calibrated Hydrolab | _____ * Field Data Logs |
| _____ * Field Notebook _____ | _____ * Waterproof Pens |
| _____ * Tags/Rubber Bands | |
| _____ * Fresh Acids- H2SO4 / HNO3 | _____ * Trip Spikes |
| _____ * Goggles, Gloves | _____ * pH Test Strips |
| _____ * 5 Gallons DI Water | |
| _____ * Bottles: 250ml- _____ 175ml- _____ 60ml- _____ 175ml- _____ 250ml/TM- _____ | |
| _____ * Bottles: One Liter Bottle for Chlorophyll | |
| _____ * Stations: Routine / As Specified: _____ | |
-
- _____ * Keys: _____

PERSONAL ITEMS

- | | | | |
|-------------------|--------------------|------------------------|----------------------------|
| _____ * Watch | _____ * Sunglasses | _____ * Drinking Water | _____ * Hat |
| _____ * Sunscreen | _____ * Raingear | _____ * Food | _____ * Mosquito Repellent |

POST TRIP PROCEDURES

- _____ * Sort Samples in Sequence and by Tag Color
- _____ * Place Samples in Refrigerator marked "incoming samples"
- _____ * Sign Chain of Custody Sheet/Get Time Stamp
- _____ * Place Chain of Custody Sheet on the "to be logged in" board
- _____ * Return Clean Bottles to Lab
- _____ * Separate Filter Holders and place in Soak Solution
- _____ * Place Dirty Bottles in Washing area
- _____ * Replace Water in Hydrolab Stand(Tap Water Only!)
- _____ * Clean Out Vehicle
- _____ * Fill Out Trip Ticket/Return Keys
- _____ * Return Field Notebooks to Office
- _____ * Report Equipment(Vehicle,Hydrolab, etc.) Problems to Supervisor

Figure 7.3
 Sample Submission Diagram For Filtered Water




Figure 7.5
 Sample Log-In to LIMS

specifyLogin DefaultInfo Samples Products Report Keymenu exit

 Sample Login Menu

Sample Definition

Working on: L5

Sample	Fld Samp	Received	Coll.Date/Time	Station Cd	QC	U	D	W	Smp
						D	S	E	Typ
L5-1	A100	24-JUL-92	23-JUL-92	1000	PELHIL	1	2	1	24

Use [Next Block] to edit product information

Extended Sample Information:

Hydrolab: Depth 2.5 Temp 23 pH 8.1 Spec Cond 2000
 DO 10 SDD 5 Salinity 20 Redox 2500
 Comments this is for documentation purposes
 Project Collector Account Permit #
 CAMB TJ 2225 WATER CONSERVATION AREA

Use [Previous Block] to edit default information

Enter in permit number

Count: *0

<OSC><DBG>

<Replace>

Figure 7.6
Sample Product List

specify login: default info Samples Products Report Keymenu exit
Sample Login Menu

Sample Definition Working on: L5

Product List For sample L5-1

Matrix	Product Code	C-Class	Price	Cost	L-Level	Location
00	Routine	CAMB				
00	Routine	ALCO3				
00	Routine	CA				
00	Routine	CL				
00	Routine	COLOR				
00	Routine	K				
00	Routine	MG				
00	Routine	NA				
00	Routine	NH4				
00	Routine	NO2				

Use [Previous Block] to return. Any changes will be POSTed.
Show children? Use [Menu] 'Products' to toggle display.

FRN-40400: Transaction complete --- 3 records posted and committed.
Count: 10 v <OSC><DBG><List><Repl>

samples analyzed by the District laboratory.

The sample preparation (digestion) records are maintained electronically in LIMS and with the analytical run. The information required is the laboratory sample number, standards, QC samples, dilution factors, person preparing the samples, and the date of preparation. The analysis logs for each type of instrument are shown in Figures 7-7 to 7-12. The microbiology logs are shown in Figures 7-13 to 7-16.

If samples are collected and sent to a contract laboratory by common carrier for analysis, the custody forms supplied by the contract laboratory are used. Only laboratories that have DEP approved comprehensive quality assurance plans are used as contract laboratories. The precleaned sample containers received from the contract laboratory are delivered to the project manager who is responsible for their secured storage.

7.2 Field Custody Protocols

The samples are tagged at the time of collection using waterproof Tyvek™ tags which have been filled out with waterproof pens. The tags are attached with rubber bands to the neck of the sample bottle.

Each project has a unique three to four character project code. The sample field numbers are generated by a combination of this project code and up to a five digit sequential project number that unequivocally links the collected sample to the time and date of collection, and sampling point. The field sample numbers are associated with the station (site) code on the Chemistry Field Data Log and in the field notebook. The field sample number can be used to trace a sample through the sample tag, the Chemistry Field Data Log, the field notebook, the laboratory information management system, and final data archival.

The field records are maintained and stored by the project manager, who is responsible for maintaining all records of the project for the period of time specified by the Florida standards for record management. The field records which are identified by the project code may include the field notebooks and the Chemistry Field Data Logs.

All physical parameter measurements obtained in the field shall be written on the Chemistry Field Data Log shown in Figure 7-4 at the time of sample collection. These measurements include temperature, pH, specific conductivity, dissolved oxygen, oxidation-reduction potential, secchi disc depth, total column depth, turbidity and salinity. Other information that must be noted on this form includes project code, collector ID, sample number, date and time of collection, station or site ID, parameter analysis, types of QC samples collected (when and where collected), depth sample was collected, requested parameters and chain of custody documentation. Beginning and ending times of any composite sampling is noted. Optional information on this form includes upstream/downstream notation, discharge information, weather, sample type, depth and a comments section.

Figure 7.8

FORM 0217
12/84

P-E 3600 I.D./W.T. FILE

I.D. / W.T. FILE: _____
DATE: _____

SAMPLE WEEK: _____
DISC. #: _____

ANALYST: _____

PAGE: _____

		PARAMETER					
		ABSORBANCE					
		DATA FILE					
	LIMS NUMBER	PROJECT CODE	SAMPLE NUMBER	DILUTIONS/RERUNS			
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							
21							
22							
23							
24							
25							

Figure 7.9



Rapid Flow Analyzer Log

Technician _____ Date _____

Parameter _____ Std. Cal _____

High or Low (Circle One)

CF Method Number _____ Instrument # _____

Tray Protocol File Name	Sample Numbers
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

Comments: _____

CALIBRATION:

CHANNEL 1:	CHANNEL 2:
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

QUALITY CONTROL:

QC1: _____
QC2: _____
C.V.: _____
STAD: _____

GAIN: _____

Figure 7.10



Calibration Log for Physical Parameters

ANALYST: _____ DATE _____
PARAMETER: _____

STANDARDS:

Known Value	Instrument Reading
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

REP NO. QUALITY CONTROL SAMPLES:

() QC1 True Value _____	QC2 True Value _____
() _____	() _____
() _____	() _____
() _____	() _____
() _____	() _____
() _____	() _____
() _____	() _____
() _____	() _____

REPEAT SAMPLE READINGS:

() _____	() _____	() _____	() _____
() _____	() _____	() _____	() _____
() _____	() _____	() _____	() _____

SAMPLES ANALYZED: _____

Figure 7.11



DEPARTMENT OF WATER RESOURCES EVALUATION

Ion Chromatography Log

Analyst: _____ Date: _____

Parameter: _____ Method: _____

Workgroup: _____

Samples Analyzed: _____

Comments: _____

Calibration: _____ CI _____ SO4 _____

Retention Time: _____

Standards:	Peak Area	Peak Area	
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

QC1 _____

QC2 _____

C.V. _____

STAD _____

GAIN _____

Figure 7.12



DEPARTMENT OF WATER RESOURCES EVALUATION

Carbon Analyzer Log

Analyst: _____ Date: _____
Parameter: _____ Method: _____
Workgroup: _____

Samples Analyzed: _____

Comments: _____

Calibration:	TC	IC	TOC/DOC
System Blank:	_____	_____	_____
Standards:			
Blank	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
QC1	_____	_____	_____
QC2	_____	_____	_____
C.V.	_____	_____	_____
STAD	_____	_____	_____
GAIN	_____	_____	_____

FECAL STREPTOCOCCUS

MEMBRANE FILTRATION (SM #910 B)

ANALYST

SAMPLE ID													
TIME & DATE COLLECTION													
TIME & DATE ANALYSIS													
DILUTION													
RAW COUNTS													
IN INCUBATION TEMP&DATE													
OUT INCUBATION TEMP&DATE													
MEDIA MANF& LOT#	Difco	Difco	Difco	Difco	Difco	Difco	Difco	Difco	Difco	Difco	Difco	Difco	Difco
FILTER MANF& LOT#	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN
PETRI DISH MANF& LOT#	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN	CELMAN
BUFFER LOT# & PH													
CONFIRMATION													
CATALASE													
MEDIA LOT# & PH BHI AGAR BHI BROTH BHI+40%BILE													
IN TIME DATE TEMP													
OUT TIME DATE TEMP													
CFU/100ML													

Bubbles indicates a + catalase (indicates presence of non streptococcus) No bubbles indicates a - catalase (streptococcus present)

Relevant field observations are noted in a bound waterproof notebook at the time of sample collection. These include sample number, station or site name, date and time, weather, flow conditions, water color, water smell, water clarity, weed conditions, number of drops of acid added to each bottle, persons other than sampling personnel at the site, type of purging and sampling equipment used with corresponding ID# (if available), field decontamination performed and if applicable, wave height, bottom conditions, algae description, use of boats and/or other fuel powered equipment, calibration information and standards used for calibration.

For monitoring wells, the following information is recorded in addition to any applicable information from the above list: depth to water, calculation for purge volume, determination of volume purged, method of purge, purging rate, date and time the well was purged, and readings taken until the well stabilized.

For wells with in place plumbing or artesian wells, the following information is recorded: plumbing or tap material, flow rate at which the well was purged, time the well was allowed to purge, and the flow rate when the sample was collected.

A monitor well database is kept for all of the wells that are sampled by the SFWMD. This includes the following information: well casing material, well diameter, type of casing, screen diameter, screen type, total depth of the well, casing depth, method of well installation, date of well installation, drillers name, latitude, longitude, measuring point elevation, and land surface elevation. In addition, new wells installed by the SFWMD also have the drilling mud type and name recorded.

For sediments, the depth the sample is taken is recorded. Drilling/boring is not used for sampling.

7.2.1 Sample Transport

Following collection of the sample, the bottles are sealed, tagged, and returned to the contract laboratory along with appropriate sample identification (i.e., sample name and number), and chain of custody form. Samples are packed in coolers with bubble-wrap or other appropriate packing material to avoid breakage. Samples for VOC analysis are placed into a separate bubble-pack bag for each station. Chain of custody forms are enclosed in a plastic bag for protection from water damage. Samples are kept at 4°C with wet ice. Coolers are taped shut using packing tape and taken to the SFWMD shipping area for pickup by common carrier. Coolers are shipped to the laboratory overnight, and the shipping receipts are retained. Chain of custody forms are returned to SFWMD with the analytical results.

7.2.2 Sample Transmittal

The Chemistry Field Data Log and examples of the sample tags are shown in Figures 7.2 through 7.4. The Chemistry Field Data Log includes the field ID number, date and time of sample collection, station (site) code, intended analyses (designated by circling desired parameters in figure 7.4), method of preservation, and limited comments about the sample or sample

container. Sample seals are not used.

7.3 Laboratory Operations

The samples are brought to the laboratory by the field sampling personnel and placed in the designated refrigerator. The person bringing the samples to the laboratory signs and stamps the Chemistry Field Data Log with the time clock. The Chemistry Field Data Logs are given to the person responsible for logging-in the sample. The information from the Chemistry Field Data Log is entered into the LIMS and a laboratory sample number is generated by LIMS consisting of a login group hyphenated with the number of the sample within the group. ex: L1-1 is the first sample of the first login group. The person logging-in the samples records the LIMS numbers on the Chemistry Field Data Log, initials and dates the form at the time the samples are logged (in the "received by" section). The LIMS numbers, date logged, and initials of the logger are also maintained in a bound laboratory notebook. A computer-generated list, verified by the data entry technician, is obtained with the respective field and LIMS numbers for each sample.

Labels for each aliquot are generated by LIMS and are manually attached to the sample field tag. As each sample is labeled, the technician checks the proper match of field and LIMS numbers, and the sample bottle for leakage, cracks, and any other obvious faults. Five percent of the samples (at least one from each batch) are checked for proper preservation. The results of this check are documented in a bound notebook.

7.3.1 Sample Rejection

Samples are rejected if the bottle is leaking or cracked, if there is no Chemistry Field Data Log submitted with the samples, if the Chemistry Field Data Log is incomplete, or if holding times or preservation protocols have been violated.

7.3.2 Sample Security, Accessibility, and Storage

The samples are accessible to the laboratory and field staff during working hours. The doors to the exterior of the building are locked when the last member of the staff leaves for the day. The door to the sample receiving area is opened at 6:00 A.M. The only persons authorized to be in the laboratory are the laboratory staff and the sampling personnel delivering samples. All visitors must be escorted by a member of the division staff.

Samples are stored in refrigerators designated for sample storage only. No VOC's are stored in the laboratory refrigerators and standards are stored in refrigerators designated for standard storage.

Sample digestates are stored in the digestion laboratory until they are analyzed. Analysis takes place within 48 hours for nutrients and within 30 days for trace metals.

7.3.3 Sample Distribution and Tracking

The analysts check the LIMS database daily for samples requiring their assigned analyses and make work groups for these samples. They are required to analyze reworks and the oldest samples on each day's first run. The supervisor receives a daily backlog report listing all samples showing incomplete analyses.

The analysts are responsible for removing the required samples from the refrigerators, analyzing them, and returning them to the refrigerators immediately following analysis. Bound notebooks are used to track digested samples. The person doing the digestion lists the samples, standards and QC samples digested in the batch.

Following completion of the analyses and review of the results, the laboratory supervisor or analyst removes the samples which are ready for disposal from the refrigerators. The status of the sample is changed from NEED to DONE in the disposal department in LIMS. Samples are disposed of by diluting and dumping into the sanitary sewer. The SFWMD does not analyze hazardous waste samples. Digestates containing hazardous materials, (i.e. the mercury in Kjeldahl nitrogen digestates), are treated as hazardous waste and are collected and stored for disposal by a commercial hazardous waste company.

7.3.4 Interlab Custody

The SFWMD has only one laboratory and does not transfer samples to any other laboratory. However, in the event of catastrophic failure, samples will be transported to a contract lab by common carrier accompanied by the log in sheets and the contract lab's chain of custody.

7.4 Electronic Data Records

A LIMS database is used by SFWMD for sample tracking, data storage, and data reduction.

7.4.1 Security System

The LIMS database has several levels of security. The Staff Programmer Analyst responsible for the operation and maintenance has the highest level of security and can access all information and programs in LIMS. The Division Director, Laboratory QA Officer, Supervising Professionals and Senior Chemists can access and modify all sample and results information. The analysts can enter analytical results but cannot modify the sample information or results. The project managers can only read and print results. User names and passwords are assigned to each person. Back up of both the file system and database is performed daily so that in the event of catastrophe, only one day of data would be lost. One week's worth of system files is stored in the vault.

7.4.2 Forms

Project Managers review weekly data results and identify tests for rework or field sample data for modification. These copies are stored by date produced and are kept according to the record management standards.

7.4.3 Electronic Data Transfer

All requests for copies of the data, electronic or hard copy, are made through the appropriate project manager. All data at SFWMD is public information. The project manager is responsible for requesting the data and submitting it to the requestor. For all data transmitted electronically, (e.g. by modem, or diskettes), hard copies are sent via mail following data transmittal.

7.4.4 Documentation and Verification

All LIMS documentation is maintained by the Staff Programmer Analyst who is responsible for maintaining the documentation records and the maintenance logs. The documentation for all instrument software is located in the laboratory at the instrument. Software problems are included in the maintenance log for the instrument. Software revisions are installed and the records maintained by the Staff Programmer Analyst.

Software is verified by comparing the results generated by the new software to the results from the old software for at least six analytical runs. The evaluation of the results must yield no significant statistical difference. The QC results are also used to show that the software is performing correctly.

Data entry is verified by comparing the results obtained by the instrument software to the results entered into the database. Manual data entry is verified by comparing the results in the database to the results on the data entry forms.

8.0 Analytical Procedures

The procedures used by SFWMD are found in Section 5.0.

8.1 Field Screening Methods

A phosphorus screening method is used in the field to determine samples which may be exceeding established limits. The screening data is used only to determine which samples should be sent to the laboratory for testing (see appendix C). The criteria used to determine which samples are to be sent to the lab are project specific.

8.2 Laboratory Glassware Cleaning and Storage Procedures

For classics and nutrients, the laboratory glassware is cleaned by washing with a hot Liquinox™ solution, rinsing with D.I. water, rinsing with 10% hydrochloric, and finally rinsing with D.I. water. Once dry, glassware is stored capped, in cabinets, in the appropriate analytical or digestion area.

Digestion tubes, beakers and other pieces of glassware are washed in a labware washer which has been programmed to follow the above procedure except for the acid rinse. The program is stopped after the first tap water rinse, the glassware is removed from the washer and rinsed with 10% hydrochloric acid and then returned to the washer for the final tap and analyte free water rinses. The glassware is then allowed to either air dry or is placed in a glassware dryer. Once dry, the glassware is stored in clean drawers in the appropriate analytical or digestion area. TKN digestion tubes are not dipped in acid.

A complete supply of glassware is dedicated for use in the metals laboratory. Pipets are soaked in 1% Liquinox™, and cleaned in a pipet washer using three volumes of D.I. water followed by soaking in 10% HNO₃, and washing with three volumes of analyte free water. Volumetric flasks and glassware dedicated to major cation analyses are rinsed three times with analyte free water after each use, air dried, and stored in cabinets in the appropriate analytical or digestion area. Glassware dedicated to trace metals is rinsed with analyte free water after each use, soaked in 20% HNO₃, rinsed with analyte free water, allowed to air dry, and stored separately from other glassware in the metals analysis area. Polyethylene trace metal bottles are soaked in 20% nitric acid for a minimum of 24 hours, analyte free water rinsed, and stored in cabinets in the metals area filled with analyte free water which has been acidified to approximately 1% with double distilled nitric acid.

8.3 Laboratory Method Modifications

The color procedure has been modified for use in the laboratory. The samples submitted to this laboratory are from natural surface and ground water sources within the boundaries of the SFWMD.

The color in the samples is due primarily to vegetative decay and not from industrial sources.

Measurement of the color at 465 nm gives results comparable to those measured visually by technicians. Use of the spectrophotometer eliminates the natural variation in color perception found in the human eye allowing the lab to consistently report results regardless of which technician performs the analysis. The procedure used is given in Appendix A.

The adopted silica method is a modified method tested by Alpkem for use in Rapid Flow Analyzers. The adopted method reduces saltwater interferences and is found in Appendix B.

Alkaline phosphatase activity is measured using the Methylumbelliferyl Phosphate (MUP) assay (Petterson and Jansson, 1978). The procedure is included as Appendix D.

Total and dissolved inorganic carbon (TIC/DIC) is determined using a high temperature combustion infrared detection system following the EPA Method 415.1 for total organic carbon. The procedure is included as Appendix F.

8.4 Laboratory Reagent Storage

The storage of the laboratory reagents and chemicals is given in Table 8-1.

Table 8-1. Reagent and Chemical Storage

<u>Chemical</u>	<u>Method of Storage</u>
Mineral acids	Stored in original containers in vented cabinet designed for acid storage. Note: each type of acid is stored in a separate cabinet.
Liquid bases	Stored in original containers in a vented cabinet designed for corrosive storage.
Organic solvents	Stored in original containers in a vented cabinet designed for flammable storage in the outside storage area.
Compressed gases	Stored in original containers in the compressed gas storage area in the outside storage area.
Dry chemicals	Stored in original containers segregated by reactivity in the dry chemical storage area.

As each chemical is received, it is dated and initialed by the person unpacking it. When a new container is opened for use, it is dated and initialed by the person who opened it.

8.5 Waste Disposal

The laboratory has a designated hazardous waste storage area outside the laboratory. The process wastes containing mercury and phenol are collected for disposal by a hazardous waste company. As each waste is generated, the volume is entered on the monthly hazardous waste report required by the SFWMD Risk Management Division. When the volume collected reaches the level set by Risk Management, the waste is removed by a commercial waste hauler and disposed of according to state and federal guidelines. The Risk Management Division is responsible for securing the services of the waste disposal company.

Any small amounts of reagents are transferred to the hazardous waste storage area as they expire or are no longer needed. The waste disposal company then picks them up for proper disposal.

Concentrated acids and bases are neutralized and put into the sanitary sewer system with copious amounts of water. Samples are disposed of by washing them into the sanitary sewer system with copious amounts of water. Nitric acid is disposed of according to local, state and federal regulations.

holders are then washed following the same procedure as inorganic sample bottles. Disposable filter holders are purchased from Poretics and are not cleaned before use. At least three 10 ml portions of sample water are passed through the filter holder and filters and discarded before filling sample containers.

Filtering syringes are washed with Liquinox™, then rinsed with hot tap water, followed by analyte free water and allowed to air dry.

Clean bottles/equipment are collected from the lab by field personnel and are not shipped. Bottles for the collection of samples sent to contract laboratories are provided by the laboratory performing the analyses. These bottles have been cleaned by the contract laboratory according to that laboratory's procedures prior to shipment to the District. The laboratory must have approved cleaning procedures in their Comprehensive Quality Assurance Plan.

6.3.7 Autosamplers

All autosampler bottles are cleaned in the laboratory using the standard bottle cleaning procedure described in Section 6.3.1 as per EPA SOP & QAM guidelines, Appendix B, Sections B.5 and B.6. Tubing is not cleaned, but rather, replaced. The sampler and sampler tubing is dedicated to the site. The sampler is programmed to rinse the tubing twice with sample before collection and then purge the tubing after the sample is taken. For the American Sigma Model 6201 autosamplers located at pump stations, the sampler and pump tubing are replaced with clean tubing every four storm events, or, at a minimum, quarterly in the absence of storm events. The pump tubing for American Sigma Models 700 and 800 autosamplers is replaced at a minimum quarterly, or preferably, monthly. The entire sampler (American Sigma Model 6201 700, or 800) is cleaned when it is removed from the site for repair or transferred to another site. In an emergency, pump tubing may be cleaned by soaking in hot water and Liquinox[®] and scrubbing outside and ends with a small bottle brush. The tubing is then rinsed with hot tap water then DI water. If only metals are required the tubing will be rinsed with 10% Nitric Acid. Finally, the tubing is rinsed with DI water and both ends capped with aluminum foil until ready for installation.

6.3.8 Atmospheric Deposition Equipment

Atmospheric Deposition dry and wetfall buckets are cleaned according to section 6.3.1. Once dry, the buckets are capped and taken to the site. Each lot of cleaned buckets is tested for proper decontamination by pouring one liter of DI water into the bucket, swirling, allowing to steep overnight, and analyzing for the parameters of interest.

6.3.9 Low Level Mercury Equipment

All low level mercury sample collection equipment (including sample bottles) is precleaned before use by the contract laboratory using approved methods. The filtration units are disposable. The

SFWMD does not clean any low level mercury equipment. Sample bottles are rinsed three times in the field with sample before the sample is collected. When using the pump, a minimum of three sample hose volumes (200 ml) of sample water are flushed through the tubing or filter before the sample is collected. The contract laboratories ship all equipment clean and in bags. All equipment is taken into the field in plastic bags.

6.3.10 Analyte Free Water Containers

The analyte free water containers are cleaned in the laboratory using the standard bottle cleaning procedures described in Section 6.3.1 as per EPA SOP & QAM guidelines, Appendix B.

6.4 Sampling Procedures


Water sampling locations in the SFWMD are frequently established at water control structures where known flow rates can be combined with chemical concentrations to determine loadings. However, in lakes, rivers, marshes, estuarine systems, storm water runoff, and agricultural point sources, sampling sites are chosen based on other criteria as described in the appropriate Quality Assurance Project Plans. In any case, the sample must be collected in such a manner as to ensure that it is representative of the water body being studied. A new pair of disposable latex/PVC gloves are used at each sampling point for all types of sampling.

6.4.1 Surface Water (except organics and autosamplers)

The surface water sampling procedures (except organics and autosamplers) are given in the flow chart in Figure 6-1. The following special considerations are observed when applicable:

1. If a boat is used, the sample is taken from the bow, and/or upwind and upstream from the motor.
2. When wading, the sample is collected upstream from the boat or body.
3. Care is taken not to disturb the sediment in the immediate sampling area.
4. Pre-preserved containers are not used as collection containers.
5. Intermediate containers are inverted, immersed to the appropriate depth, and turned upright pointed in the direction of flow, if applicable.
6. Preserve samples according to Table 6.10.

Figure 6.1 - Inorganic Surface Water Sampling Procedure



*Preservation/Filtration occurs within 15 minutes of collection

6.4.2 Surface Water - Organics

Surface water samples for organic analysis are collected using a subsurface grab sampler equipped with a precleaned glass bottle provided by the contract laboratory. All unpreserved containers (except for VOCs) are rinsed three times before sample collection. The sample containers are immersed inverted to 0.5 meter below the surface, turned upright, pointing in the direction of flow, until full.

For Volatile organic samples, the water sample is poured slowly down the edge of the 40 mL Teflon-lined septum glass vial from the precleaned glass bottle to minimize aeration. The vial is filled to the point of creating a convex meniscus. The septum is placed, teflon side down, on the meniscus and sealed with the screw cap. The vial is inverted and lightly tapped on the lid to dislodge any entrapped air bubbles. The absence of air bubbles indicates a proper seal. If air bubbles are present, the bottle is opened, additional sample is added and the vial is resealed. Additional sample is added a maximum of three times. If a seal cannot be obtained, the vial, sample, and septum are discarded, a new vial is used and the sampling procedure is repeated.

Each sample is identified by the project code and sequential sample number. During sample collection, the date, time, location, water movement, weather and site conditions are recorded. The preservation technique and holding times are shown in Table 6.10. Following collection of the sample, the bottles are sealed, tagged, and placed in wet ice.

The samples are sent to the certified contract laboratory having a DEP approved Comp QAP with appropriate sample identification and chain of custody form. The samples are packed in coolers with bubble-wrap or other appropriate packing material to avoid breakage. Samples are kept at 4°C with wet ice. Samples for VOC analysis are placed into separate bubble-pack bags for each station. At least one trip blank must be included in each separate cooler. Coolers are taped shut with shipping tape, labeled appropriately and shipped to the laboratory using common carrier overnight delivery.

6.4.3 Surface Water - Autosamplers

Water quality autosamplers are used for the collection of daily composite or discrete samples at selected sample sites. The sampling procedure is controlled by the automatic sampler programming. The automatic sampler, American Sigma Model Streamline 700 or Model 800, is programmed to rinse the sample collection tubing twice prior to sample collection and to purge the tubing following collection. The samplers are programmed depending on the project requirements. For example, the programming may include collection of daily composites at a rate of 80 mL of sample at 144 minute intervals and addition of the sample to the correct composite sample bottle for a total of 10 samples per bottle. When samples are analyzed for total Kjeldahl nitrogen, total phosphorus, ammonia and nitrate plus nitrite only, one ml of 50% sulfuric acid is pre-added to the discrete autosampler bottles before sample collection. This amount of acid has been determined to maintain pH <2 after sample collection.

After sample collection the bottles are capped, thoroughly mixed by inverting five times, labeled, placed on ice immediately and transported to the field laboratory and stored at 4°C. Within 24 hours, aliquots of the samples are transferred to 175 mL bottles, tagged, and stored at 4°C until transported to the analytical laboratory. Samples are submitted to the analytical laboratory within two days of the day the sampler was serviced. If possible, samples are poured into sample bottles at the collection site, preservation is checked and samples are immediately placed on wet ice for transport to the analytical laboratory. These procedures are only followed for long holding time parameters.

Short holding time parameters are not routinely collected using an automatic sampler. However, if short holding time parameters are collected, autosampler bottles are not preacidified. The samples are collected within 24 hours after sample uptake or compositing by the autosampler has been completed, and analyzed within holding times.

Sampler intakes at water control structures are located 1 to 2 feet below the surface and 1 foot off the wing wall.

American Sigma Model 6201 autosamplers located at pump stations are dedicated to the site and work in synchrony with each of the pumps in the pumping battery. The autosampler is activated once any of the pumps becomes operational, and is usually programmed to collect a 50 ml aliquot for every 10,000 RPM on each operating pump. The aliquot is dispensed into a refrigerated five gallon polyethylene bottle. The refrigerated autosampler is deployed for up to one week and samples picked up within 24 hours of after sample collection has been completed. Deployment times may vary depending on project requirements as specified in the project specific Quality Assurance Plan. The refrigerated composite sample is homogenized by shaking the capped bottle. An aliquot is transferred unfiltered to a 250 mL sample bottle for TKN and TPO₄ and another aliquot is filtered into a 60 mL sample bottle, for NO_x and NH₄ testing. Both aliquots are preserved as indicated in Table 6.10.

Autosamplers are dedicated to a specific location as long as they are working properly. If an autosampler fails to perform, the sampler head including the pump hose, electronics, and flow sensor are replaced as required. Sampler intake tubing is dedicated to a collection site.

6.4.4 Sediment/Soil

Sediment samples for organics are collected by hand grab using a stainless steel scoop/spoon, stainless steel petite Ponar™ dredge, or stainless steel core. The stainless steel scoop/spoon is utilized only in quiescent shallow waters. The petite Ponar™ dredge is effective over the wide range of circumstances encountered during the collection of sediment samples. The petite Ponar™ dredge can be used on soft or hard bottoms under a variety of flow and depth conditions. Sediment samples from canals are collected mid-stream using a small boat. The dredge is lowered by rope until contact with the bottom is established. Then the dredge is slowly retrieved and emptied into a stainless steel bowl/tray. Samples are not composited or sectioned into subsamples. The stainless steel core is used when collecting shallow sediment

samples. Each sample (except VOC) is thoroughly mixed in the stainless steel bowl with a precleaned stainless steel spoon. The sediment in the bowl should be scraped from the walls and bottom, rolled to the middle and initially mixed. The mass is quartered, moved to the opposite sides of the bowl, and each quarter individually mixed. Each quarter is then rolled to the center of the bowl and the entire mass is mixed again. This process is continued until the mass is as homogeneous as possible. A precleaned stainless steel spoon or small scoop is used to transfer the required amount of final homogenized material from near the center of the mass into the appropriate sample bottle. This is done quickly to prevent oxidation of metal ions or volatilization of organic compounds. The bottle is filled to its rim.

Duplicates are collected to measure the variability inherent in the sampling process. Duplicates for sediments are collected from the same sampling device. Since a true split sediment sample is almost impossible under field conditions, split soil samples are considered duplicates. If analyses for pesticides, extractable organic compounds, or VOC are to be performed, the sample containers must be glass with teflon lined lids. VOCs bottles must have a teflon-lined septum. For other classes of analytes, glass or plastic jars may be used as sample containers. The containers are completely filled with sample so that air bubbles are not trapped in the container. This is necessary in order to minimize sample oxidation that could influence certain test results. The best way to minimize trapped air is to very gently pack the sample into the container with the spoon or scoop. Preserve according to Table 6.12.

Sediment samples are collected using PVC cores if the sample will not be analyzed for organics. Intact sediment cores are obtained by driving a PVC coring tube to a depth of approximately 50 cm into the soil. Under shallow conditions (water depth < 30 cm) this is achieved using a block of wood to protect the neck of the core tube and then striking the block several times with a hammer such that the core penetrates the soil with minimal compaction. If the water is greater than 30 cm, a PVC coupling should be attached to the coring device. This attachment is comprised of a ball valve and a length of PVC pipe. This is fitted to the neck of the PVC coring tube and a closed ended piece of metal pipe with a diameter larger than the upright PVC pipe is used to pound the core into the sediment. During the insertion, the ball valve should be open.

6.4.5 Ground Water

Ground water samples are collected from the least to most contaminated wells. The order of collection is: 1) VOC, 2) Extractable Organics, 3) Total Metals, 4) Dissolved Metals, 5) Inorganics and 6) Radionucleides. The SFWMD does not use temporary well points or dedicated equipment for ground water monitoring. Samples are not collected from wells with free product showing.

A protective covering of visqueen plastic is placed on the ground at all sampling sites to reduce the potential for contamination. A fresh pair of disposable non-powdered latex gloves are worn at each sample site while purging and another fresh pair worn for collecting samples.

All samples are collected using a Teflon bailer which is suspended by a Teflon coated stainless steel lanyard and carefully lowered into the well. Disposable lanyards are not used. Prior to the collection of the sample, the bailer is rinsed three times with sample. The lanyard and bailer are not allowed to touch the ground during sampling or purging.

The bailer, filled with sample, is retrieved from the well and placed on a bailer stand. VOC samples are collected first via a controlled-flow bailer bottom and are poured slowly down the side of the sample vial to minimize sample aeration. The vial is filled to the point of creating a convex meniscus. The cap is secured with the Teflon side of the septum contacting the sample and no headspace in the sample container. The vial is inverted and gently tapped to locate bubbles. If bubbles are present, the sample in the vial is discarded and sample is recollected and checked for bubbles. Additional sample is added a maximum of three times. If unsuccessful, the sample, vial and septum are discarded, a new vial and septum are used, and the collection procedure begins again.

To filter a sample, a one piece, molded in-line disposable filter is attached to the bottom of the bailer with a Teflon adaptor. The top of the bailer is attached to a vacuum pump by a Teflon attachment so the bailer can be pressurized to force water through the filter. The filters are 0.45 micron unless a dissolved metals sample is collected for permit compliance, in which case the filter pore size is required to be 1.0 micron. The first 100 mls of sample to pass through the filter is discarded as rinse water, and the last 100 mls of sample water in the bailer is not used since it has been in contact with the air at the top of the bailer.

Bottles are supplied by the SFWMD laboratory. If the analyses are to be performed by a contract laboratory, the bottles are cleaned by the contract lab according to that laboratory's procedures prior to shipment to SFWMD. The contract laboratory must have approved cleaning procedures in their Comprehensive Quality Assurance Plan, and perform a 5% analytical check on the washing procedure. Certified disposable VOC vials are purchased from I-Chem.

Duplicates are collected by sampling from consecutive bailers.

Splits are not routinely collected in the groundwater program. When splits are required to measure the performance between two or more laboratories, the following procedure is used:

Splits are collected from the same bailer. For large volume samples that may require more than one bailer load, the first half volume of the first bailer load is poured into the first set of containers and the second half in the second set of containers. Then the first half-volume of the second bailer load is poured into the second set of containers and the second half in the first set of containers, etc., until both sets are full.

The sample bottles are labeled with the date, sample number and project name before sampling at a site begins. After collection the samples are preserved, and immediately placed on ice in a closed container. Exposure of the organic samples to sunlight is kept to a minimum. When VOCs are collected, trip blanks are kept with the organic sample bottles to make certain that the samples have not become contaminated. These trip blanks accompany the samples from the time the empty sample bottles are shipped from the contract laboratory until the samples are analyzed. Trip blanks account for at least 5% of the samples that are analyzed. At least one trip blank must be included in each separate cooler.

6.4.5.1 Wells with In Place Plumbing

These wells are purged for a minimum of 15 minutes, until three bore volumes have been removed, or until the well has chemically stabilized, whichever is greater. The sample is taken from the faucet closest to the source and before any screens, aerators, or filters. The flow rate is reduced to less than 500 ml/min when collecting samples to avoid any undue disturbance. Unfiltered samples are collected directly into the sample containers from the spigot. Filtered samples are collected by filling a bailer from the spigot. The bailer is rinsed three times with sample, then filled and pumped through a filter attached to the bottom of the bailer with a Teflon adaptor. Sample bottles are rinsed once with sample (unless sample bottles are prepreserved, collection is for Oil & Grease, bacteriologicals, Total Recoverable Petroleum Hydrocarbons or VOCs).

6.4.5.2 Flowing Artesian Wells

These wells are purged until three bore volumes have been evacuated, or until the well has chemically stabilized, whichever is greater. The flow of water from the well is adjusted to minimize the aeration and disturbance of samples. Unfiltered and filtered samples are collected directly from the discharging water into a bailer.

6.4.5.3 Monitoring Wells

The depth to water in each well, relative to a measuring point is measured using a graduated stainless steel tape and chalk or an electronic water level indicator. The depth to water is measured twice. Both values, which must be accurate to within 1/10th of a foot, are recorded in the sampler's field log. More stringent measurements may be required for specific project programs. These values will later be used to calculate water elevation relative to mean sea level.

The same measuring point is used every time the well is sampled. When the wells are surveyed, the survey point is used for the measuring point. Water level measuring devices are rinsed with dilute Liquinox[®] and DI water and wiped dry with a Kim Wipe[®] towel before measuring each well.

Prior to the collection of samples, each well is then purged by a purge pump, submersible pump or centrifugal pump. Purge pumps are gas powered and extreme care must be used when handling and locating these pumps to minimize on-site contamination (placed away and downwind from any sampling activities). The purge pump consists of a centrifugal suction pump that does not contact the water within the well. A flexible PVC suction hose connects the pump to a three foot length of rigid Teflon drop pipe. The portion of the drop pipe that comes in contact with the water column is cleaned prior to use as per the procedures outlined for other sampling equipment in section 6.3. The drop pipe is equipped with a Teflon check valve at its bottom end to prevent the back flow of purged water into the well. The drop pipe is lowered to the top of the water column so that the purging process removes all of the standing water. If the water level is drawn down during the purging process and air begins to enter the drop pipe, the drop pipe is lowered to the new water level. If the water level is drawn down too severely, the purge rate is decreased. The SFWMD does not collect samples at any well with a depth to water of more than 25 feet. Submersible pumps are used for purging only. They are decontaminated as specified in section 6.3.4.5.

The well is purged until a minimum of three standing water volumes are removed. The volume of water to be purged is calculated using the following formula:

$$\text{Minimum Purge Volume in gallons} = (D^2) * (DW - DTW) * 0.1224^1$$

D = Casing Diameter in inches

DTW = Depth to Water in feet

DW = Depth of Well in feet

(1) The constant 0.1224 is an units conversion factor

The volume of water removed from the well must be calculated to avoid excess purging. The flow rate is estimated by measuring the amount of time required to fill a bucket of known volume. The required purge volume is then divided by the estimated flow rate to find the estimated time necessary to purge the well. Flow is measured several times during purging to be certain that it does not change. Though a minimum of three water column volumes must be purged, the well is not considered purged until the well has chemically stabilized. Temperature, pH and conductivity are monitored and readings are recorded at 5 minute time intervals, or at shorter intervals if necessary. Chemical stability readings of the well are taken at intervals required to purge one half the bore volume. These measurements are made in a flow through chamber to minimize atmospheric contact with the sample. The well is considered to be chemically stable when two consecutive readings of temperature, pH and conductivity are within 5% (at least five minutes apart). The purge volume is noted in the field log. All measurements that were recorded while waiting for the well to stabilize are recorded in the field log. Only the final readings taken after the well has stabilized are input into the sample

results database. The drop pipe is slowly raised out of the well while the purge pump is still on to make certain that all of the water above the drop pipe inlet is purged. This procedure also minimizes the possibility of the back flow of water from the drop pipe or suction hose.

Wells with low hydraulic conductivity sediments that can be purged dry at one liter per minute are not sampled because the water quality is not deemed to be representative of the aquifer water quality. The SFWMD does not sample wells, dry purged or slow recovery, for which recovery cannot be accomplished in four hours or less.

All samples are collected using a Teflon bailer. When handling bailers or sample containers, disposable latex gloves are worn. Prior to the collection of the samples, a Teflon bailer is lowered into the well via a lanyard. The bailer is filled and evacuated three times. This procedure is done as a final rinse of the pre-cleaned Teflon bailer and associated Teflon parts before actual sampling begins. The bailer is suspended by a Teflon coated stainless steel lanyard that is carefully lowered into the well. The lanyard and bailer are not allowed to touch the ground during sampling or purging.

The bailer, filled with sample, is retrieved from the well and placed on a bailer stand. VOC samples are collected first via a controlled-flow bailer bottom to avoid sample aeration. The first 100 mls of sample in the bailer are discarded. The sample bottles are filled from the bailer after discarding the first rinse, according to section 6.3.4.

6.4.5.4 Porewater Wells and Peepers

Porewater is collected using porewater wells and peepers (Appendix G) for research projects only. Once collected, samples are processed and handled as routine groundwater samples.

6.4.6 Biological Tissue

6.4.6.1 Above Ground Macrophyte Biomass

Above ground vegetation is clipped at the sediment-water interface with a pair of shears or knife, and separated into subsamples by species. Each biomass subsample is placed into an individual large opaque plastic bag, labeled, and transported to the laboratory for processing.

6.4.6.2 Below Ground Macrophyte Biomass

After the above ground vegetation has been removed from the sampling quadrant, a soil core is collected from the center of the quadrant using a PVC soil corer. Each soil core is placed into an individual large opaque plastic bag, labeled, and transported to the laboratory for processing.

6.4.7 Hazardous Wastes/Drums

The SFWMD does not sample hazardous waste or drums.

6.4.8 Waste Water

The SFWMD does not sample waste water.

6.4.9 Microbiological

Microbiological samples are collected with extra care in order to prevent sample contamination. The samples are not composited. The personnel collecting the samples do not touch the rims/top of the Whirlpaks™. A Whirlpak™ must be discarded if it is suspected that the tip portion has been touched. Microbiological sampling is closely coordinated with the laboratory performing the analyses to ensure that samples are analyzed within holding times.

6.4.10 Oil and Grease

Oil and Grease samples are collected directly into the sample container without rinsing with sample. Sample containers are not pre-acidified. Sample bottles for Oil and Grease are provided by the contract laboratory.

6.4.11 Low Level Mercury

The SFWMD follows the clean hands (CH), dirty hands (DH) technique as specified in *EPA Method 1669: Sampling Ambient Water for Trace Metals at EPA Water Quality Criteria Levels* when collecting low level total and methyl mercury samples in surface water. CH, using shoulder length plastic gloves, followed by wrist length plastic or latex gloves, touches only clean surfaces (such as new gloves, new plastic bags, or surfaces that have been cleaned in an ultra-clean facility). DH, using wrist length latex or plastic gloves, touches the cooler and the sample bottle covered with two plastic bags. DH opens the outer plastic bag, CH reaches in, pulls up bottle covered with inner plastic bag, and removes the bottle from the inner bag. CH places the inner bag into the outer bag. DH seals the outer bag.

CH reaches into the water or slowly wades (so as not to suspend sediment) and collects 100 ml of surface water by immersing the sample bottle about 10 cm under the surface. The bottle is rinsed by shaking the contents (cap on) and dumping the contents away and downstream from the sampling area. This rinse is repeated two more times. The bottle is filled a final time with more than 100 mL and the cap secured tightly by hand. CH then wades out of the water, puts the bottle in the inside bag, seals it, and places the inside bag into the outside bag. DH seals the outside bag and places it back inside the container. DH records the sample bottle number.

time, site and other information on the field data log sheets, and later transfers the information to the field notebook. Relevant site conditions are recorded into the field notebook.

Samples may also be collected using a pump. The CH/DH procedure is used. DH removes the sampling train (filter holder with nitex screen, Teflon tubing and C-FLEX tubing - all connected) from the cooler and opens the outside bag. CH opens the inside bag and removes the sampling train. While CH holds the tubing, DH secures it to the sampling boom with tie wraps. DH opens the pump head, CH places the pump tubing inside, and DH closes the pump head. DH opens the ring stand clamp, CH places the pump tubing inside and DH closes the clamp.

To begin sample collection, DH positions the end of the sampling train about 10 cm below the water surface and about two meters from shore and secures it to the sampling platform (i.e. place boom end under battery handle). DH starts the pump and begins flushing. DH removes a sample bottle from the cooler and opens the outside bag. CH opens the inside bag and removes the sample bottle. CH places the inner bag into the outer bag. DH seals the outer bag. CH opens the sample bottle. The contents are dumped into a waste container if they are dilute HCL or on the ground (away from the sampling area) if contents are deionized water. After flushing with three sampling-hose volumes of sample water (approximately 400 mL), CH fills the bottle with about 50 mL from the sample water stream. Note: keep cap faced downwards. The bottle is rinsed by shaking the contents (cap on) and dumping the contents away from the sampling area. This rinse is repeated two more times. The bottle is filled a final time with more than 100 mL and the cap secured tightly by hand. CH then puts the bottle in the inside bag, seals it, and stuffs the inside bag down inside the outside bag. DH seals the outside bag and places it back inside the cooler.

To prepare a filtered surface water sample, DH stops the pump, removes a filter cartridge from a cooler and opens the outside bag. CH opens the inside bag, removes the filter and connects it to the end of the sampling train. DH starts the pump and the filtered sample is collected as above. Equipment is rinsed with sample water at the next site before collecting a sample. To control mercury contamination, CH secures a clean plastic bag to each end of the sampling train for transport between sites. The equipment and supplies are stored inside coolers and bags in a place relatively free of mercury contamination. Dust and human breath (from dental amalgams) contain large amounts of mercury and should be avoided during sample collection. Samples should not be collected in the rain.

6.4.12 Benthic Macroinvertebrates

The SFWMD follows macroinvertebrate sampling methods as detailed in the Florida Department of Environmental Protection (FDEP) Biology Section Standard Operating Procedures (Appendix E) and in EPA/600/0-90/000, *Macroinvertebrate Field and Laboratory Methods for Evaluating the Biological Integrity of Surface Waters*. The field methods used to study benthic macroinvertebrate communities relate to qualitative, semi-quantitative and quantitative sampling, and are dependent on the data quality objectives of the study. Generally, there are four categories of macrobenthic invertebrate samples collected: cores, grabs, artificial

substrates, and miscellaneous. Field equipment is cleaned prior to the field trip and after each sample according to the procedures described in Section 6.3.5.

a. Corers

Coring devices may be used at various depths in any substrate that is sufficiently compacted so that an undisturbed sample is retained; however they are best suited for sampling homogeneous soft sediments, such as silt, clay or sand. Macrofauna are sampled using single and multiple-head corers, tubular inverting samplers, open-ended stovepipe devices, hand corers and/or box corers. Devices are lowered slowly into the substrate to ensure good penetration and to prevent organism escapement. Visual inspection of each sample is necessary to ensure an adequate amount of sample is obtained.

b. Grab samplers

Grab samplers are designed to penetrate the substrate by gravity and have spring- or gravity-activated closing mechanisms. They are used to sample a unit area of the habitat. The habitat and substrate type sampled, depth of penetration, angle and completeness of jaw closure, loss of sample during retrieval, disturbance at the water-sediment interface, and effect of high flow velocities all affect the quantity and species of macroinvertebrates collected by a particular grab. Petite Ponar™ or Ekman grab samplers are typically used; however the type and size of the selected device depends on the substrate composition, water depth, and hoisting gear available. The sampler is lowered slowly to avoid a disturbance of the surface sediment and to ensure the device bites vertically. Upon tripping of the closing mechanism, the sediment is contained and the device is retrieved. The sample is then placed into a suitable container for transporting to the laboratory or is placed directly into a sieving device and processed.

c. Artificial substrate samplers

Artificial substrate samplers, such as the multiplate (modified Hester-Dendy) sampler, consist of natural or artificial materials of various composition and configuration. They are placed in the waterbody for a predetermined period of exposure and depth for the colonization of indigenous macroinvertebrates. This type of sampling is used to augment bottom substrate sampling because many of the physical variables encountered in bottom sampling are minimized, e.g. light attenuation, temperature changes, and substrate variation. The samplers are deployed in the euphotic zone of the water column to promote maximum colonization. They are exposed for a preset period of time and retrieved vertically from the water. The samplers are placed in a preservative container and transported to the laboratory for processing.

d. Miscellaneous qualitative devices

Many devices such as dip nets, bare hands, tongs, and forceps may be used to collect benthic macroinvertebrates. Dip nets are used by sweeping the net through the water or by holding the net stationary against the bottom and disturbing the substrate, causing the benthos to be swept into the net.

Samples collected by any of the aforementioned devices contain varying amounts of fine materials which can be removed by sieving immediately in the field, or in the laboratory. If

laboratory sorting can be performed within 24 hours. place samples immediately on ice. If sorting will be delayed, preserve with 10% formalin. After organisms have been removed from detritus, they should be placed into 70% ethanol.

6.4.13 Atmospheric Deposition

Atmospheric deposition samples are collected according to the National Atmospheric Deposition Program (NADP) guidelines.

The District uses Aerochem Metric Wet/Dry Precipitation collectors. The Aerochem is an aluminum structure with two 3.5 gallon PVC buckets which serve as receptacles for wet and dry deposition. An aluminum lid covers the wet side bucket under dry conditions. When the humidity sensor detects rain, the lid slides over the dry side bucket, uncovering the wet side bucket. All buckets are collected on Tuesdays only, according to NADP guidelines. The buckets are capped (while using PVC/Latex gloves) and taken to the laboratory for processing within 24 hours. Clean buckets are placed into the Aerochem collectors. Before clean buckets are placed into the collectors, one bucket from each decontamination lot is checked by pouring one liter of DI water into the bucket, allowing the water to sit overnight (capped), and processing as an equipment blank.

All wet side buckets are weighed using a top loading balance and subtracting the clean bucket tare weight (determined before deployment in the field). Live frogs, insects and nonrepresentative foreign matter are removed with a clean scoop or tweezers. Any removed contaminants are documented in the field notebook. The sample is then processed as a routine surface water sample, including filtration and preservation as required. Processing equipment blanks, splits and field blanks are submitted according to routine QA procedures.

Dry side buckets are inspected for contamination as above. Only loose contamination is removed (with clean tweezers/spatula). Approximately 600 ml of DI water is added to the dry bucket contents rinsing the sides of the bucket. A precleaned spatula is used to rub the inside bottom and sides of the bucket. The total amount of DI water added should not exceed one liter and the volume must be recorded in the field notebook. The sample is then processed as a routine surface water sample, including filtration and preservation as required. Processing equipment blanks, splits and field blanks are submitted according to routine QA procedures.

6.4.14 Duplicates/Split Samples

Duplicates are collected to measure the variability inherent in the sampling process, and should be obtained by duplicating in rapid succession during the same sampling event, the entire acquisition technique used to obtain the first sample. Duplicates for water are collected by sampling from successively collected volumes. Duplicates for soils are collected from the same sample source (i.e., soil obtained from the same soil sampling device). Split samples are collected to measure the variability between laboratories, and should be obtained as subsamples

taken from the same parent sample (a true split of soil, sediment or sludge is almost impossible to accomplish under field conditions; a split soil sample should be considered a duplicate and not a split sample). Split samples (for water) are taken from consecutive sample volumes from the same sampling device (i.e. from same bailer), or by mixing in a large intermediate vessel.

Note: for large volume samples that may require more than one bailer load, the first half volume of the first bailer load is poured into the first set of containers (the second half into the second set of containers), then first half volume of second bailer load is poured into the second set of containers (second half into first set of containers) etc. until both sets are full.

6.5 Documentation

The following is a list of the field records that are maintained.

1. Field Trip Checklist (surface water/groundwater)
2. Chemistry Field Data Log
3. Bound field notebook (project specific)

6.6 Preservation

Holding times and preservation techniques for each parameter are given in Tables 6-10, 6-11, 6-12, 6-13 and 6-14.

Samples are preserved in the field at the time of sample collection either by using preservatives provided by the SFWMD laboratory or the contract laboratory. If sample containers are received with preservative from the contract laboratory (VOCs), additional preservative is requested so it is available if needed. ACS reagent grade or better preservatives are used.

For unfiltered samples, the appropriate preservative is added to the sample bottle after the bottle is filled. The bottle is capped and shaken after which a small amount of the sample is poured onto the pH (0 - 3 pH units) test strip to ensure $\text{pH} < 2$. If pH is not < 2 , additional acid is added drop-wise, the bottle is capped and shaken, and the pH is tested again. This procedure is followed until $\text{pH} < 2$. The pH is checked on all samples requiring preservative and this amount of acid is added to the equipment blank and field blank. The amount of acid added is recorded in the field notebook.

For filtered samples, the acid is added after filtration following the procedure outlined for unfiltered samples.

Fresh preservatives are obtained from stocks weekly for all sampling trips. Preservatives are taken into the field in polyethylene dropper bottles. One mL nitric acid (70%) glass ampoules are used in groundwater sampling for the preservation of metals samples.

6.7 Sample Dispatch

Samples to be analyzed by the SFWMD laboratory are submitted to the laboratory by field personnel. The majority of samples are submitted the same day they are collected but all samples are submitted as soon as possible after collection in order to meet recommended holding times.

Samples to be analyzed by contract laboratories are shipped to the laboratory by common carrier overnight delivery the same or next day they are collected. All samples are carefully packed with appropriate material to prevent breakage. Insulated coolers are used for sample shipment and are sealed with shipping tape to avoid tampering. If samples must be kept at 4°C, wet ice is used during shipping.

6.8 Field Waste Disposal

All field generated wastes and contaminated purge waters are segregated and containerized for proper disposal including the use of the services of a commercial disposal company, if required. The SFWMD does not sample hazardous waste sites so the only field generated wastes are acids and isopropyl alcohol. The isopropyl alcohol is containerized for proper disposal. Acids are diluted and/or neutralized and flushed down the sanitary sewer. The calibration standards for field parameters are flushed into the sanitary sewer.

Table 6.10 Holding Time and Preservation for Water Inorganics¹

Parameter	Holding Time	Preservation	Container & Size
Alkaline Phosphatase	24 hours	Cool. 4°C	Plastic, 175 ml
Alkalinity	14 days	Cool. 4°C	Plastic, 175 ml
Ammonia	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Plastic, 60 ml
Biochemical Oxygen Demand	48 hours	Cool. 4°C	Plastic, 500 ml
Carbonaceous Biochemical Oxygen Demand	48 hours	Cool. 4°C	Plastic, 500 ml
Chemical Oxygen Demand	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Chloride	28 days	None required	Plastic, 60 ml
Chlorine, Total Residual	Analyze immediately	None required	Plastic/glass, 250 ml
Chromium VI	24 hours	Cool. 4°C	Plastic, 250 ml
Color	48 hours	Cool. 4°C	Plastic, 60 ml
Fluoride	28 days	None required	Plastic, 500 ml
Inorganic Carbon	14 days	Cool. 4°C	Plastic, 175 ml ²
pH	Analyze immediately	None required	Plastic, 175 ml ³
Kjeldahl nitrogen	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Mercury	28 days	pH<2(HNO ₃)	Plastic, 125 ml
Metals	6 months	pH<2(HNO ₃)	Plastic, 250 ml
Nitrate + nitrite	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Plastic, 60 ml
Nitrite	48 hours	Cool. 4°C	Plastic, 60 ml
Oil and Grease	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Glass, 1 Liter
Organic Carbon	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Orthophosphate	48 hours	Filter immediately, cool. 4°C	Plastic, 60 ml
Total phosphorus	28 days	Cool. 4°C, pH<2(H ₂ SO ₄)	Plastic, 175 ml
Residue, filterable, nonfilterable, volatile	7 days	Cool. 4°C	Plastic, 500 ml ³
Silica	28 days	Cool. 4°C	Plastic, 60 ml
Specific conductance	28 days	Cool. 4°C	Plastic, 500 ml ³
Sulfate	28 days	Cool. 4°C	Plastic, 60 ml
Turbidity	48 hours	Cool. 4°C	Plastic, 500 ml ³

(1) From 40 CFR, Part 136, Table II (7-1-90)

(2) Not listed in 40 CFR, Part 136, Table II (7-1-90) (3) 250 ml plastic bottle may be used if BOD is not requested

Table 6.11 Holding Time and Preservation for Water Organics

Parameter	Holding Time	Preservation ¹	Container & Size
Volatile (Purgeable) Organics	14 days (preserved with 1:1 HCl), 7 days (unpreserved)	Cool, 4°C	Glass, 40 ml, Teflon lined septum
Base neutral acid extractable compounds	7 days until extraction, 40 days after extraction	Cool, 4°C	Amber Glass, 1 L, Teflon lined cap
Organochlorine Pesticides and PCB's	7 days until extraction, 40 days after extraction	Cool, 4°C	Amber Glass, 1 L, Teflon lined cap

(1) If Residual chlorine is present, Sodium Thiosulfate is added to the sample vial first. The vial is then filled to almost full volume with sample, acid is added, and finally the vial is filled as per procedure.
 Note: It is not recommended to mix the two preservatives (and sample) together in an intermediate vessel.

Table 6.12 Holding Times and Preservation for Sediments, Soils and Tissues¹

Parameter	Holding Time	Preservation	Container & Size
Volatile organics	14 days	Cool, 4°C	Glass, 4 oz. widemouth with Teflon/silicone septum
Semivolatile organics	14 days until extraction, 40 days after extraction	Cool, 4°C	Glass, 8 oz. widemouth with Teflon/silicone septum (50 grams)
Total metals	6 months	Cool, 4°C	Glass or plastic, 8 oz. widemouth (200 grams)
Mercury	28 days	Cool, 4°C	Glass or plastic, 8 oz. widemouth (200 grams)
Nutrients & Inorganics ²	Not Specified	Cool, 4°C	Glass, 500 ml or plastic, 8 oz. widemouth with Teflon lined closure

(1) From Table 5, Chapter 62-160, F. A. C.

(2) From USEPA Standard Operating Procedures and Quality Assurance Manual (Appendix A).

Table 6.13 Holding Times and Preservation for Microbiologicals⁴

Parameter	Holding Time	Preservation ²	Container & Size ³
Microbiologicals ¹	6 hours	Cool. 4°C/Na ₂ S ₂ O ₃	Glass or Plastic > 125 ml

(1) Parameters included are: Fecal Coliform, Total Coliform and Fecal Streptococci.

(2) Addition of sodium thiosulfate is only required if the sample has a detectable amount of residual chlorine, as indicated by a field test using EPA Method 330.4 or 330.2 or equivalent.

(3) Presterilized Whirlpak bags (or equivalent) are typically used for sampling.

(4) From Table 8, chapter 62-160, F.A.C.

Table 6.14 Holding Times and Preservation for Biologicals

Parameter	Holding Time	Preservation	Container & Size
Chlorophyll	(i) 24 hours to filter (ii) 21 days after filtration	(i) Cool. 4°C, dark, filter with MgCO ₃ . (ii) Filter frozen (until testing)	Amber Plastic, 1 L
Macrobenthic Invertebrates	Preserved in the field with formalin, then ethanol in the laboratory or placed on ice and preserved with formalin within 8 hours, then preserved with ethanol	10% buffered formalin, then 70% ethanol	Glass or Plastic

6.9 Field Reagent and Standard Storage

Table 6-15 lists the storage procedures for reagents, standards and solvents.

Table 6-15. Reagent, Solvent and Standard Storage

<u>Chemical</u>	<u>Method of Storage</u>
Sulfuric Acid	Stored in original containers in vented acid storage cabinet. Note: Each acid is stored in a separate cabinet. Acid is transported in the field in polyethylene dropper bottles.
Nitric Acid	See above.
Phosphoric Acid	See above.
Hydrochloric Acid	See above.
Isopropyl Alcohol	Stored in original containers in vented solvent storage cabinet. No other solvents are stored in the same cabinet. Taken into the field in glass containers carried in a safety carrier.
pH Standards	Stored in cabinet designated for standard and reagent storage. Cabinet is in air conditioned laboratory. Standards taken into the field in polyethylene bottles.
Conductivity Standards	Stored in cabinet designated for standard and reagent storage. Cabinet is in air conditioned laboratory. Standards taken into the field in polyethylene bottles.
Formalin (10%)	Stored in cabinet designated for standard and reagent storage. Cabinet is in air conditioned laboratory. Taken into the field in polyethylene bottles.
Ethyl Alcohol	Stored in original containers in vented solvent storage cabinet. No other solvents are stored in the same cabinet. Taken into the field in approved non-combustible containers.

9.0 Calibration Procedures and Frequency

9.1 Instrumentation Lists

The following is a list of the laboratory instrumentation.

<u>Manufacturer</u>	<u>Model and Description</u>
Hach	Model 2100 AN Turbidimeter
Hach	Model 18900 Ratio Turbidimeter
Accumet	Model 50 Ion Analyzer
Fisher	pH Electrode
Orion	Combination Fluoride Electrode
Orion	Model 162 Conductivity Meter
Radiometer	CDM83 Conductivity Meter
Mettler	P160 Top Load Balance AE163 Analytical Balance (2)AE100 Analytical Balance
Bausch & Lomb	Spectronic 501 Visible Spectrophotometer
Dionex	4000i Ion Chromatograph DX500 Ion Chromatograph
Alpkem	(3)RFA300 Rapid Flow Analyzer with PC workstation (3)RFA500 Rapid Flow Analyzer with PC workstation
Perkin Elmer	1100B Flame Atomic Absorption Spectrophotometer with PC Workstation Z5100 Furnace Atomic Absorption Spectrophotometer with PC Workstation Optima 3000XL Sequential ICP with PC Workstation Lambda 6 uv-vis Spectrophotometer with PC workstation
Orion	Model 960 Autochemistry System
Orion	pH Electrode with ATC

<u>Manufacturer</u>	<u>Model and Description</u>
Shimadzu	Model 5050A Total Organic Carbon Analyzer
Rosemount Dohrmann	DC-190 Total Organic Carbon Analyzer
Millipore	Cytofluor 2350 Fluorescence Measurement System
Thermolyne	Type 37900 Culture Incubator
Fisher Scientific	Model 21K/R Marathon Centrifuge Model 307A BOD Incubator
Barnstead	(3) Model 2250 Autoclave
Fisher Scientific	Model SPT - IH Stereoscope
Reichert	Quebec Dark field Colony Counter
Precision	Model 6M Incubator Model 251 Coliform Water Bath
YSI	Model 59 DO Meter Model 5905 DO Probe
Hach	Model 45600 COD Reactor
Thermolyne	(3) Model 9000 Oven Furnatrol II Muffle Furnace
Thelco	GCAL Precision Scientific Oven
CEM	MDS - 2100 Microwave Digester
Tecator	(2) Digestion System 40, 1016 Digester
Powers Scientific	(2) Refrigerator
Amerikooler	Refrigerator
Jordan	Refrigerator

The following is a list of the field instrumentation.

<u>Manufacturer</u>	<u>Model and Description</u>
Hydrolab	Model 4031 pH, Conductivity, ORP, and Temperature Meter Model 4041 pH, Conductivity, DO, and Temperature Meter Model Surveyor II pH, Conductivity, DO, Temperature, Salinity, Model Surveyor III pH, Conductivity, DO, Temp., Salinity
Solomat YSI	ORP, Turbidity and Depth Meter WP803 Model 58 Dissolved Oxygen Meter Model 6000UPC Multiprobe pH, Conductivity, DO, Temperature, EH Turbidity
HF Scientific	Model 600XL Multiprobe pH, Conductivity, DO, Temperature, EH Model DRT-15CE Turbidity

9.2 Standard Receipt and Traceability

Standards are received by the Senior Chemists, Sr. Lab Techs or Sr. Scientific Techs., initialed, dated, and stored in the designated area for the particular standard.

The preparation date of in-house primary stock solutions is recorded in a log book along with the following information regarding purchased stock solutions: analyte, concentration, supplier, date opened, expiration date and date of disposal. Only one bottle of each purchased analyte stock solution may be in use at one time. Purchased stock solutions are replaced according to expiration date or sooner if the stock is depleted. Manufacturers certificates of analysis and/or records of traceability for purchased stock solutions are filed in a notebook according to analyte or analytical category. The same procedure is followed for field samples.

Primary stock solutions prepared from freshly dried, ACS reagent grade chemicals are refrigerated and prepared monthly. Records of the preparation are kept in a logbook. Working calibration standards are prepared fresh daily. The standard sources and preparation are given in Table 9-1.

The calibration procedures for laboratory instruments are given in Table 9-3, for laboratory equipment in Table 9.4, and for field instruments in Table 9-5. QC check standards from a different source than the calibration standards are used to check the initial calibration for both laboratory and field instruments.

Calibration information for the field is recorded in the Field notebook or Calibration Book. Date and time of calibration, technician, standards used, standard results and temperature, and instrument used are recorded.

Calibration information for the laboratory is recorded in the individual instrument log and includes the analyst, date of analysis, standard values, millivolt or absorbance values for the standards, correlation coefficient, results of continuing and initial calibration standards, instrument gains, precision and accuracy results, and samples analyzed.

Table 9-1 Standard Sources and Preparation

Instrument Group	Standard Sources	How Received	Source Storage	Preparation from Source	Lab Stock Storage	Preparation Frequency
Atomic Absorption/ Emission	Spex/JT Baker/ NIST/EM Scientific	Solutions of 1000 mg/L and 10,000 mg/L	Room temperature	Primary stocks (>1 mg/L) prepared from source Working stocks	0.2% HNO ₃ at room temp. NA	Weekly or as needed Daily
Continuous Flow & Ion Chromatograph	Commercial lab supplier	Dry, ACS reagent grade	Room temperature	Primary stocks, 1000 mg/L prepared from source Working stocks	Refrigerator NA	Monthly Daily
Organic Carbon/ Inorganic Carbon	Commercial lab supplier	Dry, ACS analytical grade	Room temperature	Primary stocks, 1000 mg/L prepared from source Working stocks	Refrigerator Working Stocks	Monthly Daily
Alkaline Phosphatase	Sigma Chemicals	Dry, Enzymatic grade	Room temperature	Primary stocks, 10 micro Molar in Tris Buffer Working stocks	Room temperature NA	Monthly Daily
pH Standards	Commercial lab supplier	pH 4, 7, 10 solutions	Room temperature	A	NA	Replace on expiration
Conductivity Standards	Commercial lab supplier	200, 720, 1413, 2000 uS	Room temperature	NA	NA	Replace on expiration
Turbidity	Hach	Sealed Gel Standards	Room temperature	NA	NA	Annual replacement
Color	Commercial lab supplier	500 Pt-Co units	Room temperature	Working stocks	NA	Weekly
Analytical balances	Commercial lab supplier	Class S weights	Dessicator, room temperature	NA	NA	NA

Table 9-2. Solutions Requiring Standardization

Test	Standard sources	How Received	Source Storage	Preparation from Source	Standardization procedure and Criteria	Standardization Frequency	Preparation Frequency
Alkalinity Titration	Commercial Lab Supplier	0.02 N Sulfuric Acid	Acid Storage Cabinet	Used as is from supplier (certified)	Normality checked with 0.05 N sodium carbonate, must be within 5% of expected value	Weekly and each time a new lot is used	Used as is from Supplier
Chlorine Residual Titration	Commercial Lab Supplier	1 Titer PAS titrant	Refrigerator	Used as is from supplier	Standardization and certification provided by supplier	Standardization and certification provided by supplier	Prepared monthly and as needed

Table 9-3. Laboratory Instrument Calibration

Instrument	# Standards Initial Calib.	Accept/Reject Criteria - Initial Calibration	Frequency	# Standards Continuing Calibration	Accept/Reject Criteria - Cont. Calibration	Frequency
Atomic Absorption/ Emission	2-5	Corr. Coefficient > 0.995	Daily prior to use or failure of cont. calibration	2-5	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Continuous Flow	5 - 7	Linear Regression Corr. Coefficient > 0.998	Daily prior to use or failure of cont. calibration	5-7	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Ion Chromatograph	5	Corr. Coefficient > 0.998	Daily prior to use or failure of cont. calibration	1	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Carbon Analyzer	3	Concentration within 5% of known value	Daily prior to use or failure of cont. calibration	3	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
pH Meter	3	pH 7 = 0 +/- 5 mV; pH 4 = 177 +/- 10 mV; Eff = 1.00 +/- 0.05	Daily prior to use or failure of cont. calibration	3	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Conductivity Meter	3	Concentration within 5% of known value	Daily prior to use or failure of cont. calibration	3	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Turbidimeter	3	Concentration within 5% of known value	Daily prior to use or failure of cont. calibration	3	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Visible Spectrophotometer	5	Linear Regression Corr. Coefficient > 0.998	Daily prior to use or failure of cont. calibration	1	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Fluorometer	5	Linear Regression Corr. Coefficient > 0.998	Daily prior to use or failure of cont. calibration	5	Concentration within 5% of known value (mid-range)	Initial and every 20 samples
Analytical Balance	3	Weight within 2% of known value	Monthly	3	Weight within 2% of known value	NA

(1) Cel standards are checked monthly with Formazin control solutions of the same concentration, (i.e. 1.8, 18.0 and 180 NTU), to ensure integrity.
 (2) Linearity of the dynamic range checked weekly for ICP allowing the use of 2 calibration standards.
 Note: Method Calibration requirements will be followed if more stringent than those listed in the QAP.

Table 9-4. Laboratory Equipment Calibration

Equipment	Calibration	Acceptance Criteria	Frequency
Analytical Balance	Calibrated with the 100g weight, then checked with entire set of Class S weights Maintenance contract exists for calibration	All weights within 2 % of known value NA	Monthly Semi-Annually
Autoclave	Check timer with stopwatch	Timer set for 15 min. maintains 121°C for at least 15 min. The entire cycle is completed within 45 min.	Quarterly
Incubator	Spore check Temperature recorded from a calibrated thermometer. Adjustments made as needed	No visible spores +/- 0.5 °C	Monthly Twice daily
Water Bath	Temperature recorded from a calibrated thermometer. Adjustments made as needed	+/- 0.2 °C	Twice daily

Table 9-5. Field Instrument Calibration (Hydrolabs)^{1,3}

Instrument Probe	# Standards Initial Calib.	Accept/Reject Criteria - Initial Calibration	Frequency	# Standards Continuing Calib.	Accept/Reject Criteria - Continuing Calib.	Frequency
pH	2 (1 pH 7 & 1 pH 4 or 10) ²	Reading within 0.1 pH unit	Daily prior to use or failure of cont. calibration	1	Concentration within 5% of known value	At the end of the day or within 24 H of initial calibration
Conductivity	1 in the expected Range	Concentration within 5% of known value	Daily prior to use or failure of cont. calibration	1	Concentration within 5% of known value	At the end of the day or within 24 H of initial calibration
Dissolved Oxygen	Winkler titration	Concentration within 5% of known value	Annually	Saturated air	Concentration within 5% of known value	At the end of the day or within 24 H of initial calibration
Temperature	1	Concentration within 2% of known value (NIST thermometer)	Daily	1	Concentration within 5% of known value - lab calibrated thermometer	Quarterly
Turbidity	1	Reading within 5% of known value	Daily	1	Concentration within 10% of known value	Daily prior to use
Automatic Samplers	Correct sample volume verified by using graduated cylinder	Volume within 5% of programmed volume	Daily prior to use	NA	NA	Daily prior to use

(1) The Hydrolabs all have automatic temperature compensation for pH, conductivity and DO measurements. Meters are checked daily with NIST (NBS) Certified thermometers.
 (2) Buffer pH 4 or pH 10 solution used, respectively, with acidic or basic samples.
 (3) The calibration check is conducted weekly for instruments deployed in the field

10.0 Preventative Maintenance

Preventative maintenance is a necessary part of a successful quality assurance program. Time must be allocated to clean and maintain all equipment used for the collection and analysis of a sample. Equipment which is not operating properly may give unreliable results.

10.1 Field Equipment Maintenance

Field maintenance procedures are given in Table 10-1. Maintenance that cannot be performed by SFWMD personnel is done by the manufacturer or his designee.

10.2 Maintenance Documentation

The field equipment maintenance activities are documented in bound notebooks assigned to each instrument. Service reports for repairs that cannot be done by SFWMD personnel are kept on file by the Technician Supervisor.

The laboratory equipment maintenance activities, listed in Table 10-2, are documented in a bound notebook for each instrument. Service reports for repairs that cannot be done by SFWMD personnel are kept on file in a notebook located at the instrument.

10.3 Contingency Plans

The SFWMD has replacements for most critical instruments which can be used in the event of a breakdown. Instrument service contracts are maintained on all major pieces of equipment and response time is typically 48 hours. Overflow laboratories have been contracted to provide analytical services in the event of a catastrophic failure.

Table 10-1. Field Equipment Maintenance Schedule

Instrument	Specific Activity	Frequency
Hydrolabs (all models)	DO probe membrane and electrolyte changed	Quarterly/AN
	Conductivity sensors are sanded with emery cloth	Quarterly
	pH and reference electrodes cleaned with methanol	Quarterly/AN
	pH reference electrode refilled with 3M KCl All outside surfaces cleaned and rinsed with analyte free water	Quarterly/AN Daily
Pumps (Gorman Rupp 2H Centrifugal)	Check oil and add if needed	Before use
	Drain pump of water	After use
	Wipe clean of mud and grease	After use
	Change oil & filter	Quarterly
	Change spark plugs & adjust carburator	Quarterly
Autosamplers	Check battery charge & replace as needed	Before use
	Check programming	Before use
	Magnet adjustment to seek bottle #1	After use
	Check pumping volume with a graduated cylinder	Before use
	Check indicating dessicant & change as needed	Before use
	Change pump tubing	Quarterly
	Calibrate Autosampler	Quarterly
	Distributor arm tubing replacement	Quarterly
Liquid sensor cleaned	Quarterly	
Intake tubing strainers cleaned	Quarterly	
Aerochem Collectors	Check temperature of sensor plate by touching	Before use
	Remove & cap collection buckets	Weekly
	Apply a few drops of water to sensor plate to check lid operation	Before use
	Check for snug fitting lid over collection bucket	Before use
	Check temperature of sensor plate after operation to see if warm to the touch	Before use
	Wipe top and bottom of lid & air dry	Before use
	Install clean collection buckets	Weekly

Table 10-2. Laboratory Equipment Maintenance Schedule

Instrument	Specific Activity	Frequency
Atomic Absorption	Check gases	Daily
	Service maintenance	Semiannual
	Flame: Nebulizer cleaned ultrasonically in Liquinox™ solution, rinsed with tap water, dipped in 10% HNO ₃ , rinsed with DI water.	Biweekly
	Burner head soaked in Liquinox™ solution.	After use
	Furnace: Windows inspected and cleaned with isopropyl alcohol.	Daily
	Tubes and platforms inspected and changed	AN
	Furnace decontamination as recommended by Perkin Elmer	Weekly/AN
Atomic Emission	ICP: Check pump & system tubing	Daily
	Inspect Torch & RF coil for deposits or moisture and clean if necessary	Daily
	Inspect filters	Monthly
	Clean nebulizer	Biweekly
	Flush torch with 5% HNO ₃ then DI	After use
	Pump air through spray chamber	After use

Instrument	Specific Activity	Frequency
Continuous Flow -Alpkem RFA300 -Alpkem RFA500	Inspect all tubing and fittings	Daily
	Wash manifold/flow cell	Daily
	Inspect filters	Weekly
	Replace pump tubes	Biweekly
	Clean rollers & grease	Monthly
	Service maintenance	Semiannual
Ion Chromatograph	Check tubing and fittings for leaks	Daily
	Clean columns and change bed supports	Monthly
	Preventative maintenance by manufacturer	Semiannual
Carbon Analyzer	Check/replace O-rings	Weekly/AN
	Change acid	When 2/3 empty
	Replace copper & glass wool	When discolored
	Replace injection port septa	Weekly/AN
	Inspect/replace combustion tube and catalyst	Biweekly
pH Meter	Rinse electrode with DI water	Before & after use
	Add reference solution	Daily/AN
Conductivity Meter	Rinse electrode with DI water	Before & after use
Dissolved Oxygen Meter	Rinse probe with DI water	Before & after use
	Probe membrane and electrolyte changed	Quarterly/AN
Turbidimeter	Clean cuvettes	Daily
	Adjust calibration	AN

Instrument	Specific Activity	Frequency
Visible Spectrophotometer	Clean flowcell	Daily
	Change pump tubes	Semiannual
Fluorometer	Calibration Service & inspection	Annual
Analytical Balances	Clean weighing compartment	After each use
	Clean interior/exterior	Monthly
	Calibration check against class "s" weights	Monthly/ AN
	Calibration service & inspection	Semiannual
Ovens	Check temperature	Daily
	Calibrate thermometer to NIST thermometer	Annually
Refrigerators	Check temperature	Daily
	Calibrate thermometer to NIST thermometer	Annually
Digestion blocks/ COD Reactor	Check temperature	Weekly
	Clean blocks	Monthly
	Calibrate thermometer to NIST thermometer	Annually
Centrifuge	Clean holder	After use
	Clean walls	After use
Autoclave	Check water level	Before use
	Clean interior and replace water	Before use
	Check pressure during operation	With each use
Colony Counter	Adjust focus and brightness	With each use
BOD Incubator	Check temperature	Twice daily
Incubator	Check temperature	Twice daily

Instrument	Specific Activity	Frequency
Water Bath	Check temperature	Twice daily
	Change water	Monthly
Stereoscope	Replace bulb	AN
	Wipe lense	AN

AN - As Needed

11.0 Quality Control Checks, Routines to Assess Precision and Accuracy, and Calculation of Method Detection Limits

The determination of the quality of a sample is dependent on the use of quality control samples in the field and laboratory. The definitions of each type of check can be found in DER QA-001/90, Appendix C.

11.1 Field Quality Control Checks

The field quality control checks are given in Table 11.1. The DEP QC procedures confirm the precision of the sampling techniques, that the equipment is clean, and addresses the effects of the sample handling and transport. The SFWMD routinely provides known spiked solutions to take into the field at a rate of one per quarter per project. All field blanks are preserved and transported in the same manner as the samples.

The field quality control check samples consist of the following:

- A. Field Blank - a deionized water sample poured directly into the sample container on site, preserved, and maintained open until sample collection is completed for the routine sample at that site.
- B. Trip Spike - Analyte free water spiked with a known amount of stock, prepared in the laboratory, taken into the field and handled like a routine sample (i.e. placed in coolers on ice), and returned to the lab without being opened.
- C. Replicate Sample - Samples collected at the same time from the same source.
- D. Split Sample - One sample which has been divided to make two samples which are analyzed by the laboratory.
- E. Equipment Blank - after field cleaning of sampling equipment, the final deionized water rinse is collected and analyzed as an equipment blank. Equipment blanks are also collected before sample collection begins, and at a rate of one every twenty samples. Equipment Blanks are prepared by pouring one liter of DI water into the sample collection container and through each piece of sampling equipment. For trips requiring more than one liter of water, the volume required to fill the sample bottles may be used. This volume will vary by project and must be documented in the field logbook. The Equipment Blank is filtered, preserved and handled as a routine sample.
- F. Trip Blank - Analyte free water blank (VOCs only) prepared before sample containers are transported to the field. Trip blanks remain unopened and are handled in the same manner as the samples.

Table 11.1 Field Quality Control Checks

Type	# Samples /Event	Frequency (All Parameter Groups)
Equipment Blank, Precleaned Eqpt.	>20	1 blank prior to sampling, on-site and 5% of the samples
	1-20	1 prior to sampling, on-site
Equipment Blank, Field Cleaned Eqpt.	>20	1 blank or 5% of equipment cleaned, whichever is greater
	5-10	1 blank on equipment cleaned
	<5	1 on either precleaned or field cleaned eqpt.
Trip Blank (VOC only)	1 or more	1 for each volatile organic method per cooler used to transport samples
Field Duplicate/Replicate	1 - 10	1 field duplicate/replicate
	11-20	2 field duplicates/replicates
	21-30	3 field duplicates/replicates
	>30	10% of the samples
Field Blank	1 - 20	1 field blank
	21-40	2 field blanks
	>40	5% of the samples
Trip Spike	1 or more	quarterly for each project
Field Measurements QC Check Stds. (Hydrolabs only)	1 or more	1 at the end of the day or within 24 hours of initial calibration
	(Other)	1 or more

The field quality control check samples described in A and D are included for each group of samples collected for each event at a rate of 5%. These samples are submitted to the laboratory with the routine samples for that event. Equipment blanks (E) are prepared on-site before sampling begins and at a rate of 5%.

Trip spikes are submitted to the laboratory quarterly for each project. Replicate samples are submitted at a rate of 10%. Trip blanks (F) are submitted with volatile organics (one in each cooler) and are provided by contract labs or prepared when sample containers are prepared.

Trip spikes are prepared by the QA section in the laboratory by spiking analyte free water with known stock solutions, dispensing into sample bottles and preserving. The trip spikes are taken into the field and receive the same sample handling as the samples. Project managers request the trip spikes required for their trips in advance, by entering their project name and sampling date into the project calendar. The spiked solutions are prepared using purchased or prepared stock solutions. Trip spike values are unknown to the laboratory.

11.2 Laboratory Quality Control Checks

The laboratory quality control checks are given in Table 11.2. The laboratory will meet or exceed the requirements for each method if greater than the minimum requirements. If the data from quality control check samples are not acceptable, results will be reported in the QA report to DEP.

11.3 Laboratory Quality Control Checks (Species Identification)

The SFWMD in-house specimen collections for species identification are the following:

1. Herbarium
2. Estuarine larval fish, zooplankton and benthic invertebrates
3. Phytoplankton (photographic reference record from the Caloosahatchee River)

Plant species and freshwater fish identification are done in-house. Verification of identifications by outside experts is done on an as needed basis.

Estuarine species identification is done in-house. A reference sample is made for each species collected. The reference samples and 5-10% of the samples collected, are set aside for identification by an outside expert.

Species identification for freshwater invertebrates is done inhouse. The SFWMD maintains an in-house type specimen collection. The reference samples and 5-10% of the samples collected are set aside for identification by an outside expert. Ten percent of all sorting and identification is confirmed by a different SFWMD staff member. Counts should agree within 10%.

11.4 Laboratory Quality Control Checks (Microbiology)

11.4.1 Laboratory Quality Control Checks for all Microbiology Tests

QC checks for microbiology include: (see also Table 11.3)

- Annual Water Quality Test (Suitability Test) - Standard Methods 9020
- Inhibitory Residue Test for each new lot of detergent (on file) - Standard Methods 9020
- Monthly Heterotrophic Plate Count - Standard Methods 9215
- Annual Metals Test - Standard Methods 9020
- Monthly Chlorine Residual - Standard Methods 9020
- Monthly Conductivity - Standard Methods 9020

Table 11.2. Laboratory Quality Control Checks

Type	Frequency (All parameter groups)
Method Reagent Blank	1 per sample set (batch)
Matrix Spikes (spike added prior to sample preparation)	1 or 5% of samples, whichever is greater; if more than one matrix, 1 from each matrix.
Quality Control Check Samples	Blind Performance Evaluation Samples- analyzed in duplicate semiannually*
Quality Control Check Standards	Analyzed at the beginning of each analytical run to verify standard curve. One QC is also analyzed at the end of the analytical run.
Duplicate Samples	1 or 5% of samples, whichever is greater; if more than one matrix, 1 from each matrix.
Continuing Calibration Standard	At a rate of 5% of the number of samples in an analytical set (at least one in each batch is at a concentration of 1-2 times the PQL).

* If blind QC data is not acceptable, results are reported to DEP in the QA Report.

The following QC checks are done for each microbiology test.

Table 11.3 Microbiology QC Checks

QC Check	Frequency	Purpose	Acceptance Criteria
Autoclave tape	with each sterilization batch	assure complete sterilization	Tape writing visible and dark
Spore check	monthly	assure complete sterilization	No viable cultures
Incubator temperature check	twice/day	maintain proper temperature	35.0 +/- 0.5 ° C
Water bath temperature check	twice/day	maintain proper temperature	44.5 +/- 0.2 ° C
Thermometer calibration with NIST thermometer	semi-annual	assure accurate temperature readings	< 0.2 +/- ° C correction required

All materials such as filters, plates, whirlpaks™ and media are checked for sterility upon receipt or preparation using non selective broth and inoculating for 24 hours at 35.0 +/- 0.5°C. This information is documented in a bound log book. Logs are also used for media and rinse water preparation, autoclave cycles and checks, temperature checks, equipment maintenance and sample analysis.

A control blank is run at the beginning and end, as well as every ten samples, of each filtration series or sample set. Analysis duplicates are run every ten samples. Purchases positive and negative controls are run in duplicate monthly, and with each new lot of media, for all microbiology tests. The controls are taken to the confirmation/completed stage. Ten percent of all positive samples are counted by different technicians. Counts must agree within 5%.

11.4.2 Laboratory Quality Control Checks for Membrane Filtration Procedures

A control blank is run at the beginning and end of each filtration series. In addition, a carry over blank is run every ten samples. The log sheets used for the membrane filtration procedures are found in Figures 7.13, 7.14 and 7.16.

11.4.3 Laboratory Quality Control Checks for Most Probable Number Procedure (MPN)

With each MPN analysis, Lauryl tryptose broth is inoculated with ten mL of sterile phosphate rinse buffer and used as a blank control. The completed test is performed on ten percent of all coliform positive samples and at least once per quarter.

11.5 Routine Methods Used to Assess Precision and Accuracy

11.5.1 Field Reportable Data

The formulas used to calculate the precision and accuracy of the QC checks are:

Percent Relative Standard Deviation for precision of duplicates.

$$\%RSD = \frac{141.42 * |A-B|}{(A+B)}$$

Percent Recovery of Trip spikes

$$\%Recovery = \frac{[Trip Spike Concentration]}{spike amount} * 100$$

Percent Recovery of QC Check Standards

$$\%Recovery = \frac{[Experimental]}{[Known]} * 100$$

The quality control data is kept in table format with new limits calculated annually. The field

staff is given the control limits for ready reference as samples are collected. The formulas for calculating control limits are based on the standard deviation of at least 30 measurements for each type of sample. The standard deviation is calculated according to the following formula.

$$s = \sqrt{\frac{\sum(X_i - \bar{X})^2}{(n-1)}}$$

where \bar{X} is the mean, X_i is an individual value, and n is the number of values.

$$\text{Mean} = \bar{X} = \frac{\sum X_i}{n}$$

The limits for field parameters are 2s for warning limits and 3s for control limits.

11.5.2 Lab Reportable Data

The precision and accuracy of each parameter are measured on a daily basis. The field spikes, blanks, and replicates are analyzed as routine samples.

Accuracy can be defined as the agreement between the actual obtained result and the expected result. Two QC check samples, QC1 and QC2, having a known or "true" value and standard additions also with an expected result, are used to test for the accuracy of a measurement system. Accuracy may be quantified by comparing results obtained for QC1 and QC2 to their true values and calculating a percent recovery using the following equation:

$$\text{Percent Recovery} = \%R = \frac{\text{experimental result}}{\text{"true" value}} * 100$$

The values obtained for the matrix spike are used to calculate percent recovery using the following equation:

$$\text{Percent Recovery} = \%R = \frac{[\text{matrix spike}] - [\text{sample}]}{[\text{spike}]} * 100$$

The percent recovery values may be used as an indication of bias. The control limits for accuracy are +/-2 standard deviations of the historical percent recoveries.

Precision can be defined as the agreement or closeness of two or more results and is an indication that the measurement system is operating consistently over a given time period. Because the SFWMD laboratory is highly automated and analyzes large numbers of samples in

a short period of time, one sample is chosen as the repeat (or replicate) sample for each parameter and is analyzed several times over the course of the run. Since each replicate sample is analyzed more than two times during the analytical run, it is appropriate to use the percent relative standard deviation, or the coefficient of variation, as the basis of acceptability of an analytical run. The control limits for precision are ± 2 standard deviations of the historical percent relative standard deviation.

The mean and standard deviation are calculated as shown in Section 11.4.1. The % relative standard deviation is calculated by the following equation:

$$\%RSD = (s/X) * 100$$

where X is the mean and s is the standard deviation.

On a daily basis, the results obtained for each of the quality control checks used are compared to the acceptable limits for precision and accuracy. New limits are calculated quarterly, with warning limits set at ± 1.5 standard deviations and control limits at ± 2.0 standard deviations of the most recent historical record for each type of QC check.

The procedures used to determine precision and accuracy are given in Table 11.4.

11.5 Method Detection and Practical Quantitation Limits

The method detection limits are determined by the procedures in 40 CFR Part 136, Appendix B. The method detection limits are updated annually and when necessary due to equipment or procedural changes.

The practical quantitation limit is 12 times the pooled standard deviations derived from the procedures to determine the method detection limit, and is updated annually.

Table 11.4 Procedures Used to Determine Precision and Accuracy

Method	Purpose	Concentration Level	Method References
Matrix Spike	Accuracy	Low Level Mid Level High Level	Nutrients, Trace Metals, Anions, Cations
Duplicates	Precision	Low Level Mid Level High Level	All parameters
Replicate for analytical run	Precision	Mid Level High Level	All parameters
QC Check Samples	Accuracy	Mid Level High Level	All parameters
QC Check Standards	Accuracy	Low Level Mid Level High Level	All parameters for which a sample is available
Method Reagent Blank	Accuracy	Low Level	Total Nutrients and Total Metals
Mid - Range Check Standard	Precision and Accuracy	Mid Level	Nutrients, Trace Metals, Anions and Cations

12.0 Data Reduction, Validation and Reporting

12.1 Data Reduction

12.1.1 Field Reportable Data

All the field measurement data are directly read from the instruments. These measurements include pH, specific conductance, dissolved oxygen, temperature, ORP, salinity, and depth. The data is automatically temperature compensated for pH, specific conductance, and dissolved oxygen. The cell constant for specific conductance is determined by the manufacturer. The field technician does not perform any calculations on field data.

The Technician responsible for data entry inputs all field data from the Chemistry Field Data Log into the computer.

12.1.2 Laboratory Reportable Data

The laboratory measurements which are read directly from the instruments, requiring no calculations, are pH, specific conductance, turbidity, color and fluoride. The pH and conductivity meters have automatic temperature compensation. The conductivity cells are purchased with cell constants provided by the manufacturer and are checked monthly. No calculations are performed using a single point reference or internal standard. Table 12.1 shows the formulas used to calculate specific parameters.

The analyst (Assistant Laboratory Technician, Laboratory Technician, Senior Laboratory Technician, or Chemist) is responsible for running the computer programs which provide the results in the appropriate concentrations and entering into the database.

The documentation for the results include the strip chart and chromatogram recordings and/or raw data files for the computerized calculations. The strip charts and chromatographs are identified with the date, the computer file name if applicable, and the initials of the analyst. The computer files are named using the instrument code, date, and sequential file number for the day.

12.2 Data Validation

12.2.1 Laboratory Data Integrity

The Analysts/Supervisors are responsible for checking the raw data entries and calculations for correctness. The Laboratory Quality Assurance Officer is responsible for checking sample preparation logs, and instrument/analytical logs for adherence to QC protocols and sample identification. The QA Officer is also responsible for checking the calibration integrity by comparing the present instrument responses to the historical values and the internal chain of

custody for the samples to ensure that only authorized persons have analyzed or handled a sample.

12.2.2 Field Data Integrity

The Quality Assurance Officer for each division is responsible for checking calibration integrity by checking the calibration logs and comparing present values to historical values and the sample custody integrity by checking the paperwork to ascertain that only trained personnel collected samples and that they were preserved and transported correctly. The Project Manager is responsible for checking raw data entries and calculations by reviewing the records for accuracy and use of proper formulas.

12.3 Specific Data Validation Procedures

12.3.1 Laboratory Data Validation

The analyst (Asst. Lab Technician, Lab Technician, Senior Lab Technician, or Chemist) is responsible for the first step in the validation process. It is his/her responsibility to follow the procedures correctly, perform the quality control checks, and report any discrepancies to his/her supervisor and/or the Laboratory Quality Assurance Officer. The Senior Chemists are responsible for investigating the discrepancy and determining the cause. The Senior Chemists are also responsible for the review of all data to identify obvious anomalies. The Laboratory Quality Assurance Officer is responsible for reviewing the quality control results for each run, and insuring that all QC criteria are met. The Quality Assurance Officer must update the acceptable quality control limits for all parameters quarterly.

The laboratory quality control checks are used to validate the laboratory results. Each of the QC check samples will be discussed.

Method Reagent Blank - if the result is greater than the detection limit, the run is stopped and the blank is prepared again. All samples since the last acceptable blank are reanalyzed.

Matrix Spike Sample - if the result is outside the current acceptable limits, the sample will be prepared again. If the value is outside the range again, the sample is analyzed by means of standard additions. If a matrix problem is verified, it is noted in the database that the sample exhibits matrix interference.

Quality Control Check Standards - if the result is outside the current acceptable limits, the run is stopped and the instrument recalibrated. If necessary new calibration standards are prepared and the instrument is checked for leaks, cracks in tubing, correct reaction temperature, correct wavelength or filter and correct calculation procedure in the computer.

Quality Control Check Samples - if the result for these blind samples is incorrect, the entire procedure is checked for errors. The analytical results are reported in the semiannual report.

Replicate Sample - used to assure that the same results are consistently obtained throughout the analytical run. At a minimum, it is run three times (beginning, middle, and end) during the analytical run if there are more than 20 samples in the run (two times for 20 or less samples). The %RSD should not vary more than the acceptable limit for that parameter.

Duplicate samples - at least 5% of the samples received by the laboratory are done in duplicate. The sample will not achieve "complete" status until the data is reviewed for precision.

Continuing calibration standard - used to confirm that the calibration curve remains constant throughout the analytical run. The value must remain within 5% of the initial value throughout the run.

12.3.2 Field Data Validation

The field sampling personnel (Scientific Technician, Senior Scientific Technician, Technician Supervisor, or Staff Environmental Scientist) is responsible for following the sampling procedures, reviewing the Chemistry Field Data Logs, and filling out all forms correctly and completely. The Technician is responsible for reviewing field data submitted to the laboratory for data entry for accuracy, initial data review following analysis, and review of the field quality control results for adherence to established standards. The Staff Environmentalist or Senior Environmentalist is responsible for reviewing all data for his/her project(s) to assure that the data quality objectives for the project(s) are being met.

The field quality control standards are used to validate the field collected data and the sample collection process.

Quality Control Check Standards - used to check the calibration of the instruments on a continuing basis. Results must be within the established acceptable limits.

Equipment Blanks - used to measure the effectiveness of sample equipment decontamination.

Field Blanks - used to measure the amount of environmental contamination in the samples.

Field Duplicates - used to measure the precision of the collection process.

Trip Spikes - used to measure sample shipping/handling contaminants and problems.

Trip Blanks - (for VOC only) used to measure shipping/handling contaminants and problems.

12.3.3 Project Data Validation

The Project Manager is responsible for the final review of data and its release to requestors. The Laboratory Quality Assurance Officer is responsible for review of the laboratory QC data and the

field quality control data. The Project Manager is responsible for review of all supporting documentation and the review of the data for anomalous results.

12.4 Data Reporting

Following the data validation, the results are entered in the LIMS data base by the Analysts (Analysts, Senior Lab Analyst, Staff Lab Analyst or Chemist). The Senior Chemists check entered data, by means of automated computer programs, to ensure detection of aberrant data, e.g., $\text{NO}_2^- > \text{NO}_3^-$, to avoid its inclusion into final reports. All reports are generated electronically from this data base by the Water Quality Monitoring Division Programmer Analyst. Printouts are routinely generated for all electronically transmitted data.

All requests for data reports must be made through Data Management. The Project Manager is responsible for reviewing the data before it is released to the requestor. An Example of the Final Report is given in Figure 12-1.

12.5 Data Storage

The records that will be retained are the strip charts, chromatograms, data files, Chemistry Field Data Logs, manual data entry records, daily QC reports, instrumentation logs and LIMS back up tapes.

The hard copy information is retained in the laboratory for no more than one year. After that it is sent to warehouse records storage and will be microfilmed if storage for longer than three years is needed.

The magnetic tapes used for weekly back-up of the laboratory information management system are maintained in two sets. One set is kept in an in-house vault for immediate access and one set is sent to secure storage outside the SFWMD for recovery in the event of a catastrophic event.

Archived records are indexed based on date for hard copy laboratory records, and by project for field notebooks and results in the water quality database. The Project Manager is responsible for the storage of all project for at least three years.

Table 12.1 Formulas Used for Calculations

Parameter	Formula
Alkaline Phosphatase	Computer, linear regression
Alkalinity	Computer, ml of titrant X N X 50,000/ml sample
Ammonia	Computer, linear regression
Ammonia, unionized	Ammonia value X chart value (pH & Temp)
Biochemical Oxygen Demand/ Carbonaceous Biochemical Oxygen Demand	$(DO_{(5 \text{ days})} - DO_{(\text{Initial})} - \text{Seed Correction})/P$ P = decimal volumetric fraction of sample used Seed Correction = DO loss in seed control X f f = ratio of seed in sample to seed in control
Calcium, dissolved	Computer, linear regression
Chemical Oxygen Demand	Computer, linear regression
Chloride	Computer, quadratic regression
Chlorine Residual, total	1 mL of FAS titrant = mg/L Cl
Inorganic Carbon, total and diss.	Computer, mean of 2 Repeats, multiple point calibration
Iron, total and diss.	Computer, linear regression
Nitrogen, organic	Computer, TKN - Ammonia
Nitrogen, total dissolved	Computer, TDKN + (Nitrate + Nitrite)
Nitrogen, total	Computer, TKN + (Nitrate + Nitrite)
Nitrogen, total Kjel.	Computer, linear regression
Magnesium, dissolved	Computer, linear regression
Nitrate + nitrite	Computer, linear regression
Nitrite	Computer, linear regression
Nitrate	Computer, (Nitrate+nitrite)-Nitrite
Total phosphorus	Computer, linear regression
Organic Carbon, total and diss.	Computer, Total - Inorganic, mean of 2 Repeats, multiple point calibration
Orthophosphate	Computer, linear regression

South Florida Water Management District, DHRs Certification E46077
 (Project KPA)

Lab Sample:	L777-1	L777-2	L777-3	L777-4	L777-5	L777-6	L777-7	L777-8	L777-9	L777-10	L777-11	L777-12
Field No.:	0946	0947	0948	0949	0950	0951	0952	0953	0954	0955	0956	0957
Station:	OKF-820	OKF-820	OKS-8291	OKS82501	OKF-84	OKF-84	OKF-84	OKF-84	OKF-84	OKF-84	OKF-84S1	OKS-94M1
FQC:	FB	OKF	OKS	OKS	SS	SS	SS	FB	FSB	RS	OKF	OKS
Date Collected:	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93	12-APR-93
Product:	ALCO3	CA	CL	COLOR	F	K	LCOMD	LPH	MG	NA	NH4	NH2
Units:	mg/l	mg/l	mg/l	units	mg/l	mg/l	mg/l	mg/l	mg/l	mg/l	mg/l	mg/l
Value:	<5	153.789	40.618	25.223	274.343	14.889	245.968	311.94	10.745	174.104	299.881	35.16
Value:	<5	289.24	11.844	13.292	14.889	14.755	14.755	132.271	<5	112.544	14.686	22.088
Value:	<.01	.366	<.004	<.004	.677	.358	.368	.647	<.01	.356	.358	2.505
Value:	<.004	<.004	.026	.005	<.004	.018	<.004	<.004	<.004	<.004	<.004	<.004
Value:	<.004	.007	.007	.052	.093	.09	.142	.142	.004	.119	.092	3.628
Value:	<2	175.793	<2	<2	<2	<2	58.725	<2	58.176	<2	<2	<2
Value:	<.05	<.05	.508	.545	1.29	1.205	1.372	<.05	.302	1.258	.329	.298
Value:	<.004	.028	.091	.126	.132	.251	<.004	.158	.123	.092	.092	.092
Value:	<.004	.028	.091	.126	.132	.251	<.004	.158	.123	.092	.092	.092
Value:	<.004	.028	.091	.126	.132	.251	<.004	.158	.123	.092	.092	.092

13.0 Corrective Action

Corrective action is required in those cases when the criteria levels for the quality control measures are not met. The specific corrective action for each type of quality control measure is given in Tables 13.1 and 13.2.

The analyst (Asst. Lab Technician, Lab Technician, Senior Lab Technician or Chemist) or field sample collection technician (Asst. Scientific Technician, Scientific Technician, Senior Scientific Technician) is responsible for assessing each QC measure and initiating corrective action according to Tables 13.1 and 13.2, respectively. The Supervising Professionals, Senior Chemists, Senior Supervising Technicians-Scientific, and the divisional Quality Assurance Officers are responsible for approving the corrective action taken or for initiating further steps to solve the problem.

External sources which may initiate corrective action are performance evaluation results, performance audits, system audits, split sample results, and laboratory/field comparison studies.

The problem and corrective action are documented in detail in one of the following: analysis logbooks, digestion logbooks, or instrument maintenance logs depending on the nature of the problem and how it was solved. The supervisor will report the problem to the appropriate Quality Assurance Officer who has the responsibility for determining if the solution is acceptable and if not, what further steps should be taken.

DEP recommended corrective action will be initiated as a result of systems or performance audits, split samples or data validation review.

Table 13.1 Corrective Actions for the Laboratory

QC Activity	Acceptance Criteria	Recommended Corrective Action
Initial Instrument Blank	Instrument response < MDL response	Prepare new blank, if same response determine cause of contamination: reagents, environment, equipment failure, etc.
Initial Calibration Standards	Coefficient of correlation >0.995 for AA/AE, >0.998 Visible spectrometer, fluorometer	Reanalyze standards, if same response, reoptimize instrument, if same response, prepare new standards
QC Check Standards	Value within +/- 2 standard deviations of the historical value	Reanalyze QC check standard, if same response, prepare new QC check standard, if same response, prepare new primary and calibration standards
Continuing Calibration Standards	+/- 5% of expected value	Reanalyze standard, if same response, recalibrate and reanalyze run from last continuing calibration standard
Replicate Sample	Coefficient of variation within established limit	Determine cause: baseline drift, carryover, etc. Reanalyze all samples if correctable cause not found
Duplicate Sample	Coefficient of variation within established limit	Reanalyze duplicates, reanalyze all samples between duplicates
Matrix Spikes	Value within +/- 2 standard deviations of the historical value	Remake spike and reanalyze, if acceptable reanalyze affected portions of run, if still not acceptable, spike a different sample. If second sample spike is acceptable, analyze first sample by standard addition. If second sample is not acceptable, spike all samples in that LIMS group in order to check for matrix interference.
Microbiology Control Blanks	Control Blank < 1	Sample resubmitted

Table 13.2 Corrective Actions for the Field

QC Activity	Acceptance Criteria	Recommended Corrective Action
Initial Calibration Standards	Value within +/- 5% of expected value	Reanalyze standards, if same response, optimize instrument, if same response, use new standards.
QC Check Standards	Value within +/- 5% of expected value	Reanalyze QC check standard, if same response, prepare new QC check standard, if same response, recalibrate

Table 13.3 Corrective Actions resulting from Field Quality Control Checks

QC Activity	Acceptance Criteria	Recommended Corrective Action
Equipment/Trip/Field Blank	Value < 2 MDL	Laboratory should reanalyze blanks: If same response, project manager should check recorded cleaning procedures and mark sample trip results for affected and related parameters questionable or invalidate data, as required. If second analysis acceptable, reanalyze affected samples in first run.
Trip Samples	Value within ± 2 standard deviations of the historical mean	Laboratory should reanalyze Trip Spikes: If Same Response, Sample Prep should be checked as well as other projects containing the same solution. If no other projects were affected, results for affected and related parameters should be invalidated or marked questionable, as required. If second analysis acceptable, reanalyze affected samples in first run.
Duplicate Samples	%RSD < 10	Laboratory should reanalyze duplicates: If same response, mark sample trip results for affected and related parameters questionable or invalidate data as required. If reanalysis shows Field Collection to be acceptable, reanalyze all samples analyzed with the Field samples the first time.

14.0 Performance and Systems Audits

Audits are an essential part of the quality assurance program for both laboratory and field operations. A systems audit is conducted to measure compliance with the comprehensive and project quality assurance plans. A performance audit is conducted to evaluate the quality of the data outputs with respect to mandatory limits or the laboratory's own performance standards.

14.1 Systems Audits

The systems audit is used to evaluate the entire measurement system, field and laboratory. It is a detailed review of each component of the sample collection process from equipment cleaning through submission of the samples to the laboratory and the laboratory process from sample log-in to archival of the results. Each element must be evaluated for conformance to appropriate methodology, approved procedures and the appropriate Quality Assurance Plans. A list of deficiencies must be made and addressed to correct, improve, or modify the system as necessary.

14.1.1 Internal Systems Audits

Systems audits are performed using the forms shown in Figure 14.1. The audit is conducted Semiannually by one of the following persons: the Division Quality Assurance Officers, Water Quality Monitoring Senior Technician, Senior Scientific Technician Supervisor or Environmental Scientist.

14.1.2 External Systems Audits

There are no regularly scheduled external systems audit. The SFWMD will submit to audits conducted by the DEP Quality Assurance Section.

14.2 Performance Audits

A performance audit is used to evaluate the routine quality control program of the laboratory.

14.2.1 Internal Performance Audits

Internal performance audits are conducted semiannually by the Laboratory Quality Assurance Officer. The audit may consist of any or all of the following: resubmission of previously analyzed samples under a different LIMS number, preparation of additional QC samples, samples split with another laboratory, and submission of spike samples, all of which are blind to the analysts.

The audit may be conducted at any time on a suspect parameter in addition to the semiannual audits for the entire laboratory. A performance audit on a parameter is mandatory when 10% of the analytical runs for that parameter fail one or more quality control criteria. Blind samples are immediately submitted.

The results of the performance audits are included in the quarterly quality assurance reports issued by the Laboratory Quality Assurance Officer.

14.2.2 External Performance Audits

The laboratory participates in three external performance audit programs. They are:

1. Florida Environmental Laboratory Certification Program administered by the Department of Health and Rehabilitative Services, semiannual
2. United States Geological Survey, Denver, annual
3. Florida Department of Environmental Protection Phosphorus Round Robin, semiannual

FIGURE 14.1 - SYSTEMS AUDIT CHECKLIST

Project Name: _____ Project No: _____

Project Manager: _____ Date: _____

Project Code (field and lab use): _____

Field Auditor: _____ Lab Auditor: _____

Signature: _____ Signature: _____

YES NO

Planning and Preparation:

- | | | |
|---|-------|-------|
| 1. Was QA Project Plan prepared for this project? | _____ | _____ |
| 2. Was briefing held with project participants, both field and lab? Date: _____ | _____ | _____ |
| 3. Were Additional instructions given to participants (ie., changes in project plan)? | _____ | _____ |
| 4. Was there a written list of sampling locations and descriptions? | _____ | _____ |
| 5. Was there a map of sampling locations available to field personnel? | _____ | _____ |
| 6. Was sampling scheduled with field technician supervisor in advance (minimum one week)? | _____ | _____ |
| 7. Were analyses scheduled with laboratory in advance (minimum one week)? | _____ | _____ |

Comments : _____

<u>General Sampling Procedures:</u>	YES	NO
1. Were sampling locations properly selected?	_____	_____
2. For what analyses were the samples collected? _____	_____	_____
3. Was sampling equipment protected from possible contamination prior to sample collection?	_____	_____
4. Were clean disposable latex or vinyl gloves worn during sampling?	_____	_____
5. Were gloves changed for each sample station?	_____	_____
6. If equipment was cleaned in the field, were proper procedures used?	_____	_____
7. What field instruments were used? _____ _____ _____		
8. Were calibration procedures documented in the field notes?	_____	_____
9. Was calibration check conducted within 24 hours or at the end of the sampling day for Hydrolabs and every 4 hours and at the end of the day for other field instruments?		_____
10. Were samples chemically preserved at time of collection?	_____	_____
11. Was the preservative amount recorded in the field notes?	_____	_____
12. Was the pH of preserved samples checked to insure proper preservation?	_____	_____
13. Were the proper sample bottles used?	_____	_____
14. Were samples iced at the time of collection?	_____	_____
15. Were sample bottles rinsed with sample before filling?	_____	_____
16. Were field conditions recorded in the field notes?	_____	_____

General Sampling Procedures: (continued) YES NO

- 17. Was filtering equipment pre-rinsed with sample? _____
- 18. Were samples processed on a clean tray/surface? _____
- 19. Was waste material containerized and maintained separate from samples and equipment? _____

Comments: _____

Surface Water Sampling: YES NO

- 1. What procedures were used to collect the surface water samples? _____

- 2. Was a Niskin bottle used for sample collection? _____
- 3. Was a pump used for sample collection? _____
- 4. If pump was used, what type? _____
- 5. Was the pump properly cleaned before/between sites? _____
- 6. Was the pump tubing properly cleaned before/between sites? _____
- 7. Was sample collected in polyethylene bucket? _____
- 8. Did samplers wade in stream during sample collection? _____
- 9. Were samples collected upstream from the sampler or vehicle? _____

Surface Water Sampling (continued)

YES NO

10. Were autosamplers used? _____
11. Were samples composited? _____
12. How were samples composited? _____

13. Were autosampler bottle(s) capped and shaken well before compositing or processing? _____
14. Was tubing changed within the last three months? _____
15. Were Preventative maintenance protocols followed and documented? _____

Comments: _____

Well Sampling:

YES NO

1. Was depth of well determined? _____
2. Was depth to water determined? _____
3. Was measuring tape properly decontaminated between wells? _____
4. How was the volume of water originally present in each well determined? _____
5. Was the volume determined correctly? _____

6. How was completeness of purging determined:

Volume measure _____
Time _____
Flow rate _____
Scond./pH/T _____

YES NO

7. Was a sufficient volume purged? _____

8. Was the well over-purged? _____

9. Was the disposal of purge water handled properly? _____

10. Was a dedicated (in-place) pump used? _____

If no: Describe method of purging: _____

11. How were the samples collected?

Bailer? _____

Pump? _____

Other? _____

12. Construction material of bailer? _____

13. If a pump was used, describe how it was cleaned
before and/or between wells: _____

14. Were the samples properly transferred from bailer
to sample bottles (i.e., was the purgeable sample
agitated, etc.)? _____

15. Was the rope or line allowed to touch the ground? _____

16. Was a teflon coated stainless steel cord used? _____

Comments: _____

Sediment Sampling:

YES NO

1. What procedures were used to collect the samples?

2. Were the samples well mixed prior to placing the sample in the sample container?

3. Were the samples composited?

Comments:

Atmospheric Deposition:

1. Were buckets maintained capped before deployment? _____

2. Were previously deployed buckets sealed immediately and lids taped? _____

3. Were field conditions documented in the field log? _____

4. Were sample buckets secured and kept clean in vehicle? _____

5. Was preventative maintenance and cleaning performed as specified in the SOP/CompQAP? _____

6. Was balance calibration performed? _____

7. Was tare weight legible? _____

8. Was tare weight determined within the last quarter? _____

9. Was Bucket Equipment Blank performed? _____

Comments:

Other Sampling:

1. What other types of samples were collected during this investigation? _____

2. What procedures were used for the collection of these samples? _____

Comments: _____

Field Quality Control:

YES NO

1. Were QC samples specified in the QA Project Plan? _____ _____
2. Were the QC samples collected in accordance with the QA project plan? _____ _____
3. Did sampling personnel utilize any trip blanks? _____ _____
4. Did sampling personnel utilize any field blanks? _____ _____
5. Were any equipment blanks collected? _____ _____
6. Were any duplicate samples collected? _____ _____
7. Were any trip spikes utilized? _____ _____
8. Check method used to collect split sample
 _____ Filled one large container and then transferred portions
 _____ Sequentially filled bottles
9. Were chain of custody records completed for all samples? _____ _____

<u>Field Quality Control (continued):</u>	YES	NO
10. Were all samples identified with appropriate tags?	_____	_____
11. Were sample I.D. tags filled out properly?	_____	_____
12. Did information on sample I.D. tags and Chemistry Field Data Log match?	_____	_____
13. Were samples kept in a secure place after collection?	_____	_____
14. Was Chemistry field Data Log signed by sampling personnel?	_____	_____
15. Were amendments to the project plan documented (on the project plan itself, in a project logbook, elsewhere)?	_____	_____

Comments: _____

<u>General Laboratory Procedures:</u>	YES	NO
1. Have unique sequential laboratory numbers been assigned to each sample?	_____	_____
2. Has the data from the Chemistry Field Data Log been input to the computer directly?	_____	_____
3. Have samples been stored in an appropriate secure area?	_____	_____
4. Has sample custody been maintained by the laboratory?	_____	_____
5. Has the proper bar code label been attached to each sample I.D. tag?	_____	_____
6. Were the samples aliquoted properly?	_____	_____

Comments: _____

<u>Analytical Methods:</u>	YES	NO
1. Have approved analytical methods or procedures been followed?	_____	_____
2. Does the project plan include copies of any non-standard methods without appropriate quality assurance results for validation of the method?	_____	_____
3. Does use of the analytical methods specified result in data of adequate detection limit, accuracy, and precision to meet the requirements of the project?	_____	_____

Comments: _____

<u>Laboratory Quality Control:</u>	YES	NO
1. Have approved sample holding times been observed?	_____	_____
2. Have replicate analyses been performed on at least one sample?	_____	_____
3. Have spike analysis been performed on at least one sample?	_____	_____
4. Have the quality control reporting forms been properly filled out?	_____	_____
5. Are current instrument calibration curves used for all methods?	_____	_____
6. Did the spiking procedures follow acceptable protocols for quantity and concentration?	_____	_____

Laboratory Quality Control (continued):

YES NO

- | | | |
|--|-----|-----|
| 7. Are quality control charts used to track QC precision and accuracy? | ___ | ___ |
| 8. Are QC charts kept up to date? | ___ | ___ |
| 9. Is the precision of the data presented within acceptable limits? | ___ | ___ |
| 10. Is the accuracy of the data presented within acceptable limits? | ___ | ___ |
| 11. Are recent (one year or less) performance audit results available? | ___ | ___ |
| 12. Has the laboratory followed the preventative maintenance procedures outlined in the QA plan? | ___ | ___ |
| 13. Is the completeness of the data acceptable? | ___ | ___ |

Comments: _____

Data Validation and Reporting:

YES NO

- | | | |
|--|-----|-----|
| 1. Were all the steps in the data validation procedure outlined in the QA plan followed? | ___ | ___ |
| 2. Was the data reported in the proper format with the proper units? | ___ | ___ |
| 3. Was the laboratory I.D. number included on each page of the data? | ___ | ___ |

Comments: _____

15.0 Quality Assurance Reports

The Divisional Quality Assurance Officers are responsible for preparing a quarterly internal quality assurance report. These reports are for internal use and are not submitted to DEP. These reports to the Division Director include the following:

1. An assessment of data accuracy, precision, and method detection limits
2. Results of performance and systems audits
3. Significant quality assurance/quality control problems and the recommended solution
4. Outcome of any corrective action.

Quality assurance reports to DEP for SWIM projects are made on a variable schedule based on the frequency of sampling since these are on-going monitoring projects. For weekly, biweekly, or monthly sampling schedules, QA reports are submitted semiannually. For bimonthly, quarterly or semiannual sampling schedules, the QA reports are submitted annually. The reports to DEP are written by the Divisional Quality Assurance Officers, Laboratory Quality Assurance Officer and the Project Managers. The Project Manager is responsible for submitting the report to DEP.

If no project audits are performed and no significant quality assurance/quality control problems occur for a specific project, a letter stating these facts will be sent to DEP in lieu of the quality assurance report.

The quality assurance reports must include the following for performance audits:

1. Date of the audit
2. System tested
3. Persons performing/administering the audit
4. Parameters analyzed
5. Reported results
6. True values of the samples (if applicable)
7. If any deficiencies or failures occurred, a summary of the problem and the corrective action taken
8. Copies of documentation.





The quality assurance reports must include the following for systems audits:

1. Date of the audit
2. System tested
3. Who performed/administered the audit
4. Parameters analyzed
5. Results of tests
6. Parameters for which results were unacceptable
7. Explanation of the unacceptable results including probable reasons and the

- corrective action taken
- 8. Copies of documentation.

For significant quality assurance/quality control problems, the following information must be included in the report:

1. Identify the problem and the date it was found
2. Identify the individual who reported the problem
3. Identify the source of the problem
4. Discuss the solutions and corrective actions taken to eliminate the problem.

METHOD # 1100.1	METHOD NAME COLOR	REVISION 1.1	REFERENCE SM 2120B	STATUS EFFECTIVE
PAGE 1 OF 5				DATE 01/18/94
APPROVALS/DATE				
DIVISION DIRECTOR 	LAB SUPERVISOR  07/22/94	QUALITY ASSURANCE  8-23-94	METHODS  8/22/94	

SECTION 1.0 - TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
2.0 - METHOD DESCRIPTION/HISTORY	1
3.0 - SAFETY PRACTICES	1
4.0 - LIST OF EQUIPMENT/INSTRUMENTATION	2
5.0 - REAGENTS	2
6.0 - STANDARDS	2
7.0 - QUALITY CONTROL	2
8.0 - STEP-BY-STEP PROCEDURE	3
9.0 - DATA HANDLING	4
10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING	5
11.0 - REFERENCES	5

SECTION 2.0 - METHOD DESCRIPTION/HISTORY

2.1 Color in water may result from the presence of metallic ions (iron and manganese), humus and peat materials, plankton, weeds, and industrial waste.

2.2 The term "color" here is termed to mean true color, that is, the color of the water from which turbidity has been removed. In our case the turbidity is removed by filtration (0.45 um membrane filter). It should be noted that the filtration procedure may remove some of the color from the sample.

2.3 The color of the sample is determined by spectroscopic comparison to platinum-cobalt color solutions at 465 nm in a 1 cm quartz flowcell. One unit of color corresponds to 1 mg/L of platinum in the form of the chloroplatinate ion. The method is applicable to waters which are colored naturally due to vegetative decay and the presence of metallic ions.

2.4 This method version 1.1 was adopted on January 18, 1994 and is an internally developed (SFWMD) method and is based on Standard Methods SM 2120B. The previously used method was the same in principle, but this method reflects a new format and includes changes due to the implementation of a new LIMS system in October of 1992.

SECTION 3.0 - SAFETY PRACTICES

3.1 Wear safety glasses and a full-length, long-sleeved laboratory coat.

3.2 Latex or polyethylene gloves (non-powdered) may be worn when handling the samples.

METHOD # 1100.1	METHOD NAME COLOR	REVISION 1.1	REFERENCE SM 2120B	STATUS EFFECTIVE
PAGE 2 OF 5				DATE 01/18/94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 3/4/94	LAB SUPERVISOR <i>[Signature]</i>	QUALITY ASSURANCE 1-26-94 <i>[Signature]</i>	METHODS <i>[Signature]</i> 1/18/94	

SECTION 3.0 - SAFETY PRACTICES (CONT'D)

3.3 All personnel conducting this method should be familiar with the SFWMD Chemical Hygiene Plan and should have reviewed any pertinent Material Safety Data Sheets. Note: pay special attention to the MSDS for Platinum Cobalt Solution, as it is a suspected carcinogen (Note: Prolonged exposure to Platinum-Cobalt solutions causes degradation of tooth enamel).

3.4 No hazardous wastes are generated by this procedure; liquid waste may be disposed of in the sink.

SECTION 4.0 - LIST OF EQUIPMENT/INSTRUMENTATION

- 4.1 Visible Spectrometer, Model Spectronics 501
- 4.2 Sample introduction system (flowcell, transfer tubing, and pump)
- 4.3 Beaker 1 L, polypropylene, for waste collection (Fisher #02-586-6H)

SECTION 5.0 - REAGENTS

5.1 Cell cleaning solution (5% Potassium Hydroxide in Methanol).

SECTION 6.0 - STANDARDS

- 6.1 Platinum Cobalt Color Stock Solution/Standard 500 c.u. (500 mg/L) - (Fisher Scientific #SO-P-120).
- 6.2 Platinum Cobalt Color Standard 300 c.u. (300 mg/L) - measure 30 mL of stock solution with a class A volumetric pipet and dilute to 50 mL in a class A volumetric flask.
- 6.3 Platinum Cobalt Color Standard 100 c.u. (100 mg/L) - measure 10 mL of stock solution with a class A volumetric pipet and dilute to 50 ml in a class A volumetric flask.
- 6.4 Platinum Cobalt Color Standard 50 c.u. (50 mg/L) - measure 5 mL of stock solution with a class A volumetric pipet and dilute to 50 mL in a class A volumetric flask.

SECTION 7.0 - QUALITY CONTROL

7.1 The Platinum Cobalt Color standards are checked by the analyst before conducting the analyses. The results are recorded on the physical parameters log.

7.2 QC1 and QC2 will be prepared as required by the QA unit. QC1 and QC2 are analyzed at beginning of each set of analyses, by the analyst, and are repeated after each 20 samples.

METHOD # 1100.1	METHOD NAME COLOR	REVISION 1.1	REFERENCE SM 2120B	STATUS EFFECTIVE
PAGE 3 OF 5				DATE 01/18/94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 3/6/94	LAB SUPERVISOR <i>[Signature]</i>	QUALITY ASSURANCE 1-26-94 <i>[Signature]</i>	METHODS <i>[Signature]</i> 1/18/94	

SECTION 7.0 - QUALITY CONTROL (CONT'D)

7.3 A repeat analysis should be conducted for every 20 samples analyzed. These results are recorded on the QC sheet and physical parameters log. The mean and coefficient of variation of the replicate set is determined and noted on the QC sheet and physical parameters log.

7.4 All quality control data must be within the current established limits before entering sample data into the LIMS system.

7.5 Samples should be shaken gently prior to analyses; excessive shaking will entrain air and result in erroneous readings.

SECTION 8.0 - STEP-BY-STEP PROCEDURE

8.1 Create a workgroup for the analyses to be conducted. Refer to the ACS LIMS users Guide for instructions on "Creation of a Workgroup". Remove the samples to be tested from the refrigerator and allow them to warm to room temperature (this may require several hours).

8.2 Turn on the spectrometer and allow it to stabilize for at least 30 minutes. Select 465nm wavelength if required (press "4", "6", "5", <SECOND FUNCTION>, <YES>).

8.3 Log on to the LIMS terminal adjacent to the spectrometer and type "color" at the UNIX prompt. Enter the workgroup number and the system will prompt you for the result of the first sample.

8.4 Add 4 drops of cell cleaning solution into 20 mL of fresh D.I. water in a glass graduated cylinder. This mixture should be used to clean-up the cell before conducting any further work.

8.5 Insert the inlet tubing into the cleaning solution prepared in 8.3 and turn on the pump and allow all of the solution to pass into the flow cell.

8.6 Enter the Absorbance mode on the spectrometer by pressing the <%T/A/C> button on the spectrometer.

8.7 Rinse a beaker with D.I. water and fill it with fresh D.I. water. Place the inlet tube into the beaker of distilled water and turn the pump on. Allow the water to pass into the flow cell until the reading on the spectrometer stabilizes. When stable zero the spectrometer by pressing the <SECOND FUNCTION> button and the <100 %T/ Zero A> button.

8.8 Next pump the 500 c.u. color standard into the flow cell until the reading stabilizes, check the absorbance. If the reading is in the acceptable range (0.130 - 0.145) continue to the next step. If the reading is not acceptable repeat steps 8.4 to 8.6. If an acceptable reading cannot be obtained, contact your shift supervisor or the QA officer for assistance.

METHOD # 1100.1	METHOD NAME COLOR	REVISION 1.1	REFERENCE SM 2120B	STATUS EFFECTIVE
PAGE 4 OF 5				DATE 01/18/94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 3/4/94	LAB SUPERVISOR <i>[Signature]</i>	QUALITY ASSURANCE <i>[Signature]</i> 1-26-94	METHODS <i>[Signature]</i> 1/18/94	

SECTION 8.0 - STEP-BY-STEP PROCEDURE (CONT'D)

- 8.9 Enter the concentration mode by pressing the <%T/A/C> button. Remove the inlet tube from the water and clean the inlet tube with a Kimwipe. Place the inlet tube into the beaker of D.I. water and zero the instrument as in step 8.7. Place the inlet tubing into the beaker of 500 c.u. standard and allow the reading to stabilize. Calibrate the instrument by pressing "5", "0", "0", and <ENTER>.
- 8.10 Clean the inlet tube with a Kimwipe prior to its insertion into the next solution. Pump the 300 c.u. standard into the flow cell allow the reading to stabilize and record the reading on the physical parameters log sheet. The reading should be 300 +/- 5 c.u.
- 8.11 Repeat step 8.10 for the remaining standards. The 100 c.u. and 50 c.u. standards should read within +/- 3 c.u.
- 8.12 If the standards are within limits repeat step 8.8 for the QC1 and QC2 solutions; they should read within the current limits for the QC solution.
- 8.13 The inlet tube should be cleaned with a Kimwipe after removing it from any standard, QC, or sample solution to prevent cross contamination.
- 8.14 If the QC samples are within limits proceed to analyze the samples following step 8.10. When the reading has stabilized press the "send" button on the spectrometer to send the result to the LIMS system. Proceed to analyze each sample as prompted by the LIMS.

SECTION 9.0 - DATA HANDLING

- 9.1 The detailed data handling procedures are covered in section 8.0, however 9.2 through 9.8 summarize the data handling process.
- 9.2 Sign onto the workstation
- 9.3 Go into the ACS LIMS via the command `sl` and create a workgroup for the test COLOR and exit the SEEDPAK1 Main Menu (see Creating a Workgroup and Printing a Workgroup in the ACS LIMS Users Guide).
- 9.4 At the UNDX prompt type "color" .
- 9.5 The system will prompt you for the workgroup number you just created for the samples. Enter the workgroup number (i.e. WG19) and hit "Enter".

METHOD # 1100.1	METHOD NAME COLOR	REVISION 1.1	REFERENCE SM 2120B	STATUS EFFECTIVE
PAGE 5 OF 5				DATE 01/18/94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 3/4/94	LAB SUPERVISOR <i>[Signature]</i>	QUALITY ASSURANCE <i>[Signature]</i> 1-26-94	METHODS <i>[Signature]</i> 1/18/94	

SECTION 9.0 - DATA HANDLING (CONT'D)

9.6 The system will then prompt you, one by one, for each sample in the workgroup. Hit the send button on the spectrometer and the system will send the result and advance to the next sample. Hitting "Enter" with no test result will cause the program to skip that sample and go on to the next sample. To get back to a sample that you skipped, hit "Enter" until the program asks for the workgroup number again. Once you enter the workgroup number again, you will be prompted with the samples that are incomplete.

9.7 To end the program, hit "Enter" until you see "Done..." and are back at the UNIX prompt.

9.8 Pick up the data entry report for the samples you have just entered from the system printer.

SECTION 10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING

10.1 The spectrometer should be wiped clean after each use.

10.2 The flow cell should be cleaned after each use with the cell cleaning solution by pumping it into the cell and allowing it to stand for a few minutes. Pump D.I. water through the cell after the cleaning solution for a few minutes to completely remove it from the flow cell.

10.3 Leave the cell filled with D.I. water when not in use.

10.4 Release the clamps from the pump tubes when not in use.

10.5 Consult your supervisor before making any major changes, adjustments, and/or repairs to the instrumentation.

SECTION 11.0 - REFERENCES

11.1 EPA Methods for Analysis of Water and Wastes, EPA-600/4-79-020, March 1979.

11.2 Standard Methods for the Examination of Water and Wastewater, 18th Edition, 1992.

11.3 ACS LIMS Users Guide, version 1.0, 1992.

11.4 SFWMD Comprehensive Quality Assurance Manual, current version.

11.5 Spectronics 501 Operators Manual

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 1 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i>	LAB SUPERVISOR <i>[Signature]</i> 10-01-94	QUALITY ASSURANCE 11-1-94 <i>[Signature]</i>	METHODS <i>[Signature]</i>	

SECTION 1.0 - TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
2.0 - METHOD DESCRIPTION/HISTORY	1
3.0 - SAFETY PRACTICES	2
4.0 - LIST OF EQUIPMENT/INSTRUMENTATION	2
5.0 - REAGENTS	3
6.0 - STANDARDS	4
7.0 - QUALITY CONTROL	4-5
8.0 - STEP-BY-STEP PROCEDURE	5-7
9.0 - DATA HANDLING	7-8
10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING	8
11.0 - REFERENCES	9

SECTION 2.0 - METHOD DESCRIPTION/HISTORY

2.1 The reaction of silicate with molybdate forms B-molybdosilicic acid at a pH of 1.0-1.8. The B-molybdosilicic acid is reduced by tin (II) to form molybdenum blue, which is measured at 820 nm.

2.2 Interference from orthophosphate and tannin is eliminated by the use of tartaric acid. Color absorbing at the analytical wavelength will interfere. Samples for analysis are filtered in the field through a 0.45 µm filter.

2.3 The analyses are conducted in a highly automated instrument called a Rapid Flow Analyzer (RFA). This instrument is equipped with an autosampler for sample introduction, a peristaltic pump, a mixing manifold, and a photometer for colorimetric measurement. The analog output of the photometer is relayed to a personal computer equipped with ALPKEM data collection software. The same software is used to calculate sample concentrations.

2.4 This modified method version is based on APHA Standard Method 4500SiD, and is a revision of SOP dated November, 1990. This method reflects a new format and includes changes due to change in instrumentation and the implementation of a new LIMS system in October of 1992.

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 2 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11-1-94	LAB SUPERVISOR <i>[Signature]</i> 12-26-94	QUALITY ASSURANCE <i>[Signature]</i> 11-1-94	METHODS <i>[Signature]</i>	

SECTION 3.0 - SAFETY PRACTICES

- 3.1 Wear safety glasses and a full-length, long-sleeved laboratory coat.
- 3.2 Latex or polyethylene gloves (non-powdered) may be worn when handling the samples.
- 3.3 All personnel conducting this method should be familiar with the SFWMD Chemical Hygiene Plan and should have reviewed any pertinent Material Safety Data Sheets.
- 3.4 The disposal of samples can be done in the sink, flushing with ample amounts of tap water.
- 3.5 Preparation of reagents containing hydrochloric acid, chloroform and stannous chloride solutions should be conducted in a fume hood. The reagents should be prepared by slow addition of concentrated hydrochloric acid to D.I. water. Use an acid resistant bottle carrier when carrying glass containers of concentrated hydrochloric acid and chloroform.
- 3.6 Before starting any run, all lines connecting the instrument to the reagents should be checked and tightened if necessary. In case of a leak onto an electrical system, the power should be disconnected before conducting any repairs.
- 3.7 The electrical power should be disconnected before conducting any repairs inside the instrument on controllers, electrical wiring or any other components near sources of electricity.
- 3.8 In case of spills of concentrated hydrochloric acid the spill should be first treated with an appropriate spill kit and the contaminated absorbent should be collected and placed into adequate storage containers for disposal.

SECTION 4.0 - LIST OF EQUIPMENT/INSTRUMENTATION

- 4.1 ALPKEM™ Rapid Flow Analyzer (RFA), Model 300 with XYZ autosampler
- 4.2 Personal Computer equipped with A/D converter, printer and ALPKEM data reduction software.
- 4.3 Class A volumetric glassware (pipets and volumetric flasks)
- 4.4 Clean Nalgene plastic containers.

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 3 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 11-26-94	QUALITY ASSURANCE <i>[Signature]</i> 11-7-94	METHODS <i>[Signature]</i>	

SECTION 5.0 - REAGENTS

5.1 Sodium lauryl sulfate solution - Dissolve 5.0 grams of dodecyl sodium sulfate in 90 mL of deionized water contained in a 250 mL Erlenmeyer flask. It may be necessary to warm the solution to obtain complete dissolution. Transfer the solution to a 100 mL volumetric flask and dilute to mark with deionized water. Transfer to a small plastic dropping bottle for daily use. (Safety PP = A)

5.2 Tartaric acid 10% w/v - Dissolve 100 grams of tartaric acid in approximately 800 mL of deionized water. Dilute the solution to 1 liter. Transfer the solution to a liter plastic container. Add 2 drops of chloroform and shake well. Store the reagent at 2-6 °C. (Safety PP = A)

5.3 Hydrochloric acid 1.2 N - In a fume hood, cautiously add (with stirring) 100 mL of concentrated (12N) hydrochloric acid to approximately 800 mL of deionized water contained in a 1 liter volumetric flask. When the solution has returned to room temperature, dilute to 1 liter and mix well. Store the reagent in a 1 liter plastic container. (Safety PP = F, C)

5.4 Stock stannous chloride - In a fume hood, cautiously add (with stirring) 10 mL of concentrated (12N) hydrochloric acid to 10 mL of deionized water in a 50 ml pyrex beaker. Dissolve 10 grams of stannous chloride in the acidic solution. Heating may be required to obtain complete dissolution. Store the stock solution in a tightly closed plastic container and refrigerate at 2-6 °C. (Safety PP = F, C)

5.5 Working stannous chloride reagent - In a 60 mL plastic container, mix together 50 mL of 1.2 N hydrochloric acid and 0.5 mL stock stannous chloride. This reagent should be prepared fresh daily. (Safety PP = A)

5.6 Ammonium molybdate reagent - Dissolve 1.080 grams of ammonium molybdate in approximately 80 mL of deionized water. Add 0.3 mL (9 drops) of 50% sulfuric acid. Add 1mL of 5% sodium lauryl sulfate, dilute to 100 mL with deionized water, and mix well. Transfer the solution to a clean plastic 175 mL container. This reagent should be prepared fresh daily. (Safety PP = A)

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 4 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 11-25-94	QUALITY ASSURANCE <i>[Signature]</i> 11-1-94	METHODS <i>[Signature]</i>	

SECTION 6.0 - STANDARDS

NOTE: Immediately after mixing, standards should be poured into labelled, clean 175 mL plastic containers.

6.1 Silica stock standard - NBS/NIST standard solution or traceable stock, 10 mg/mL Si or 21.3 mg/mL SiO₂.

6.2 Solution A (1065 mg/L or 1.065 mg/mL) - In a 100 mL class A volumetric flask, pipet 5 mL of stock solution and dilute to the mark with deionized water. Cap the flask and mix well. Transfer the solution to a clean 175 mL plastic container.

6.3 Standard 1 (21.3 mg/L) - In a 100 mL class A volumetric flask, pipet 2 mL of solution A and dilute to the mark with deionized water. Cap, mix well and transfer the solution to a clean 175 mL plastic container.

6.4 Standard 2 (10.65 mg/L) - In a 100 mL class A volumetric flask, pipet 1 mL of solution A and dilute to the mark with deionized water. Cap, mix well and transfer the solution to a clean 175 mL plastic container.

6.5 Standard 3 (5.33 mg/L) - In a 200 mL class A volumetric flask, pipet 1 mL of solution A and dilute to the mark with deionized water. Cap and mix well. After making standard 4, transfer 100 mL to a clean 175 mL plastic container.

6.6 Standard 4 (2.67 mg/L) - In a 100 mL class A volumetric flask, pipet 50 mL of standard 3 and dilute to the mark with deionized water. Cap, mix well and transfer to a clean 175 mL plastic bottle.

6.7 Standard 5 (0.0 mg/L) - Deionized water only.

SECTION 7.0 - QUALITY CONTROL

7.1 QC1 and QC2 are prepared fresh monthly or as needed by the QA unit. QC1 and QC2 are analyzed at beginning of each set of analyses. QC2 is repeated at the end of the analytical run.

7.2 Spikes are prepared from samples selected at random (1 for every 20 samples analyzed), and are made by adding 0.1 mL of the solution A (with a 1.0 mL Tensette) to a 10 mL volumetric flask and diluting to the mark with the sample. Calculate the STAD recovery as:

$$\text{STAD Recovery} = \frac{\text{STAD Conc.} - [\text{Sample Conc.} \times 0.99]}{100} \times 100$$

Enter the %STAD Recovery on the QC and RFA log forms.

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 5 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 11-26-94	QUALITY ASSURANCE <i>[Signature]</i> 11-1-94	METHODS <i>[Signature]</i>	

SECTION 7.0 - QUALITY CONTROL (CON'T)

7.3 A repeat analysis (of a matrix sample chosen at random) must be run for every 20 samples analyzed. Enter each repeat values, the calculated Mean, Standard Deviation, and Coefficient of Variation on the QC sheet.

7.4 All quality control data must be within the current established limits, and the run must be checked and initialed by the supervisor or the QA officer before sending the sample data into the LIMS system.

SECTION 8.0 - STEP-BY-STEP PROCEDURE

- 8.1 Sign on to LIMS and create a workgroup (refer to section 9.0 - Data Handling).
- 8.2 Turn 314 regulated power ON.
- 8.3 Turn on light source.
- 8.4 Place all reagent lines in DI H₂O containing 5 drops of sodium lauryl sulfate.
- 8.5 Latch platens and turn on 302 pump module.
- 8.6 After flow has stabilized, verify smooth and consistent bubble pattern throughout the manifold.
- 8.7 Set photometer parameters as defined by flow diagram.
- 8.8 On photometer, set center knob to "sample" position.
- 8.9 Slowly turn the sample fine adjust knob to set the LCD display to 5.00 volts.
- 8.10 On photometer, set center knob to "reference" position.
- 8.11 Slowly turn the reference fine-adjust knob to set the LCD display to 5.00 volts.
- 8.12 Place center knob in "absorbance" position. Using the reference fine adjust knob, set the LCD display to 0.20 (± 0.01).

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 6 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 11-07-94	QUALITY ASSURANCE <i>[Signature]</i> 11-1-94	METHODS <i>[Signature]</i>	

SECTION 8.0 - STEP-BY-STEP PROCEDURE (CON'T)

8.13 Place reagent lines in their respective containers. Allow 2 to 3 minutes before placing stannous chloride line in the reagent bottle.

8.14 After 15 minutes, reagent flow should be stabilized. Repeat steps 8.8 through 8.12.

8.15 While reagents are stabilizing, prepare sample table.

8.16 Turn on main power switch of computer. A blue menu screen appears after a screen of credits.

8.17 Build a sample table by selecting F4 and entering standard, sample, and quality control identifiers with dilution factors (if applicable).

8.18 When sample table is finished, press ALT S to save the table. Sample tables are named using A 3 letter parameter code, a 3 digit Julian date, and a letter corresponding to the order of analysis; example: SIO136A - 1st run; SIO136B - 2nd run.

8.19 Press ALT P to print a hard copy of the table file.

8.20 Press ESC to return to main menu.

8.21 Press F5 to get to data collection mode. Use the "space" bar to turn on the channel to be viewed. Press F3 to monitor baseline.

8.22 While baseline is being monitored, begin pouring the analytical run according to the sample table.

8.23 Manually activate the sampler and set the stop count to the cup number corresponding to the end of the first set of standards.

8.24 Observe that the standards appear linear and that peaks do not have spikes or any unusual shape to them. Press reset on sampler.

8.25 If standards appear normal, press escape on computer. To begin data collection, press ALT 1 (simultaneously).

8.26 Pour all samples, including any dilutions.

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 7 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i>	QUALITY ASSURANCE 11-1-94 <i>[Signature]</i>	METHODS <i>[Signature]</i>	

SECTION 8.0 - STEP-BY-STEP PROCEDURE

8.27 When the run is complete, the cup position on the computer screen will read PS (Pause). After the last peak is shown on the screen, the baseline is collected until the timer in the right corner of the screen stops. At this point, press escape to return to the blue data collection screen. Press ALT 1 S to stop the analytical run.

8.28 To calculate data, refer to sections 9.3 - 9.13.

SECTION 9.0 - DATA HANDLING

9.1 Sign onto the workstation.

9.2 Go into the ACS LIMS via the command `s1` and create a workgroup for SIO2 and exit the SEEDPAK1 Main Menu (see Creating a Workgroup and Printing a Workgroup in the ACS LIMS Users Guide).

9.3 At the main menu on the computer, press **F8** (calculation sub-menu).

9.4 Make sure the sample table to be calculated is the current one loaded. Otherwise, go to **F4** and recall the right sample table file (ALT 1, L, Filename)

9.5 Press **F8** again; select: from raw data. Press "ENTER".

9.6 At filename prompt, make sure extension corresponds to the channel name to be calculated; example:
*.D11 = channel 1; *.D12 = channel 2; *.D13 = channel 3. (* = filename)

9.7 Review peak finding parameters and press "enter" at execute Y prompt.

9.8 Observe that peak markers appear on each peak. Press **F9** to continue.

9.9 After observing standard curve, press **F9** until calculation values appear on screen.

9.10 Check each quality control value to be sure it falls within acceptable limits.

9.11 If run is acceptable, press **F2** (LPT1), then press **Y** to print a hard copy of the file.

9.12 When hard copy is finished printing, press **F2** again. Backspace over LPT1. Type in filename followed by .IN extension. Example: SIO136A.IN. Press **N** at Y/N prompt.

9.13 Press **F9** to save and exit.





METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 8 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i>	LAB SUPERVISOR <i>[Signature]</i> 10-01-94	QUALITY ASSURANCE 11-1-94 <i>[Signature]</i>	METHODS <i>[Signature]</i>	

SECTION 9.0 - DATA HANDLING (CON'T)

- 9.14 To transfer a file to LIMS, hit ESC to escape to main menu screen. Press F10, then press S (shell).
- 9.15 At the C:> prompt, type CD Data.
- 9.16 At the C:\Data > prompt, type: S2SEND Filename.IN 2
- 9.17 After file is sent, pick up report at the printer.
- 9.18 Fill out all log books and QC sheets daily for each analytical run.

SECTION 10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING

- 10.1 Remove stannous chloride reagent line 2 to 3 minutes before other reagent lines. Pre-rinse all lines in clean DI water.
- 10.2 Place all lines in clean DI water and flush instrument for 20 - 30 minutes.
- 10.3 Turn pump module off and unlatch platens.
- 10.4 Turn off light source and power module.
- 10.5 Dispose of all sample cups, clean the work area, and rinse and store all glassware.
- 10.6 If troubleshooting is necessary, refer to the RFA manual - troubleshooting section.
- 10.7 Consult the Supervisor or Quality Assurance Officer before making any major changes, adjustments or repairs to the instrument.

METHOD # 3120.1	METHOD NAME SILICA (SiO ₂)	REVISION 0.0	REFERENCE SM4500SID (Modified)	STATUS EFFECTIVE
PAGE 9 OF 9				DATE 9-26-94
APPROVALS/DATE				
DIVISION DIRECTOR  11/2/94	LAB SUPERVISOR  10-31-94	QUALITY ASSURANCE  11-94	METHODS 	

SECTION 11.0 - REFERENCES

- 11.1 EPA Methods for Analysis of Water and Wastes, EPA-600/4-79-020, March 1979.
- 11.2 Standard Method for the Examination of Water and Wastewater, 17th Edition.
- 11.3 ACS LIMS Users Guide, version 1.0, 1992
- 11.4 SFWMD Comprehensive Quality Assurance Manual, current version.
- 11.5 ALPKEM RFA 300 series Operator's Manual
- 11.6 ALPKEM RFA Software Manual

July 14, 1995

Page 1 of 2

Soluble Reactive Phosphate

Reference: Stannous Chloride - APHA Standard Methods, 15th ed., p. 417. Method 424E (1980).

Equipment: Spectrophotometer, Hach 2000

Reagents: The reagents for use in this procedure are purchased as part of a test kit, K-8513, Phosphate (M-Blue), 0-4 ppm, from Chemetrics, Inc., Route 28, Calverton, Virginia, 22016-0214. The ammonium molybdate (R-8513) is supplied in evacuated 13mm diameter glass ampoules. The stannous chloride (A-8500) is supplied in a plastic dropper bottle. Each kit contains enough reagent for testing 30 samples.

Standards: Stock solution = 1000 mg/l $\text{PO}_4\text{-P}$ = 4.394 grams potassium phosphate, monobasic (KH_2PO_4) dissolved in 1 liter deionized (18 megohm) water. This solution is prepared monthly in West Palm Beach by the Lab QA Officer.

Working standard solutions are prepared from the stock and preserved with 10 drops of 50% sulfuric acid before dilution to 1000 ml.

4.0 mg/L = 40 ml stock diluted to 1000 ml with di H_2O
2.0 mg/L = 20 ml stock diluted to 1000 ml with di H_2O
1.5 mg/L = 15 ml stock diluted to 1000 ml with di H_2O
1.0 mg/L = 10 ml stock diluted to 1000 ml with di H_2O
0.5 mg/L = 5 ml stock diluted to 1000 ml with di H_2O
Blank = 1000 ml di H_2O

Standard Additions (Spikes): 5.0 ml of the 2.0 mg/L working standard is added to 20 ml of sample. Standard addition concentration is equal to 0.80 times the samples concentration plus 0.5 mg/L.

QC Check Solutions (Known): Prepared monthly by the Laboratory Quality Assurance Officer in West Palm Beach.

Procedure: Sample and standard treatment: 25 mls of sample or standard is poured into a beaker. Two drops of stannous chloride (A-8500) is added and mixed well. The tip of the evacuated ampoule containing the ammonium molybdate solution is broken off under the surface of the sample. The ampoules then fill automatically. The ampoule is inverted several times to mix the sample and solution and the color is allowed to develop for at least 10 minutes but less than 30 minutes. The absorbance is measured by inserting the ampoule into the spectrophotometer.

Calibration: The spectrophotometer is calibrated by treating the working standards and deionized water blank by the procedure described above. The absorbance is measured at 690 nm and a calibration curve is generated.

Quality Control Procedures:

1. Calibration standards are run every 20 samples.
2. A standard addition is run every 10 samples.
3. A QC check samples is run every 10 samples.
4. A repeat (duplicate) is run every 10 samples.

The results from all quality control samples must fall within the current acceptable limit ranges.

General Description: This method measures reactive (ortho) phosphate on an unfiltered sample. Suspended sediment is not found to produce significant interference due to its settling out during the 10 - 30 minute color development time. Similarly, at this wavelength, sample color does not interfere. The sample values range between the laboratory values for total phosphorus and orthophosphate. The sensitivity of this procedure is 0.02 mg/L.

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 1 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>Mell</i> 11/12/94	LAB SUPERVISOR <i>J...</i> 11-1-94	QUALITY ASSURANCE <i>M...</i> 11-1-94	METHODS <i>DB Duerrhoff</i>	

SECTION 1.0 - TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
2.0 - METHOD DESCRIPTION/HISTORY	1
3.0 - SAFETY PRACTICES	2
4.0 - LIST OF EQUIPMENT/INSTRUMENTATION	2
5.0 - REAGENTS	2
6.0 - STANDARDS	2
7.0 - QUALITY CONTROL	2
8.0 - STEP-BY-STEP PROCEDURE	3
9.0 - DATA HANDLING	5-7
10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING	7
11.0 - REFERENCES	7

SECTION 2.0 - METHOD DESCRIPTION/HISTORY

2.1 Phosphatases can be associated with algal and bacterial cell wall and released into water by disintegrating algal cells. These enzymes catalyze the hydrolysis of phosphomonoesters to orthophosphate and an alcohol:

$$R-PO_4 + H_2O \rightleftharpoons ROH + H_2PO_4$$

2.2 Phosphatases are classified as either acid or alkaline, depending on the pH of the environment in which they exist. The determination of acid phosphatase activity (APA) and alkaline phosphatase activity is conducted at this native pH by adjusting the pH of the buffer solution.

2.3 The substrate used in this assay is methylumbelliferyl phosphate (MUP), which has a low background fluorescence, thus allowing assay of wide variety of concentration with very high sensitivity. The amount of substrate added is determined by preparing increasing amount of substrate solution. V_{max} is calculated as the optimum amount of substrate for enzymatic hydrolysis.

2.4 Basically, MUP is prepared in a pH adjusted buffer and added into the sample. The phosphatase enzyme that maybe present in the sample will hydrolyze MUP into methylumbelliferone and phosphate. Methylumbelliferone fluoresces at a specific wavelength when excited with UV light and can be quantified by a spectrophotometer or a fluorometer. A computer aided Cytofluor, a fluorescence plate scanner, is used in our laboratory to perform the analysis.

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 2 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/1/94	LAB SUPERVISOR <i>[Signature]</i> 10-11-94	QUALITY ASSURANCE <i>[Signature]</i> 11-1-94	METHODS <i>[Signature]</i>	

SECTION 3.0 - SAFETY PRACTICES

- 3.1 All personnel conducting this method should be familiar with the SFWMD Chemical Hygiene Plan and should have reviewed any pertinent Material Safety Data Sheets.
- 3.2 No hazardous wastes are generated by this procedure; liquid waste may be disposed of in the sink.
- 3.3 Follow the Personal Safety Protection Codes below during analysis:
- A = Lab coat, glasses, gloves
 - B = Lab coat, apron, acid resistant gloves, face shield+goggles+respirator or full-face respirator
 - C = Lab coat, acid resistant gloves, goggles, face shield, apron
 - D = Lab coat, acid resistant gloves, glasses
 - E = Lab coat, glasses, apron, gloves
 - F = Rubber acid carrier
 - G = Flush sink drain with ample amount of tap water

SECTION 4.0 - LIST OF EQUIPMENT/INSTRUMENTATION

- 4.1 Millipore CytoFluor 2350 Multiwell Fluorescence Plate Reader
- 4.2 IBM-compatible PC, with Windows 3.1 and Excel 3.0, and a printer
- 4.3 Cytofluor computer interface
- 4.4 Multiwell Low fluorescence plates (24 wells), opaqued
- 4.5 Multiwell pipetter, 0-150 μ L capacity
- 4.6 Eppendorf Micropipettor, adjustable 250 μ L-1000 μ L
- 4.7 Analytical balance, 0.1 mg sensitivity
- 4.8 Volumetric flasks, 10 and 1000 mL capacity
- 4.9 Graduated pipets, 1, 2, 5, 10 mL capacity
- 4.10 Freezer
- 4.11 Incubator, ambient to 40°C

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 3 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 10-11-94	QUALITY ASSURANCE 11-1-94 <i>[Signature]</i>	METHODS <i>[Signature]</i>	

SECTION 5.0 - REAGENTS

5.1 Tris stock buffer: Prepare a 0.1M solution of Tris (base) buffer. Add 12.11 g Tris to 1000 mL volumetric flask and add enough Millipore water to bring it into solution. Add 0.2037 g of anhydrous $MgSO_4$, to get a final concentration of 0.01M. (Safety PP = A)

5.2 Working buffer A: Measure 100 mL of Tris stock buffer into 1L volumetric flask. Bring to volume using sterile Millipore water. Adjust the pH to pH 8.00 for alkaline phosphatase, by slowly adding 1N HCl while stirring. The final concentration of this solution is 0.01M Tris/0.001M $MgSO_4$. (Safety PP = A)

5.3 Working buffer B: Measure 100 mL of Tris stock buffer into 1L volumetric flask. Bring to volume using sterile Millipore water. Adjust the pH 6.5 to for acid phosphatase, by slowly adding 1N HCl while stirring. The final concentration of this solution is 0.01M Tris/0.001M $MgSO_4$. (Safety PP = A)

5.4 Substrate : Methylumbelliferyl phosphate (MUP) (FW=256.2). Weigh 0.128 g of MUP and bring to 250 mL volume using the TRIS stock buffer. (Safety PP = A)

SECTION 6.0 - STANDARDS

6.1 Stock standard (1000 μ M MU): Dry approximately 1g of methylumbelliferone (MU) overnight at 105 °C. Weigh out 0.1982 g of oven-dried MU into a 1L volumetric flask and dilute to volume with working Tris buffer A or B. Keep at room temperature, in a dark container.

6.2 Secondary stock standard (10 μ M MU): Pipet 100 μ L of stock MU solution into a 10 mL volumetric flask and add working Tris buffer A or B to volume.

6.3 Working Standards

Std.	Methylumbelliferone conc.,		Vol. of Stock MU, mL	Final Volume, mL (Add Tris Buffer C to volume)
	μ M	nM		
S1	0.1	100	0.1	10
S2	0.3	300	0.3	10
S3	0.5	500	0.5	10
S4	1.0	1000	1.0	10
S5	2.0	2000	2.0	10

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 4 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 10-2-94	QUALITY ASSURANCE <i>[Signature]</i> 1-11-94	METHODS <i>[Signature]</i> DB Duanoff	

SECTION 7.0 - QUALITY CONTROL

- 7.1 Run a duplicate set of each working standard for each analysis. Record the fluorescence value on the instrument log.
- 7.2 Run a QC solution after a set of working standard and at the end of each analysis.
- 7.3 A repeat analysis should be conducted for every 20 samples analyzed. These results are recorded on QC sheet and instrument log. The mean and coefficient of variation of the replicate set is determined and noted on the QC sheet and physical parameters log.
- 7.4 All quality control data must be within the current established limits before entering sample data into the LIMS system. Consult the supervisor or QA officer if unable to obtain acceptable QC result.
- 7.5 Samples should be mixed thoroughly each time when taking an aliquot.

SECTION 8.0 - STEP-BY-STEP PROCEDURE

- 8.1 Create a workgroup for APA (Product). Refer to the ACS LIMS users Guide for instructions on "Creation of a Workgroup" and to Section 9.2.1 and Section 9.2.1. Remove the samples to be tested from the refrigerator and allow them to warm to room temperature (this may require several hours).
- 8.2 Turn on the Cytofluor and allow it to warm up for at least 15 minutes. Select the CYTOCALC Program (double click with the mouse).
- 8.3 Determine the pH of the water samples by referring to the hydrolab data or by laboratory measurement with a pH meter. Prepare a fresh set of working Tris buffer A or B, based on the pH range of the samples to analyze, and by following the procedure in Section 5.0. (Therefore, TRIS buffer A should be prepared and used if sample pH is greater than 7, and TRIS buffer B should be used when the sample pH is less than 7).
- 8.4 Prepare a fresh set of working standards, as described in Section 6.0.
- 8.5 The well plate template is displayed automatically when entering CYTOCALC program. If the 24-well plate template is not displayed, open the file 24well.CFL (file, open, 24 well.CFL), or the most recent run file with 24 well. Edit the plate protocol by entering the last four digits of the sample number. Use one plate for every 10 samples. It is important to have a duplicate of each standard for the CYTOCALC to function.

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 5 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/2/94	LAB SUPERVISOR <i>[Signature]</i> 10-11-94	QUALITY ASSURANCE 1-11-94 <i>[Signature]</i>	METHODS <i>[Signature]</i>	

SECTION 8.0 (CON'T)

8.6 Alternatively, a new protocol can be created by following these steps:

- a. On the well displayed on the screen, click the well you want to use.
- b. On the upper right hand corner of the screen, click the well assignment you need (Bl=Blank, Un=Sample, St=Standard, Co=Control, Em=Empty)
- c. On the edit bar located at the upper portion of the screen, edit the correct I.D. You will need to enter the concentration values of each standard and control.

8.7 Measure 2 mL of blank, standard and QC into each assigned well.

8.8 Pipet 1.8 mL of samples according to the tray protocol. Add 200 μ L of MUP into each sample well. **DO NOT ADD MUP SUBSTRATE TO STANDARD, BLANK AND QC WELLS.**

8.9 Immediately place the well plate in the door transport, ensure proper seating of the plate and that the first well is on the top right position. Scan within 10 seconds (Click RUN, COLLECT DATA, then enter WORKGROUP number under file name). Record the exact time of the start of the incubation. The system will automatically scan for time 0 reading, then a message "50 % completed will be displayed" on the scan window. (This means that a second scan will resume after 1800 seconds or 30 minutes).

8.10 At the end of the second scan, the system will prompt a message, "100 % completed", and will give you an option to CLOSE DOOR or IGNORE DOOR. Choose IGNORE DOOR if another plate is to be scanned, and CLOSE DOOR if no more plate is to be scanned.

8.11 At this stage, a message may appear "Cytofluor setting has been changed". Hit OK to display the plate. Note that the fluorescence data will be displayed on the plate layout.

8.12 On the top right corner of the screen, depress the LINK button. Note that any associated well (standard, blanks, QC and sample well) is highlighted (black background). Double click each well that is not highlighted. **IMPORTANT:** Only highlighted well (blackened background) will be calculated with the standard curve.

8.13 After highlighting each well; depress the DATA button on the right hand corner of the screen. The system will prompt "ACCEPT CURRENT LINK?". If this is your final link, select YES to proceed. Otherwise, depress the plate button on the upper right hand corner of the screen and make any ID corrections. (At this stage, any outlier standard or sample can be hidden to exclude from calculations. Press the HI button).

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 6 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>Mell</i> 11/12/94	LAB SUPERVISOR <i>J. [unclear]</i> 11/08/94	QUALITY ASSURANCE <i>Man Smith</i> 11/08/94	METHODS <i>OPDuanoff</i>	

SECTION 9.0 - DATA HANDLING

9.1 Cytofluor PC-Data Handling

8.14 (After LINK command from Steps 8.12 and 8.13) On the menu bar, select DATA, FORMAT. Specify SCAN 1, then press CALCULATE (Ensure that Calculate and Statistics are marked with X on the screen).

8.15 The computer will automatically go to Excel Program, and display the worksheet with the calculated values. Enter the complete LIMS ID for each sample and QC.

8.16 Select FILE, PRINT to print the time 0 worksheet. Select FILE, SAVE AS, and give file name (.csv) to save file

8.17 Select FILE, END DATA DISPLAY to return to Calculate window. Specify Scan 2. Repeat steps 8.14 to 8.16 to calculate, print and save results of SCAN 2. Save as Filename.csv.

9.2 APA Calculation (PC Windows-Excel)

9.2.1 To calculate final APA values in nM/min-mL, the equation below is used:

$$\text{APA (nM/min-mL)} = \frac{\text{MU conc. at time 0 } (\mu\text{M}) - \text{MU conc. after 30 minutes } (\mu\text{M})}{(30 \text{ minutes} \times 1.8 \text{ mL})} \times 1000$$

This equation is entered into Excel worksheet, with filename APA.XLM. Calculation can be done in these worksheet by opening this file, then the time 0 and time 30 minutes files. Copy and paste the time 0 values into APA.XLM worksheet. Repeat to copy time 30 values into the calculation worksheet. Calculation will be automatic. Review the QC and blank values to ensure accurate calculations.

9.2.2 The calculated APA values can be sent to LIMS following Step 9.2.2.

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 7 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 10-1-94	QUALITY ASSURANCE <i>[Signature]</i> 1-11-94	METHODS <i>[Signature]</i> Q. B. Duanoff	

SECTION 9.0 - DATA HANDLING (CON'T)

9.2 LIMS Database

9.2.1 Prior to running samples: Go into the ACS LIMS via the command `sl` and create a workgroup for the test APA and exit the SEEDPAK1 Main Menu (see Creating a Workgroup and Printing a Workgroup in the ACS LIMS Users Guide).

9.2.2 After running the samples: Examine the results for acceptability of QC samples and repeats. either by manual entry or by transfer from disk through the LIMS Database Analyst. To enter manually, go the WORKSTAT, MANUAL ENTRY, APA. Enter the workgroup number. This will prompt the sample numbers one by one. Enter the concentration results.

9.2.2 Pick up the LIMS data entry report for the samples you have just entered from the system printer and examine against the raw data report. Submit any necessary corrections to the DATA UNIT by filling a LIMS Database Correction form.

SECTION 10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING

10.1 Wipe dry the surfaces of the plate incubator, the Cytofluor and the work areas. Use D.I. water to clean any spill on the instrument.

10.2 Rinse the plate thoroughly with dilute Liquinox™ and D. I. water. Place the plate upside down on a paper towel, to dry. Discard the plate if any visible crack or scratches is observed.

SECTION 11.0 - REFERENCES

11.1 Pettersson, K. and M. Jansson. 1978. Determination of phosphatase activity in lake water-a study of methods. Verh. Internat. Verein. Limnol. 20:1226-1230.

11.2 Prof. Robert G. Wetzel. 1994. Personal Communication. Department of Biological Sciences, The University of Alabama, Tuscaloosa, Alabama 35487-0344, USA.

11.3 Cytofluor 2300 Manual. 1992. Millipore Corporation, Bedford, MA.

METHOD 3160.1	METHOD NAME ACID AND ALKALINE PHOSPHATASE ACTIVITY (APA)	REVISION 00	REFERENCE Pettersson and Jansson, 1978	STATUS EFFECTIVE
PAGE 8 OF 8				DATE 10-01-94
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i> 11/12/94	LAB SUPERVISOR <i>[Signature]</i> 10-28-94	QUALITY ASSURANCE 11-11-94 <i>[Signature]</i>	METHODS <i>[Signature]</i>	

FIGURE 1.0 SAMPLE MULTIWELL PROTOCOL (24-well plate)

(A1) S1 (2000 nM)	(A2) S2 (1000 nM)	(A3) S3 (500 nM)	(A4) S4 (300 nM)	(A5) S5 (100 nM)	(A6) BLANK
(B1) S1 (2000 nM)	(B2) S2 (1000 nM)	(B3) S3 (500 nM)	(B4) S4 (300 nM)	(B5) S5 (100 nM)	(B6) BLANK
(C1) QC1	(C2) QC2	(C3) Sample	(C4) Sample	(C5) Sample	(C6) Sample
(D1) Sample	(D2) Sample	(D3) Sample	(D4) Sample	(D5) Sample	(D6) RPT

DER Biology Section

v. 1-7/6/93

SOP #BA-7
Benthic Macroinvertebrate
Dip Net Sample Collection

(based on Plafkin, et al, 1989, *Rapid bioassessment protocols for use in streams and rivers: benthic macroinvertebrates and fish*, EPA/444/4-89-001)

STEPS**COMMENTS****Materials**

1. *Field Physical / Chemical Characterization Data Sheet*
2. *Habitat Assessment Sheet*
3. Dip Net with No. 30 mesh
4. 4-liter wide-mouth plastic jugs
5. 100% formalin

Prepare per SOP #BA-2.1

Methods

1. Visually examine the area or reach to be sampled. You must either walk or boat throughout the aquatic system, paying close attention to its physical and habitat characteristics.
2. Fill out *Field Physical / Chemical Characterization Data Sheet* and *Habitat Assessment Sheet*. The percent coverage of substrate type refers to how much of each habitat type is actually present at the sampling site.
3. Determine the number of sweeps to perform in each habitat type out of the 20 total sweeps per station. This requires a two step process. First, select the "major" or "most productive" habitats for the stream type. Use the following formula to calculate the number of sweeps in each habitat type:

$$\text{Number of sweeps per Major Habitat} = \frac{20}{(\text{Number of Major Habitats}) + 1}$$

The result is rounded to the nearest integer. The remaining number of sweeps (to make a total of 20) is evenly divided among the minor habitats (such as sand, mud, or muck in most cases).

In fairly small (1st to 4th order) streams, the length of a discrete station should consist of a 100 m stretch of stream, and the width should be from bank to bank. In very large systems it may be necessary to establish more than one station to adequately characterize the biota.

See SOP's #BA-17 and BA-18 for instructions on filling out these forms.

It is important to accurately determine the spatial extent of each substrate type (in a 3 dimensional context) for habitat scoring procedures.

Generally, the most (to least) productive habitat types are as follows: snags, aquatic vegetation, leaf packs, roots, undercut banks, rocky outcrops, muck, and sand. All but the last two can be considered "major" or "productive".

Example: If 3 major habitat types are present, perform 5 sweeps in each of these habitats and divide the remaining 5 sweeps up among the other non-major habitats, so that a total of 20 sweeps are performed. If 4 major habitats are present, perform 4 sweeps in each of these, then 4 divided up among the remaining types. For 5 major habitats do 3 in each major habitat and divide the remaining 5 sweeps up among the other non-major types.

Proper interpretation of benthic collections requires that samples be collected from multiple habitats that are representative of the site. If possible, the same habitats should be sampled at reference and test sites the same number of times to isolate the effects of water quality on the benthic community.

DER Biology Section

STEPS

4. Perform 20 discrete 0.5 meter sweeps with the dip net. Sample the available substrates as determined by the above procedures.
 - a. In streams with sufficient water velocity, the most effective way to capture invertebrates is to place the bottom rim of the dip net downstream of the area to be sampled. Disturb, agitate, or dislodge organisms (with hands and/or feet) from substrates (snags, etc.) at a distance of 0.5 m upstream of the net.
 - b. For areas without flow, disturb an area of substrate that is one dip net width wide and approximately 0.5 m long, and sweep the net over the area a few times to ensure the capture of organisms which were living there.
 - c. For heavily vegetated areas (some streams, lake margins, or wetlands) jab the net into the base of the vegetation, digging down to the substrate, and dislodge organisms using a one-half meter sweeping motion with the net.
 - d. Sample leaf packs (if present) by disturbing leaf pack areas with hands or feet before scooping one-half meter worth of material into the net.
 - e. Sand, muck, mud, and silt (non-major habitats) can be sampled by taking 0.5 meter sweeps with the net while digging into the bottom approximately 1 cm.
5. Record the number of sweeps for each habitat on the *Field Physical/Chemical Characterization Data Sheet*.
6. Reduce the sample volume after each discrete sample by dislodging organisms from larger debris (but retaining invertebrates in the net or sieve) and discarding the debris. Save finer debris plus organism mixture in large wide mouth jugs. Try to reduce enough of the sample volume in the field so that no more than 2 gallons of material are collected. If this is not possible, put the material into additional jugs. Sample reduction is easier in the laboratory.

v. 1—7/6/93

COMMENTS

When performing an upstream/downstream type of study, sample the downstream station first to prevent upstream invertebrates from drifting into a location they were not originally inhabiting.

Catch organisms by allowing them to flow into the net and also by sweeping the net towards disturbed material.

Several sweeps over the same 0.5 meter area are recommended to make sure all organisms are captured. This sampling effort in a discrete 0.5 meter spot is considered as 1 sweep.

Where a continuous half meter sweep is impossible, take 2 quarter meter sweeps in the same area to attain a full 0.5 meter sweep.

If the net is pushed too deep in coarse sand, very little of the sand will be washed through the net resulting in a sample that contains few organisms and is hard to process.

Ideally, control and test sites will be sampled the same number of times in the same habitats.

The relative proportions of the organisms collected must be maintained intact to calculate many community metrics. Some field picking of delicate organisms is acceptable as long as community composition is not altered.

Indicate on the label how many jugs the entire sample is contained in, e.g., "1 of 2", "2 of 2".

DER Biology Section

STEPS

- 7. If laboratory sorting can be performed within 24 hours, place samples immediately on ice. Cold temperatures should slow organisms enough to prevent predation (and subsequent alteration of community structure). If sorting will be delayed, preserve with 10% formalin (do this by adding one part of 100% formalin to the jug with 9 parts ambient water). After organisms have been removed from detritus, they should be placed into 70% ethanol.

v. 1-7/6/93

COMMENTS

If organisms are too active during sorting, pour some carbonated water, clove oil, or other relaxing agent into the sample.

Samples that will not be sorted within 2 days should be preserved in formalin. Ethanol alone will not prevent the vegetative debris from decomposing, resulting in a sample that is very unpleasant to sort.

DER Biology Section

v. 1—7/6/93

SOP #BA-8
Benthic Macroinvertebrate Qualitative (Dip Net)
Sample Handling

STEPS

COMMENTS

Materials

- | | |
|---|--------------------------|
| <ol style="list-style-type: none">1. Waterproof paper and permanent marker2. U.S. 30 mesh sieve3. U.S. 10 mesh sieve4. Ethanol filled squeeze bottle (80%)5. White enamel pan, marked with a grid of 5 cm squares6. List of random numbers7. 250-mL glass jar8. Dissecting microscope9. 100 x 15 mm petri dish10. Forceps11. Vials for picked organisms (1 or 2 dram)12. Laboratory counter13. <i>Benthic Macroinvertebrate Bench Sheet</i> | <p>For making labels</p> |
|---|--------------------------|

Methods

- | | |
|--|--|
| <ol style="list-style-type: none">1. Check labels so you know which sample you are dealing with (control site, test site, etc.). Make a very clear label to go into the bottle of picked bugs with station identification, date sampled, replicate number, and your initials.2. Place a portion of the contents of the sample (fist sized) into a U.S. 10 mesh sieve with a U.S. 30 mesh sieve underneath.3. Rinse with tap water (a small hose attached to the faucet works best), spraying organisms and small detritus down into the U.S. 30 mesh sieve. Visually inspect large debris (leaves, plants, twigs) held in the U.S. 10 mesh sieve for animals before discarding. Wash fine debris (silt, mud) through the bottom (U.S. 30 mesh) sieve. Repeat procedures #2 and #3 until all the sample from all the jugs has been processed.4. Place sample in gridded pan. Each 5 cm grid should have a pre-assigned number. Liquid present in the sample should be sufficiently reduced to prevent material from shifting among grids during the sorting process. | <p>Make sure that you know how many containers in which the particular sample is stored (there may be several jugs). The entire sample must be included in this reduction and homogenization process.</p> <p>This inspection is best accomplished by placing the debris in a white pan and observing it with the Luxo[®] lighted magnifier. Organisms found (generally the ones too large to pass through the U.S. 10 mesh sieve) should be placed into the U.S. 30 mesh sieve with the rest of the unpicked sample.</p> <p>There are 24 total 5 cm grids in a standard white enamel pan.</p> |
|--|--|

DER Biology Section

STEPS

v. 1—7/6/93

COMMENTS

- | | |
|--|--|
| 5. Thoroughly mix the sample so that a homogenous distribution of organisms is achieved in the detrital matrix. | |
| 6. Select a grid using the random number table. Remove the contents of the entire grid and place in a glass jar. | Use a ruler to delineate the edges of a grid while removing the sample. |
| 7. Take a small amount of this detritus plus organism mixture and place it in the bottom portion of a petri dish. Add enough liquid (usually ethanol) so that the material is wet enough to move around easily with forceps. | The top portion of the petri dish may later be placed over the sample to prevent desiccation if you must leave the sample overnight. |
| 8. Place the petri dish under a dissecting scope set at low power (approximately 7x or 10x). In a deliberate, systematic manner scan back and forth or up and down picking each and every organism from the aliquot and placing it into an alcohol-filled vial (clearly identified as per step #1). | Using a pair of forceps in each hand enables you to better tease organisms out of fibrous detritus. Forceps should be sharp and properly aligned.

Use the laboratory counter to keep a running total of the number of organisms picked. |
| 9. After scanning dish in one pattern (e.g., up and down), go back through using a different pattern (back and forth) to assure that all organisms have been removed from the aliquot. | Picking accuracy should be checked by a co-worker in 10 % of the samples. |
| 10. Continue steps 6, 7, 8, and 9, until you obtain a minimum of 100 organisms. Once a grid is selected, its entire contents must be sorted. | If an obvious organism is observed but its grid number was not selected and no examples of that organism were present in grids which were selected, that organism may be noted as qualitatively observed. The organism should NOT be included in the analysis. |
| 11. Record the information requested on the <i>Macroinvertebrate Lab Bench Sheet</i> which includes site, laboratory sample number, STORET station number, sample type, replicate number, and date collected. Include the initials of the persons who collected and sorted the sample. Record the number of grids selected (e.g., "4 of 24") to enable conversion to total abundance present in the original sample. | Failure to record the number of grids selected (out of the total grids possible) seriously compromises the usefulness of the data. |

DER Biology Section

v. 1—7/6/93

SOP #BA-8.1
Preparation of 80% Ethanol**STEPS****COMMENTS****Materials**

1. 100% Ethanol (HPLC grade)
2. D.I. water
3. 4000-mL graduated cylinder
4. 4-L glass amber jug
5. Plastic funnel to fit into 4-L jug

Methods

1. Fill graduated cylinder with 3200 mL of 100 % ethanol.
2. Add D.I. water to the graduated cylinder until the total volume is 4000 mL.
3. Using the funnel, transfer the dilute ethanol to the 4-L jug which should be properly labeled as 80% ethanol.
4. Rinse graduated cylinder and funnel with D.I. water and return them to the shelf above the sink.

This operation is easier with two people as the ethanol is purchased in 5-gallon cans. One person holds the cylinder while the other pours the ethanol.

DER Biology Section

v. 1—7/6/93

SOP #BA-9
Benthic Macroinvertebrate
Grab Sample Collection
(Modified from *Standard Methods 10500B.3*)

STEPS

COMMENTS

Materials

1. Ekman or Petite Ponar dredge
2. U.S. 30 mesh box sieve
3. White enamel pan
4. Plastic squeeze bulb ("turkey baster")
5. Small bucket
6. Wide mouth plastic sample containers
7. Tape and permanent markers

For making labels

Methods

1. Use of the Ekman dredge is restricted to sampling soft substrates (silt, muck) in areas with little current. The Ponar dredge may be used for sampling under these conditions and also in areas with a harder substrate (rocks, shells, sand).
2. When sampling from a boat, dredge samples are collected from the rear and downstream of the vessel to avoid contamination of other types of samples with disturbed sediments. Rinse the box sieve with ambient water and tie it to the side of the boat where samples will be collected.
3. Ekman: Open the spring-loaded jaws and attach the chains to the pegs at the top of the sampler. Lower the dredge to the bottom, making sure it settles flat. Holding the line taught, send down the messenger to close the jaws of the Ekman dredge. Pull the sampler to the surface and place it immediately into the box sieve. Carefully open the jaws and disgorge the contents into the sieve, rinsing to assure complete sample purging.

The number of replicates collected is dependent upon several factors including the area sampled by the device, the purpose of the study, and the degree of patchiness in the distribution of the organisms at the site. Routinely, take 3 dredges. All replicates are placed in separate sample containers (for statistical analyses). If it has been determined that you are sampling in an exceptionally depauperate area, additional replicates may be required (pilot study needed). In that case the number of replicates sampled at the group of stations you wish to compare should be equal.

The box sieve is constructed of fiberglass-coated wood and U.S. 30 mesh screen. When placed in the water, it will float at the surface. If a box sieve is unavailable, the dredged material may be washed in the dip net providing it is fitted with a U.S. #30 mesh sieve material. The disadvantage of using the dip net is that it requires 2 people (1 to hold the net, 1 to manipulate the dredge).

The spring-loaded Ekman is dangerous. Hold the dredge firmly above the hinges, and be very careful that no body parts get pinched by the snapping jaws, which could produce serious injury.

Check to make sure the jaws are fully closed and that no sample was lost while lifting the dredge.

DER Biology Section

STEPS

3. **Ponar:** Open the jaws and place the cross bar into the proper notch. Lower the dredge to the bottom, making sure it settles flat. When tension is removed from the line, the cross bar will drop, enabling the dredge to close as the line is pulled upward during retrieval of the dredge. Pull the Ponar to the surface and place it immediately into the box sieve. Carefully open the jaws and discharge the contents into the sieve, rinsing to assure complete sample purging.
4. Swirl the box sieve in the water with a back and forth motion to wash the fine sediments through. Concentrate the remaining sample into one corner of the sieve.
5. Fill the small bucket with ambient water and use this water to fill the squeeze bulb. Using the squeeze bulb, rinse the sample from the sieve to the enamel pan.
6. Transfer the sample from the enamel pan into the pre-marked wide mouth jug (again, using the squeeze bulb), making sure the location, date, and replicate number is accurate.
7. Preserve the sample with 10% formalin by adding a 10 to 1 ratio of water to 100% formalin. If laboratory processing is possible within 8 hours, the samples may be stored on ice, without addition of formalin.

v. 1—7/6/93

COMMENTS

Check to make sure the jaws are fully closed and that no sample was lost while lifting the dredge.

If a sediment type is especially clayey or mucky, it may be necessary to use a hand to break up clumps and agitate the sample to reduce it. Make sure you rinse any detritus from your hand back into the sieve.

Take care to rinse the entire contents of the sample into the pan. Some organisms may stick to the screen.

Rose bengal dye may be added to the sample, as a picking aid, if desired.

DER Biology Section

v. 1—7/6/93

SOP #BA-10
Benthic Macroinvertebrate Grab
Sample Handling
 (Modified from *Standard Methods 10500C*)

STEPS**COMMENTS****Materials**

- | | |
|---|---|
| <ol style="list-style-type: none"> 1. Waterproof paper and permanent marker 2. U.S. 30 mesh sieve 3. Ethanol filled squeeze bottle (80%) 4. Glass jars 5. Dissecting microscope 6. 100 x 15 mm petri dish 7. Forceps 8. Vials for picked organisms (1 or 2 dram) 9. <i>Benthic Macroinvertebrate Bench Sheet</i> | <p>For making labels</p> <p>Prepare per SOP #BA-8.1</p> |
|---|---|

Methods

- | | |
|--|--|
| <ol style="list-style-type: none"> 1. Check labels so you know which sample you are dealing with (control site, test site, etc.). Make a very clear label to go into the bottle of picked bugs with station identification, date sampled, replicate number, and your initials. 2. Pour the contents of the wide mouth jug over a U.S. 30 mesh sieve. Rinse the jug with tap water to make sure all organisms are put into the sieve. 3. Rinse with tap water (a small hose attached to the faucet works best). Wash fine debris (silt, mud) through the sieve. Any large debris (leaves, shells, etc.) present should be brushed clean of organisms and discarded. Rinse the organism plus detritus mixture to one small area of the sieve. 4. Using an ethanol-filled (80%) squeeze bottle, rinse the organism plus detritus matrix into the smallest practical container (usually a 100 mL to 250 mL glass jar). Put the label inside the jar. 5. Record the information requested on the <i>Macroinvertebrate Lab Bench Sheet</i>, which includes site, laboratory sample number, STORET station number, sample type, replicate number, and date collected. Include the initials of the persons who collected and sorted the sample. | <p>It is critical that no errors be made during this step. If samples are mixed up, the entire study could be rendered invalid.</p> <p>Place jars on the sample shelf so that samples for a given study are organized together and clearly marked.</p> |
|--|--|

DER Biology Section

STEPS

6. Take a small amount of detritus plus organism mixture and place it in the bottom half of a petri dish. Add enough liquid (usually ethanol) so that the material is wet enough to move around easily with forceps.
7. Place the petri dish under a dissecting scope set at low power (approximately 7x or 10x). In a deliberate, systematic manner scan back and forth or up and down, picking each and every organism from the aliquot and placing it into an alcohol-filled vial (clearly marked as per step #1).
8. After scanning dish in one pattern (e.g., up and down), go back through using a different pattern (back and forth) to assure that all organisms have been removed from the aliquot.
9. Continue steps 6, 7, and 8 until the sample is finished.

v. 1-7/6/93

COMMENTS

The top portion of the petri dish may later be placed over the sample to prevent desiccation if you must leave the sample overnight.

Do not overload the petri dish with too much sample this results in sloppy work.

Using a pair of forceps in each hand enables you to better tease organisms out of fibrous detritus. Forceps should be sharp and properly aligned.

Picking accuracy should be checked by a co-worker in 10% of the samples.

DER Biology Section

v. 1-7/16/93

SOP #BA-11
Benthic Macroinvertebrate Core
Sample Collection
(Modified from *Standard Methods 10500B.3*)

STEPS

COMMENTS

Materials

1. Coring Device
2. U.S. 30 mesh box sieve
3. White enamel pan
4. Plastic squeeze bulb
5. Small bucket
6. Wide mouth plastic sample containers
7. Tape and permanent markers

For making labels

Methods

1. Use of coring devices is restricted to sampling fairly soft substrates (silt, muck, with only small amounts of sand or shell) usually in marine systems. The Biology Section uses two sizes of coring devices.
2. When sampling from a boat, use the 4 inch diameter coring device that is attached to a long pole, and has a valve near the top. Core samples are collected from the rear and downstream of the vessel to avoid contamination of other types of samples with disturbed sediments. Rinse the box sieve with ambient water and tie it to the side of the boat where samples will be collected.
3. Lower the coring device to the bottom with the valve open. After quickly pushing the device into the sediments, close the valve. The resulting vacuum will keep the material in the tube as it is raised up to the boat.
4. When collecting samples in wadable waters, a smaller coring device (2 inch diameter) can be used. This corer utilizes a flapper-valve equipped stopper which is inserted into the top of the pipe. Vacuum inside the pipe holds the material until the stopper is removed.

The number of replicates collected is dependent upon several factors, including the area sampled by the device, the purpose of the study, and the degree of patchiness in the distribution of the organisms at the site. Routinely, take enough cores so that an area equivalent to 3 Ponar dredges is collected (approximately 675 cm²). With our large (4 inch diameter) coring device, collect 8 replicates to achieve this. All replicates are routinely placed in separate sample containers (for statistical analyses). Depending on the study objectives, replicates may also be composited, as long as the number of replicates is equal for each station and clearly recorded so that the number of organism per square meter can be calculated.

The box sieve is constructed of fiberglass-coated wood and U.S. 30 mesh screen. When placed in the water, it will float at the surface. If a box sieve is unavailable, the dredged material may be washed in a dip net, providing it is fitted with a U.S. #30 mesh sieve material. The disadvantage of using the dip net is that it requires two people (one to hold the net, and one to manipulate the coring device).

Many clean water organisms are somewhat motile and may elude capture if you are not quick during sampling.

This small corer should be used primarily for non-biological sediment sampling (grain size, metals, etc.), as it is thought to be too small to effectively capture many organisms (e.g., crustaceans or tubicolous worms which are generally large in size) considered useful in impact determination.

DER Biology Section

STEPS

5. Pull the sampler to the surface, open the valve or remove the stopper, and place it immediately into the box sieve. Disgorge the contents into the sieve, rinsing to assure complete sample purging.
6. Swirl the box sieve in the water with a back and forth motion to wash the fine sediments through. Concentrate the remaining sample into one corner of the sieve.
7. Fill the small bucket with ambient water and use this water to fill the squeeze bulb. Using the squeeze bulb, rinse the sample from the sieve to the enamel pan.
8. Transfer the sample from the enamel pan into the pre-marked wide mouth jug (again, using the squeeze bulb), making sure the location, date, and replicate number is accurate.
9. Preserve the sample with 10% formalin by adding a 10 to 1 ratio of water to 100% formalin. If laboratory processing is possible within 8 hours, the samples may be stored on ice, without addition of formalin.

v. 1—7/6/93

COMMENTS

If a sediment type is especially clayey or mucky, it may be necessary to use a hand to break up clumps and agitate the sample to reduce it. Make sure you rinse any detritus from your hand back into the sieve.

Take care to rinse the entire contents of the sample into the pan. Some organisms may stick to the screen.

Rose bengal dye (use a very small amount) may be added to the sample, as a picking aid, if desired.

DER Biology Section

v. 1-7/6/93

STEPS**COMMENTS**

- | | |
|--|--|
| <p>6. Take a small amount of detritus plus organism mixture and place it in the bottom half of a petri dish. Add enough liquid (usually ethanol) so that the material is wet enough to move around easily with forceps.</p> | <p>The top portion of the petri dish may later be placed over the sample to prevent desiccation if you must leave the sample overnight.</p> <p>Do not overload the petri dish with too much sample. this results in sloppy work.</p> |
| <p>7. Place the petri dish under a dissecting scope set at low power (approximately 7x or 10x). In a deliberate, systematic manner, scan back and forth or up and down, picking each and every organism from the aliquot and placing it into an alcohol-filled vial (clearly marked as per step #1).</p> | <p>Using a pair of forceps in each hand enables you to better tease organisms out of fibrous detritus. Forceps should be sharp and properly aligned.</p> |
| <p>8. After scanning dish in one pattern (e.g., up and down), go back through using a different pattern (back and forth) to assure that all organisms have been removed from the aliquot.</p> | |
| <p>9. Continue steps 6, 7, and 8 until the sample is finished.</p> | <p>Picking accuracy should be checked by a co-worker in 10% of the samples.</p> |

DER Biology Section

v. 1—7/6/93

SOP #BA-13
Benthic Macroinvertebrate Hester-Dendy
Sample Collection
(Modified from *Standard Methods 10500B.5*)

STEPS

COMMENTS

Materials

1. 3 or 4 Hester-Dendy artificial substrates
2. Customized Hester-Dendy block, with coupling nuts for attachment of HD samplers and eye bolts for attachment of cable
3. Stainless steel cable
4. Nico-press[®] tool with fasteners
5. Whirl-pak[®] bags
6. Permanent marker

For labeling bags

Methods

1. Attach three HD's to the HD block, and place the block at a depth of one meter (or the deepest spot available if shallower than one meter). Take care to place control and test site blocks in areas of similar flow and habitat type. Space for 4 H-D's has been provided on the block, for use in studies requiring additional replication.
2. Attach cable to a point on the bank sufficiently high to enable recovery even if the water level increases.
3. After a 28 day incubation period, recover the HD samplers. Approach the block carefully, without disturbance, from the downstream position. In a deliberate, gentle manner lift the block straight up from the bottom and immediately place on a flat surface.
4. Quickly place the Whirl-pak bags over all the HD's, and unscrew them from the block. If an organism is observed crawling off a HD, capture it and put it in the appropriate Whirl-pak. Fill the Whirl-paks with ambient water (so that all the plates are wet), secure them (twirl three times and twist the ends), and place on ice.

Knowledge of the system's hydrologic regime is important to make sure samplers will not go dry during the 28 day incubation period. For example, if it is flood stage and you expect the water to drop 2 meters in the next few weeks, place sampler so that it will be 1 meter deep at the end of incubation.

In shifting sand substrates place the block so that existing snags will deflect sand from being deposited on the samplers. This can be determined by close examination of the bottom topography.

Wrap the cable around the base of a tree on the bank and use the Nico-press[®] tool and fasteners to secure the block. If vandalism is a potential problem, attempt to conceal the cable so that no one but you can find it. If the Nico-press[®] tool is unavailable, the fasteners may be cramped by hammering (two hammers are needed).

Wade or use a boat. DO NOT pull the block up from the shore.

Whirl-paks should be pre-labeled with the station, sample date, and replicate number, using the permanent marker.

Samples should NEVER be preserved until after organisms are scraped from the Hester-Dendy plates. Preservatives will poison the plates, preventing them from being used again.

DER Biology Section

v. 1—7/6/93

SOP #BA-14
Benthic Macroinvertebrate Hester-Dendy
Sample Handling

(Modified from *Standard Methods 10500C*)

STEPS**COMMENTS****Materials**

- | | |
|--|---|
| <ol style="list-style-type: none"> 1. Waterproof paper and permanent marker 2. U.S. 30 mesh sieve 3. Wrench for dismantling HD 4. Ethanol filled squeeze bottle (80%) 5. Glass jars 6. Dissecting microscope 7. 100 x 15 mm petri dish 8. Forceps 9. Vials for picked organisms (1 or 2 dram) 10. <i>Benthic Macroinvertebrate Bench Sheet</i> | <p>For making labels</p> <p>Prepare per SOP #BA-8.1</p> |
|--|---|

Methods

- | | |
|---|---|
| <ol style="list-style-type: none"> 1. Check labels so you know which sample you are dealing with (control site, test site, etc.). Make a very clear label to go into the bottle of picked bugs with station identification, date sampled, replicate number, and your initials. | <p>It is critical that no errors be made during this step. If samples are mixed up, the entire study could be rendered invalid.</p> |
| <ol style="list-style-type: none"> 2. Pour the contents of the Whirl-pak bag (an assembled Hester-Dendy and associated detritus) over a U.S. 30 mesh sieve. Using a wrench to remove the bottom nut, dismantle the Hester-Dendy. | <p>Rinse the Whirl-pak bag with tap water to make sure all organisms are put into the sieve.</p> |
| <ol style="list-style-type: none"> 3. Rinse with tap water (a small hose attached to the faucet works best). Scrape and simultaneously rinse organisms off HD plates with fingers (or a soft brush), using care not to damage the organisms. Wash fine debris (silt, mud) through the sieve. Any large debris (leaves, twigs) present should be brushed clean of organisms and discarded. Rinse the organism plus detritus mixture to one small area of the sieve. | <p>Save the HD plates and hardware, and place them into the drying oven.</p> |
| <ol style="list-style-type: none"> 4. Using an ethanol-filled (80%) squeeze bottle, rinse the organism plus detritus matrix into the smallest practical container (usually a 100 mL to 250 mL glass jar). Put the label inside the jar. | <p>Place jars on the sample shelf so that samples for a given study are organized together and clearly marked.</p> |

July 14, 1995
Page 17 of 18

DER Biology Section

STEPS

v. 1—7/6/93
COMMENTS

5. Record the information requested on the *Macroinvertebrate Lab Bench Sheet*, which includes site, laboratory sample number, STORET station number, sample type, replicate number, and date collected. Include the initials of the persons who collected and sorted the sample.
6. Take a small amount of detritus plus organism mixture and place it in the bottom half of a petri dish. Add enough liquid (usually ethanol) so that the material is wet enough to move around easily with forceps.

The top portion of the petri dish may later be placed over the sample to prevent desiccation if you must leave the sample overnight.


Do not overload the petri dish with too much sample, as this results in sloppy work.
7. Place the petri dish under a dissecting scope set at low power (approximately 7x or 10x). In a deliberate, systematic manner, scan back and forth or up and down, picking each and every organism from the aliquot and placing it into an alcohol-filled vial (clearly marked as per step #1).

Using a pair of forceps in each hand enables you to better tease organisms out of fibrous detritus. Forceps should be sharp and properly aligned.
8. After scanning dish in one pattern (e.g., up and down), go back through using a different pattern (back and forth) to assure that all organisms have been removed from the aliquot.
9. Continue steps 6, 7, and 8 until the sample is finished.

Picking accuracy should be checked by a co-worker in 10% of the samples.

Macrobenthic Invertebrate Keys

- Merritt, R.W. and K. W. Cummins, eds., 1984. An Introduction to the Aquatic Insects of North America. Second Edition. Kendall/Hunt Publishing Company, Dubuque, IA. 722 pp.
- Heard, W. H. 1979. Identification Manual of the Freshwater Clams of Florida. State of Florida Department of Environmental Regulation; Technical Series Vol 4, No. 2. 83 pp.
- Brigham, A.R., M.V. Brigham, and A. Gnilka. 1982. Aquatic Insects and Oligochaetes of North and South Carolina. Midwest Aquatic Enterprises, Mahomet, IL.
- Berner, L. and M.L. Pescador. 1988. The Mayflies of Florida (Revised edition). Univ. of Fla. Press, Gainesville, FL. 267 pp.
- Thompson, F.G. 1984. The Freshwater Snails of Florida, A Manual for Identification. Univ. of Fla. Press, Gainesville, FL. 94 pp.
- Needham, J.G. and M.J. Westfall, Jr. 1955. A Manual of the Dragonflies of North America (Anisoptera). University of California Press, Berkeley. 615 pp.
- Wiggins, G.B. 1977. Larvae of the North American Caddisfly Genera (Trichoptera). University of Toronto Press, Toronto. 401 pp.
- Young, F.N. 1954. The Water Beetles of Florida. Univ. of Fla. Press, Gainesville. 238 pp.
- Epler, J.H. 1992. Identification Manual for the Larval Chironomidae (Diptera) of Florida. Final Report for Florida Department of Environmental Regulation, Contract Number SP251.

METHOD # 3150.2	METHOD NAME TOTAL INORGANIC CARBON (TIC) DISSOLVED INORGANIC CARBON (DIC) <i>APP F</i>	REVISION 1.0	REFERENCE EPA 415.1 EPA-600 4-79-020	STATUS EFFECTIVE
PAGE 1 OF 6				DATE 7/14/95
APPROVALS/DATE				
DIVISION DIRECTOR 	LAB SUPERVISOR <i>DB Swansoff</i>	QUALITY ASSURANCE <i>7-15-95</i> <i>Alvin Swansoff</i>	METHODS <i>DB Swansoff</i>	

SECTION 1.0 - TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
2.0 - METHOD DESCRIPTION/HISTORY	1
3.0 - SAFETY PRACTICES	1-2
4.0 - LIST OF EQUIPMENT/INSTRUMENTATION	2
5.0 - REAGENTS	3
6.0 - STANDARDS	3
7.0 - QUALITY CONTROL	4
8.0 - STEP-BY-STEP PROCEDURE	4-5
9.0 - DATA HANDLING	6
10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING	6
11.0 - REFERENCES	6

SECTION 2.0 - METHOD DESCRIPTION/HISTORY

2.1 Inorganic carbon species in natural waters include gaseous and dissolved CO₂, carbonic acid (H₂CO₃), bicarbonate (HCO₃⁻) and carbonate (CO₃²⁻), occurring in either free forms or as conjugate base (example CaCO₃, NaCO₃, MgCO₃, NaHCO₃, etc). The nature and abundance of individual species depends on the source of water and the underlying soil/sediment material.

2.2 The carbonate species determine the capacity of water to neutralize acid or base (alkalinity or acidity, respectively). Carbon dioxide is a participant in the biological processes of respiration (CO₂ produced), and biosynthesis of autotrophs or photosynthetic organisms (CO₂ consumed).

2.3 Inorganic carbon is determined directly using a high temperature combustion-infrared detection system. Analysis is done on a non-acidified sample, **unfiltered for TIC and filtered** (non-carbonate filter) **for DIC**. In the absence of a set holding time for inorganic carbon, the **holding time for alkalinity (14 days)** is followed.

2.4 Calibration is done using 200 mg/L C (in the form of Na₂CO₃) and concentration of samples is expressed as mg/L C.

SECTION 3.0 - SAFETY PRACTICES

3.1 Wear safety glasses and a full-length, long-sleeved laboratory coat.

3.2 Use acid resistant gloves when handling concentrated phosphoric acid.

3.3 All personnel conducting this method should be familiar with the SFWMD Chemical Hygiene Plan and should have reviewed Material Safety Data Sheets for phosphoric acid and sodium carbonate.

METHOD # 3150.2	METHOD NAME <i>APP F</i> TOTAL INORGANIC CARBON (TIC) DISSOLVED INORGANIC CARBON (DIC)	REVISION 1.0	REFERENCE EPA 415.1 EPA-600 4-79-020	STATUS EFFECTIVE
PAGE 2 OF 6				DATE 7/14/95
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i>	LAB SUPERVISOR <i>DB Swaroff</i>	QUALITY ASSURANCE <i>7-15-95</i> <i>[Signature]</i>	METHODS <i>DB Swaroff</i>	

SECTION 3.0 - SAFETY PRACTICES (CON'T.)


3.4 The electrical power should be disconnected before conducting any repairs inside the instrument on controllers, electrical wiring or any other components near sources of electricity.

3.5 In case of spills of concentrated H_3PO_4 , treat first with an appropriate spill kit and collect the contaminated absorbent and place into adequate storage containers for disposal.

3.6 When changing compressed air bottles, take extra precaution in transporting the bottle to and from the room. It is advisable to request delivery by the vendor to the desired location. Release the pressure slowly in the old tank until the pressure gauge displays 0 psi. Disconnect the regulator from the bottle and cap the bottle tightly with the provided cap. If transporting is required, carefully load the bottle on the hand truck bottle carrier and secure tightly with the strap. Use the same carrier to transport a new compressed air bottle. Secure the new bottle with the strap and connect the regulator. Use teflon tape to avoid any leak. After the regulator is secured, use the SNOOP™ to check for any leaks.

SECTION 4.0 - LIST OF EQUIPMENT/INSTRUMENTATION

- 4.1 DOHRMANN High Temperature TOC Analyzer, Model DC-190, equipped with autosampler and output printer
- 4.2 Volumetric Flasks (Class A): 1000 mL, 100 mL capacity
- 4.3 Volumetric pipets (Class A): 20, 10, 5, and 2 mL
- 4.4 Amber bottle, 1 L
- 4.5 Tensette pipet, 0.1-1.0 mL capacity
- 4.6 8 mL glass vials


METHOD # 3150.2	METHOD NAME <i>UAP F</i> TOTAL INORGANIC CARBON (TIC) DISSOLVED INORGANIC CARBON (DIC)	REVISION 1.0	REFERENCE EPA 415.1 EPA-600 4-79-020	STATUS EFFECTIVE
PAGE 3 OF 6				DATE 7/14/95
APPROVALS/DATE				
DIVISION DIRECTOR 	LAB SUPERVISOR <i>AB Ivanoff</i>	QUALITY ASSURANCE <i>7-15-95</i> <i>Mark Kowalski-Meyer</i>	METHODS <i>AB Ivanoff</i>	

SECTION 5.0 - REAGENTS

- 5.1 Phosphoric acid, concentrated, reagent grade
- 5.2 Carbon-free water (distilled water, preferably double distilled)
- 5.3 Working Phosphoric acid solution (20%): Add 20 mL to 80 mL reagent water. Transfer to the amber borosilicate acid reservoir (located on the right side of the analyzer).
- 5.4 Oxygen 99.99% pure, or purified air (<1 ppm CO₂ or hydrocarbons)

SECTION 6.0 - STANDARDS

- 6.1 Stock Inorganic Carbon Solution : Dissolve 0.883 g of Sodium Carbonate (Na₂CO₃) in carbon-free water and bring to 100 mL final volume, in a Class A volumetric flask. (Prepared by QA Unit monthly).
- 6.2 Calibration Blank (Deionized water)
- 6.3 Calibration standards: 200, 100, 20 mg C/L
 - 200 mg C/L: Dilute 20 mL of stock IC solution and dilute to 100 mL using a Class A volumetric flask.
 - 100 mg C/L: Dilute 10 mL of stock IC solution and dilute to 100 mL using a Class A volumetric flask.
 - 20 mg C/L : Dilute 2 mL of stock IC solution and dilute to 100 mL using a Class A volumetric flask.

METHOD # 3150.2	METHOD NAME <i>APPF</i> TOTAL INORGANIC CARBON (TIC) DISSOLVED INORGANIC CARBON (DIC)	REVISION 1.0	REFERENCE EPA 415.1 EPA-600 4-79-020	STATUS EFFECTIVE
PAGE 4 OF 6				DATE 7/14/95
APPROVALS/DATE				
DIVISION DIRECTOR 	LAB SUPERVISOR <i>DB Ivanoff</i>	QUALITY ASSURANCE <i>Marie K... 7-15-95</i>	METHODS <i>DB Ivanoff</i>	

SECTION 7.0 - QUALITY CONTROL

7.1 Calibration standard recovery must be within 5% of true value.

7.2 QC1 and QC2 will be prepared by the QA unit. QC1 and QC2 are analyzed at beginning of each set of analyses, by the analyst, and are repeated after each 20 samples. Results must be within current QA acceptance limits. In case of unacceptable recoveries, see you supervisor or QA unit staff.

7.3 A repeat sample, selected at random should be analyzed for every 20 samples analyzed. These results are recorded on the QC sheet. The mean and coefficient of variation of the replicate set is determined and noted on the QC sheet. If recoveries are not within the current QA acceptance limits, see your supervisor or QA unit staff.

7.4 A spiked sample must be analyzed for every 20 samples. Add 0.5 mL of stock IC solution, by using a micropipetter, into a 10 mL volumetric flask. Bring to volume with the sample being spiked and mix well. This results in a standard addition value of 50 mg/L C.

7.5 All quality control data must be within the current established limits before entering sample data into the LIMS system. Complete and submit the QC result form.

SECTION 8.0 STEP-BY-STEP PROCEDURE




8.1 Daily Start-up

8.1.1 Check utility supply. Ensure that there is enough gas for a day's operation. Change the air tank when pressure goes below 100 psi.

8.1.2 Check acid reservoir level (located on the right side of the analyzer). Maintain the reservoir at least 1/3 full with 20% Phosphoric Acid Reagent.

8.1.3 Prime IC chamber. With the gas off, fill the IC chamber by using the "Acid to IC chamber" function (Press MAIN, 2, 5). Allow this process to finish, until the flashing message on the screen stops.

8.1.4 Turn on the gas. Ensure that the main valve and the line valves are turned on. Press CARRIER and ensure that the green light is displayed next to the CARRIER button.

METHOD # 3150.2	METHOD NAME TOTAL INORGANIC CARBON (TIC) DISSOLVED INORGANIC CARBON (DIC)	REVISION 1.0	REFERENCE EPA 415.1 EPA-600 4-79-020	STATUS EFFECTIVE
PAGE 5 OF 6				DATE 7/14/95
APPROVALS/DATE				
DIVISION DIRECTOR 	LAB SUPERVISOR DB Swanoff	QUALITY ASSURANCE 7-15-95 	METHODS 	

SECTION 8.0 STEP-BY-STEP PROCEDURE (CON'T)

8.1.5 Press MAIN, 1, to view the status menu. The following default status should be displayed:

Flow rate = 180 - 220 cc/min.
 Dryer temperature = 0 - 10 °C
 Furnace Temperature = 680 °C (Furnace light is green)
 Furnace set point = 680 °C
 Set up number = 1

8.1.6 Check set up. Press IC and make any necessary changes at this point. Ensure that the following settings are entered:

Number of repeats : 2
 Calibration standard: 200 mg/L C calibrant
 Number of rinse: 1 (Wash) and 1 (Sample)

8.2 Tray Protocol

Fill up the glass vials with standards, QC solutions and samples, according to the tray protocol below:

- 1 = D.I. water
- 2 = Blank (D.I. water)
- 3 = Calibration Standard (200 mg/L)
- 4 = Check standard (100 mg/L)
- 5 = Check standard (20 mg/L)
- 6 = QC1
- 7 = QC2
- 8 = Sample #1
- 9 = STAD
- 10-30 = Samples #2-21
- 31 = RPT
- 32 = QC2

Dilute (with D.I. water) any sample with concentrations greater than 200 mg/L.

- 5x : 2 mL sample to 10 mL
- 10x: 1 mL sample to 10 mL

8.3 Put a peg on the outside ring opposite vial #2, and a peg on the inside ring opposite vial #3. Press START when tray is ready. Sampling should automatically initiate and will stop automatically at vial #32. If less than 32 places are used, place a peg on the outside ring opposite the first empty hole.

METHOD # 3150.2	METHOD NAME <i>APP F</i> TOTAL INORGANIC CARBON (TIC) DISSOLVED INORGANIC CARBON (DIC)	REVISION 1.0	REFERENCE EPA 415.1 EPA-600 4-79-020	STATUS EFFECTIVE
PAGE 6 OF 6				DATE 7/14/95
APPROVALS/DATE				
DIVISION DIRECTOR <i>[Signature]</i>	LAB SUPERVISOR <i>AB Ivanoff</i>	QUALITY ASSURANCE <i>[Signature]</i>	METHODS <i>AB Ivanoff</i>	

SECTION 9.0 DATA HANDLING

- 9.1 Follow the instructions given in the LIMS manual to create a workgroup for TIC or DIC.
- 9.2 Review the data printout and complete the INSTRUMENT LOG, QC RESULT and MAINTENANCE LOG form, and ensure that the all the QC parameters are within the acceptable range.
- 9.3 Log into the LIMS, and at ISOTOPE%, type DENTRY. Enter TIC or DIC under parameter, then enter the Workgroup #. Type the value of each data corresponding to the sample number and hit ENTER, after each entry.

SECTION 10.0 - INSTRUMENT CLEANUP/SHUTDOWN/TROUBLESHOOTING

- 10.1 At the end of each run, allow sampling of distilled water twice, to clean the clean the system.
- 10.2 Check the RUN status. The unit should not be in the run mode.
- 10.3 Ensure that the gas is shut off. Carrier gas flow should automatically shut off at the end of each run (Green light next to the CARRIER button is off). Shut off the line valve if the instrument will not be used for the day. Shut off the main valve on the gas tank, if the instrument will not be used for a prolonged period.
- 10.4 Leave the furnace and the NDIR on, unless the unit will not be used for a prolonged period.
- 10.5 Remove used vials, empty the contents on the sink, and dispose of glass vials in designated glass disposal receptacle. Wipe up any spill around the work area and on the instrument.

SECTION 11.0 - REFERENCES

- 11.1 U. S. Environmental Protection Agency. March 1979. Methods for Chemical Analysis of Water and Wastes. EPA-600/4-79-020.
- 11.2 American Public Health Assoc., 1989. Standard Method for the Examination of Water and Wastewater, 17th Edition.
- 11.3 Snoeyink, V.L. and Jenkins, D. Water Chemistry. John Wiley and Sons, NY. 1980.
- 11.4 Dohrmann DC-190 Operator's Manual. 1991.

TIC/DIC ANALYSIS ID-WT

Analyst: _____ WG _____ Date: _____

#	LIMS #	Dilution	Comments	#	LIMS #	Dilution	Comments
1	Wash		D.I.	17			
2	Blank		Calibrate	18			
3	200 mg/L		Calibrate	19			
4	100 mg/L			20			
5	20 mg/L			21			
6			QC1	22			
7			QC2	23			
8				24			
9				25			
10				26			
11				27			
12				28			
13				29			
14				30			
15				31			RPT
16				32			QC2

EVERGLADES SYSTEMS RESEARCH DIVISION STANDARD OPERATING PROCEDURES	Version: 3.0 Revision: 1.1
Everglades Nutrient Removal Project - SOP #2 ENR Sediment/Porewater Monitoring Plan	Date: 08/03/95 FileName: ENRSED01.new

1.0 INTRODUCTION

1.0 PURPOSE AND SCOPE

The Everglades Nutrient Removal (ENR) Project (Figure 1) is a demonstration-scale constructed wetland (1500 ha) designed to reduce total phosphorus (TP) loads in surface water entering the Loxahatchee National Wildlife Refuge. The immobilization of phosphorus (P) within the sediments will be monitored by measuring the relative concentrations of different P forms over time. The field sampling program designed to collect these later data, the ENR Sediment/Porewater Monitoring Plan (SPMP), is described in this Standard Operating Procedure.

1.2 SAMPLING LOCATIONS

A total of 36 sediment/porewater sampling sites will be distributed throughout the ENR Project (Figure 1). Site names reflect which Cell the sites are located in, east/west orientation of the site within the Cell, the sequence number of the site within the transect and, in Cell 3, the type of vegetation the site is located in. Sites will be oriented along the water quality gradient that is anticipated to develop when the wetland is in operation as follows:

1. *Cell 1* - Based on its physical layout and close proximity to the inflow pump station, this Cell will be used to evaluate the peat accumulation and build-up of P in the sediments as a function of distance from the inflow. Porewater wells and core samples will be collected along two transects located 250, 500, 1000, 2000 and 4000 m downstream of the inflow culverts (i.e., a total of 10 sites). Wells are placed at each site at 0-5, 5-10 and 10-30 cm depths. Porewater collection will be quarterly from these wells. Three cores will be collected bi-annually at each site and cut into 0-5, 5-10, and 10-30 cm depth sections. Core sections from the same depth interval at each site will be composited for chemical analyses;
2. *Cells 2 and 4* - The complex shape of Cell 2 and the small size and open-water environment of Cell 4 do not lend themselves to the sampling design specified for Cell 1. Instead, two sites near the inflow and two sites near the outflow will be established for Cells 2 and 4 (i.e., a total of four sites). Wells are placed at each site at 0-5, 5-10 and 10-30 cm depths. Porewater collection will be quarterly from these wells. Three cores will be collected bi-annually at each site and cut into 0-5, 5-10 and 10-30 cm depth sections. Core sections from the same depth interval at each site will be composited for chemical analyses;
3. *Cell 3* - The complex vegetation community established in Cell 3 allows for the opportunity to examine sediment biogeochemistry as affected by different macrophyte species. The macrophyte communities that will be examined are *Eleocharis interstitia*, *Pontederia cordata* and *Sagittaria latifolia*. Three replicate sites will be established at both the inflow and outflow within each of the planting blocks of these species of vegetation (i.e., a total of 18 sites). Wells are

Project Manager Approval: _____

Date: _____

placed at each site at 0-5, 5-10 and 10-30 cm depths. Porewater collection will be quarterly from these wells. Three cores will be collected bi-annually at each site and cut into 0-5, 5-10 and 10-30 cm depth sections. Unlike cores from the other cells, each depth section will be analyzed separately for chemical analysis rather than be composited. This will permit us to perform statistical comparisons of P retention and other sediment parameters among the different vegetation types.

1.3 SAMPLING FREQUENCY

The ENR Sediment/Porewater Monitoring Plan consists of two distinct sampling programs: porewater and soil core analysis. Porewater will be collected quarterly (April, July, October, January). Soil cores will be collected bi-annually (April, July). P fractionation analysis will be done on the soil cores collected in July.

2.0 FIELD AND LABORATORY EQUIPMENT

2.1 ITEM CHECKLIST

2.1.1 Sediment Cores

PVC coring tube (3 in) with spares
 PVC coring attachment
 rubber stoppers
 hammer
 pipe pounder
 wooden pounding block
 knife
 sample labels and rubber bands
 compass
 metric measuring tape
 clipboard
 pencils/pens
 water-proof marking pens
 deionized water
 SPMP chain-of-custody and field data sheets
 field notebook
 ENR Sediment/Porewater Monitoring Plan SOP

2.1.2 Porewater Wells

Peristaltic pump
 marine battery
 sample labels for APA analysis and rubber bands
 60 mL SFWMD bottles
 Contract lab bottles
 0.45µm groundwater filters, 25 mm plain
 light-tight ice chests w/ice
 preservation chemicals (H₂SO₄, HNO₃, NaOH, Zn Acetate)
 deionized water in 5 gallon carboy
 drinking water
 latex disposable gloves
 pH test strips (<2, >9)
 spare parts for wells
 clip board
 pencils/pens/water-proof marking pens
 waders
 SPMP chain-of-custody and field data sheets
 field log book
 ENR Sediment/Porewater Monitoring Plan SOP

2.2 EQUIPMENT LOCATION

1. All field and laboratory equipment is stored in the ESRD equipment lockers located in the Division's offices in B-50 or in the B-113 storage area. Staff should be familiar with the location of all equipment and supplies required for performing this procedure.
2. Deionized water, APA labels, APA bottles, HNO₃ and SPMP chain-of-custody forms are obtained from WRED.
3. After being notified of the field sampling schedule, the contract laboratory will ship sampling bottles and acids in coolers. These coolers will be stored in building B113.

2.3 SPECIAL INSTRUCTIONS

1. Work orders should be delivered to WRE at least two weeks prior to initiating field sampling.
2. Both the contract laboratory and SFWMD laboratory should be notified at least two weeks prior to initiating field sampling and informed of the sampling schedule. Request the contract laboratory to send sampling bottles, coolers and acids (H₂SO₄, NaOH, Zn Acetate) and request the SFWMD laboratory to prepare HNO₃ when sampling porewater.
3. Filters should be ordered at least one month prior to field sampling when collecting porewater.
4. PVC piping and fittings for building the soil coring apparatuses should be ordered at least two months before sampling soil cores. The apparatuses should be assembled at least one month prior to sampling.
5. Arrange for space within the walk-in refrigerator one month prior to sampling soil cores.
6. Arrange for over-night mail pick-up of sample coolers one day prior to sampling porewater.
7. ESRD equipment must be checked out in advance on the equipment sign-out board. All equipment should be tested prior to usage. Report equipment damage or malfunction to the ESRD equipment coordinator.

ENR Sediment/Porewater Monitoring Plan - ENR SOP #2

Version:

3.0

Revision:

1.1

Date:

01/13/95

Page 6 of 24

2.4 EQUIPMENT CALIBRATION

1. Two pH buffers which bracket the analytical range of interest (e.g., 4 and 7) should be used to calibrate the pH meter following the manufacturer's instructions.
2. Meter calibration must be rechecked every four hours while in the field and at the end of the day before leaving the last sampling site.
3. All field calibration records must be kept in the field notebook.

3.0 FIELD AND LABORATORY PROCEDURES

3.1 SEDIMENT CORE COLLECTION

The number of cores that will be collected are as follows:

Cell1 (10 sites):	3 cores/site (Nutrient analysis) 4 cores/site (Nutrient analysis and Phosphorus Fractionation)
Cell2 (4 Sites):	3 cores/site (Nutrient analysis) 4 cores/site (Nutrient analysis and Phosphorus Fractionation)
Cell3 (6 Sites):	9 cores/site (Nutrient analysis) 12 cores/site (Nutrient analysis and Phosphorus Fractionation)
Cell4 (4 Sites):	3 cores/site (Nutrient analysis) 4 cores/site (Nutrient analysis and Phosphorus Fractionation)

3.1.1 Sediment Core Sampling Procedure

- 1.? Immediately upon arriving at the sampling site, insert ? redox rods at ? depth. Allow rods to equilibrate at least 15 minutes before taking a reading. Take redox readings and record on headersheets.
2. Intact sediment cores will be obtained by driving a PVC coring tube (Figure 2) to a depth of approximately 50 cm into the peat. Prior to driving the core tube into the sediment, a knife should be used to cut through any surface roots, plant material, or other detritus at the bottom of the tube. This will help minimize compaction of the core during sampling.
3. Under shallow conditions (i.e., water depth < 30 cm) this will be achieved using a block of wood to protect the neck of the core tube and then striking the block several times with a hammer such that the core penetrates the soil with minimal compaction.
4. If the water depth is greater than 30 cm, a PVC coupling should be attached to the coring device (Figure 2). This attachment is comprised of a ball valve and a length of PVC pipe. This is fitted to the neck of the PVC coring tube and a closed ended piece of metal pipe with a diameter larger than the upright PVC pipe is used to pound the core into the sediment. During the insertion, the ball valve should be open. Prior to pulling the core out of the sediment, the ball valve should be closed.
- 5.? The sediment core, with the overlying water, will then be labeled, capped at both

ENR Sediment/Porewater Monitoring Plan - ENR SOP #2				
Version:	3.0	Revision:	1.1	Date: 01/13/95
				Page 8 of 24

ends using rubber stoppers, stored out of direct sunlight, and transported to the laboratory for analysis within 12?? hrs. Compaction ratios will be recorded on the headersheets. Each sample label will contain the following information:

- a. Station designation,
 - b. Sample replicate information if appropriate
 - c. Date,
 - d. Time, and
 - e. Initials of field personnel collecting sample.
6. If the sediment surface within the core is at a significantly different depth than the adjacent soil, this suggests that compaction has occurred. The core should then be discarded and another one collected.
 - 7.? Triplicate cores will be collected at six randomly selected sites in Cell 3 to determine within-site variability for three different species of vegetation.
 8. Field personnel will complete all appropriate sections of the SPMP sample chain-of-custody form (see Appendix) immediately after collecting a sediment core sample at each station.
 9. **Both crew members will verify that all necessary samples have been collected and the information recorded on the sample label and chain-of-custody form is correct before leaving a station.**
 10. All samples should be delivered to WRED and placed in the walk-in refrigerator. Arrangements for sample pick-up should be made with the contract lab prior to sampling .

3.1.2 Sediment Core Laboratory Preparation

Preparation of soil cores for nutrient analysis:

A. Cell1: There are 10 sites in Cell1. At each of these sites, three cores will be collected. The laboratory thus receives a total of 30 cores from this cell. At the laboratory, each core is cut into 0-5, 5-10 and 10-30 cm depth sections. Core sections from the same depth interval at each site are then composited for chemical analysis (i.e. a total of 30 samples: 10 at 0-5 cm depth, 10 at 5-10 cm depth, 10 at 10-30 cm depth).

B. Cells 2 and 4: Both of these cells have 4 sites each. At each of these sites, three cores are collected. The laboratory thus receives 12 cores from each cell. At the laboratory, each core is cut into 0-5, 5-10 and 10-30 cm depth sections. Core sections from the same depth interval at each site are then composited for chemical

ENR Sediment/Porewater Monitoring Plan - ENR SOP #2						
Version:	3.0	Revision:	1.1	Date:	01/13/95	Page 9 of 24

analysis (i.e. 12 samples for each cell: 4 at 0-5 cm depth, 4 at 5-10 cm depth, 4 at 10-30 cm depth).

C. Cell3: There are 6 sites in Cell 3. At each of these sites, three replicate samples consisting of 3 cores each are collected (a total of 9 cores per site). The laboratory thus receives 54 cores from this cell. At the laboratory, each core is cut into 0-5, 5-10 and 10-30 cm depth sections. Sections from the same depth interval for each replicate at a site are then composited for chemical analysis (i.e. a total of 54 samples: 18 at 0-5 cm depth, 18 at 5-10 cm depth, 18 at 10-30 cm depth).

3.2 POREWATER COLLECTION

Porewater will be collected using porewater wells which will have screens installed to isolate one of the following depths: 0-5 cm, 5-10 cm, or 10-30 cm. Three wells, one for each depth increment, will be required at each site. The choice of this collection technique was based on relatively volume of sample needed for chemical analysis.

Each porewater well is composed of a PVC pipe with a well screen (Figure 3). The size and position of the screen on each well will depend on the desired sampling depth (e.g., 0-5 cm, 5-10 cm, or 10-30 cm depth increments). The end of the well that enters the soil is capped with a pointed PVC fitting while the sampling end is closed off using a rubber stopper. A clear acrylic (or similar inert plastic) tube extends through the rubber stopper to the base of the well. The acrylic tube is closed to the atmosphere by the attachment of a small length of impermeable hose which is closed with an inert stopper. This hose is also used to attach the well to a peristaltic pump. The wells are inserted into the soil and left to equilibrate for two weeks.

3.2.1 Porewater Sampling Procedure

1. Each well is evacuated using a peristaltic pump. The well is then allowed to recharge for approximately 15 minutes or until there is enough volume to meet sample needs. Run the pump at half speed when sampling the 10-30 cm wells. Run pump at full speed when sampling the 0-5 and 5-10 cm wells.
2. After the well has been evacuated and before samples are collected, the pH of the porewater is measured using a compact pH meter.
3. Plastic gloves will be worn when collecting and processing the sample. Because the

analytes of highest priority come from filtered samples, the filtered samples will be collected first. Samples will be collected with a peristaltic pump and will be deposited directly into labeled sample bottles. In-line 0.45 micron filters will be used for filtered samples.

4. Each sample is processed immediately (i.e., filtered and preserved, if necessary), labeled, and stored on ice in a light-tight cooler. The sample preparation, filtering, and handling protocols used for this study will follow procedures listed in the SFWMD's Laboratory Comprehensive Quality Assurance Plan (Section ___) and the PPB Environmental Laboratories sampling instruction sheet. Samples will be transported back to the laboratory within 24 hrs from the time of field collection. Each sample label will contain the following information:
 - a. Station designation,
 - b. Sample replicate information if appropriate
 - c. Date
5. The Equipment Blank (EB) samples for QA analysis will be obtained at the start of each sampling event, after every 20 samples and at the end of the sampling event. To do this, the sampling system will first be purged with deionized water. The sample bottles will then be filled with deionized water and processed in the same manner as the actual samples (i.e., filtered and preserved, if necessary).
6. The split sample (SS) will be collected concurrently along with the sample. This will be done by alternating between the sample and the split sample bottles when filling the bottles. Refer to section 6.2.1 for further information concerning QA/QC frequency.
7. When collecting the field blank (FB) samples, the lids will be kept off the sample bottles until all bottles for that sample are filled. Samples will be acidified immediately before putting on the lids. Refer to section 6.2.1 for further information concerning QA/QC frequency.
8. Field personnel will complete all appropriate sections of the SPMP chain-of-custody and field data form (see Appendix) immediately after collecting a porewater well sample at each station.
9. Both crew members will verify that all necessary samples have been collected and the information recorded on the sample label and chain-of-custody form is correct before leaving a station.
10. Samples analyzed for APA should be delivered to WRED as soon as possible. Samples going to the contract laboratory are placed in a cooler and completely iced. Copies of the SPMP chain-of-custody and field data form are sealed in a zip-lock plastic bag and placed in the cooler. The cooler is then securely taped, addressed

and put in the proper location for over-night pickup.

ENR Sediment/Porewater Monitoring Plan - ENR SOP #2						
Version:	3.0	Revision:	1.1	Date:	01/13/95	Page 12 of 24

4.0 SAMPLE SUBMISSION PROCEDURES

Sediment and porewater samples will be collected by ESRD personnel and subsequently delivered to the contract laboratory for analysis. The sample transfer mechanism will be coordinated with WRED and the contract laboratory. The chain-of-custody sheets will accompany the samples at all times. These sheets will document the transfer of samples between departments within the District and between the District and the contract laboratory.

5.0 HEALTH AND SAFETY

ESRD personnel will observe all standard District health safety policies that pertain to conducting field studies while collecting and analyzing ENR sediment and porewater samples. These policies are detailed in the ENR Safety Manual. Copies of this manual are available from the Project Manager and are stored in the ESRD offices located in B-50. Health and safety items of particular importance for this project include the following:

1. The maximum speed limit is 35 mph on paved roads and 30 mph on unpaved roads and levees. The narrow roadways on the ENR levees would dictate that maximum speeds be kept under 30 mph, especially when towing a boat.
2. Field work shall generally be conducted by crews consisting of at least two people. The only exception to this rule will be road surveys (e.g., reading staff gauges, bird censuses, etc.) which are conducted from a vehicle.
3. No field work will be performed during a thunderstorm or impending lightning. If threatening weather develops, seek shelter in a building or vehicle.
4. Personnel will wear a life jacket at all times while working from or riding in a boat and wading in the marsh.
5. Personnel shall have access to a cellular phone, radio, or other forms of emergency communication at all times while working in the field.
6. Avoid contact with wildlife, especially snakes and alligators.
7. The District's drug-free and smoke-free workplace regulations in buildings and vehicles will be adhered to at all times. Do not smoke in the vicinity of boat gas tanks.
8. All personnel must have received Hazard Communication Training and have read and be familiar with Material Safety Data Sheets for H_2SO_4 , HNO_3 , HCl, pH buffer solutions, and pH electrolyte filling solution.

ENR Sediment/Porewater Monitoring Plan - ENR SOP #2						
Version:	3.0	Revision:	1.1	Date:	01/13/95	Page 14 of 24

6.0 QUALITY ASSURANCE/QUALITY CONTROL

The Project Manager will issue all field and laboratory personnel a copy of the SPMP Standard Operating Procedures and will verify that each staff member has read the SOP. A Quality Control program will be established for the analyses of sediment and porewater. Results of the QC analyses will be reported along with the regular data by the contract laboratory to the Project Manager. Any discrepancies in field or laboratory procedures will be documented by completing a deviation report (see Appendix). Completed deviation reports will be submitted to the Project Manager for review and appropriate action. Sample chain-of-custody forms (see Appendix) will document the transfer of samples between departments within the District and between the District and the contract laboratory.

The Quality Control Program will consist of the following additional samples collected and analyzed for each sampling trip.

6.1 SEDIMENT CORE COLLECTION

1. Triplicate cores collected at six randomly selected sites and three vegetation types in Cell 3.
- 2.???? One known spike prepared and submitted to the contract laboratory by the WRED QA laboratory unit.

6.2 POREWATER COLLECTION

6.2.1 Porewater Wells

1. Equipment blank (EB) samples are generated by passing at least 1 L of deionized water through all field sampling equipment. One EB will be collected prior to field sampling and additional EBs collected (a) after every 20 field samples and (b) at the end of sampling if the total number of field samples is not an even multiple of 20, e.g., if 19 samples are collected then two EBs are required; if 22 samples are collected then three EBs are required.
2. Split samples (SS) are generated by dividing a single field sample into two separate samples. One SS will be collected (a) after every 20 field samples and (b) at the end of sampling if the total number of field samples is not an even multiple of 20, e.g., if 19 samples are collected then one SS is required; if 22 samples are collected then two SSs are required.
3. Field duplicate (FD) samples are generated by collecting an additional sample immediately after the routine field sample using the same procedures. One FD will be collected (a) after every 10 field samples and (b) at the end of sampling if the total number of field samples is not an even multiple of 10, e.g., if 9 samples are collected then one FD is required; if 12 samples are collected then two FDs are required.

4. Field blank (FB) samples are generated by filling a sample bottle with deionized water and allowing that bottle to remain open while collecting the routine field samples. One FB will be collected (a) after every 20 field samples and (b) at the end of sampling if the total number of field samples is not an even multiple of 20, e.g., if 19 samples are collected then one FB is required; if 22 samples are collected then two FBs are required.

Table 1. Chemical constituents to be measured in ENR sediment and porewater samples.

POREWATER (Wells)

alkalinity
 alkaline phosphatase activity
 chloride
 conductivity
 dissolved organic carbon
 total inorganic carbon
 total organic carbon
 total dissolved phosphorus
 ortho-phosphorus
 Al - total
 Ca - dissolved and total
 Fe - dissolved and total
 K - dissolved and total
 Mg - dissolved and total
 Mn - dissolved and total
 Na - dissolved and total
 NH₄
 NO₂
 NO₂+NO₃
 P - hydrolyzable
 pH
 silica
 sulfate
 sulfide
 total dissolved Kjeldahl nitrogen

SEDIMENT NUTRIENT PROFILES (Soil Cores)

%ash
 %moisture
 Al - exchangeable and total
 alkalinity
 bulk density
 Ca - exchangeable and total
 cation exchange capacity
 chloride

Table 1. (continued) Chemical constituents to be measured in ENR sediment and porewater samples.

SEDIMENT NUTRIENT PROFILES (cont.)

- conductivity
- Fe - exchangeable and total
- K - exchangeable and total
- Mg - exchangeable and total
- Mn - exchangeable and total
- Na - exchangeable and total
- pH
- silica
- sulphur
- total carbon
- total organic carbon
- total nitrogen
- total phosphorus

PHOSPHORUS FRACTIONATION (Soil Cores)

- Inorganic Fractions:
 - 1 Molar KCl SRP
 - 0.1 Molar NaOH SRP
 - 0.1 Molar NaOH TP
 - 0.5 Molar HCl SRP
 - Residue TP

- Organic Fractions:
 - NaHCO₃ SRP
 - NaHCO₃ TP
 - Chloroform Treated NaHCO₃ TP




Figure 1. ENR site map showing locations of sampling stations.




Figure 2. Schematic of PVC sediment core device used to collect core samples in the ENR.




Figure 3. Schematic of PVC well used to collect porewater samples in the ENR.

APPENDICES

**EVERGLADES SYSTEMS RESEARCH DIVISION
PROJECT DEVIATION REPORT**

Project Name: _____

Sampling Program: _____

Report Submitted By: _____ Date: _____

Description of Deviation: _____

Probable Impact of Deviation on Data: _____

Corrective Action to Prevent Deviation From Reoccurring in the Future: _____

Project Scientist: _____ Date: _____

Project Manager: _____ Date: _____